VOCO**(GB)****Admira****CE**⁰⁴⁸²**In accordance with DIN EN ISO 4049****Instructions for use****Product description**

Admira is a light-curing restorative material based on Ormocer[®]s and the proven composite technology. **Admira** contains 3-dimensionally linked inorganic-organic co-polymers (Ormocer[®]s) and additive aliphatic and aromatic dimethacrylates. Ormocer[®]s as a material class are protected by patent.

Admira contains 78 % inorganic filler material (= 56 % microfillers by volume, ca. 0,7 µm).

Admira is available in practical syringes and as **Admira Caps** (Composite Application System) for direct application.

Admira is available in 11 shades:	
red-brown shades:	A1, A2, A3, A3.5, A4
yellow shades:	B2, B3
grey shade:	C2
incisal:	I
opaque:	OA2
bleach shade:	BL

Indications:

- Class I to V fillings
- Reconstruction of traumatically affected anteriors
- Veneering of discoloured anteriors
- Correction of shape and shade for improved aesthetic appearance
- Locking, splinting of loose anteriors
- Repair of veneers
- Core build-up under crowns
- Composite inlays

Directions:**Preparation / Shade selection:**

Clean the teeth to be treated with a fluoride-free cleaning paste. Mark occlusal contact points (posterior area); a minimal separation facilitates the design of the approximal contact and placement of the matrix.

Moisten the **Admira** shade guide and match the shade to the moist and clean tooth in daylight before anaesthesia.

Cavity preparation:

In principle, minimal preparation should be carried out according to the rules of the adhesive filling technique in order to conserve as much healthy tooth substance as possible.

Cleaning / Drying:

Cleaning and drying of the cavity is carried out with special preparations or a hydrogen peroxide solution. The work area must be dried before **Admira** is applied. Avoid contamination of the cavity with blood or saliva.

Lining:

A calcium hydroxide lining should be applied in the proximity of the pulp. Deep cavities should be lined. Due to their fluoride release, glass polyalkenoate (glass ionomer) cements are most suitable for this purpose.

Etching / Bonding:

Apply **Admira**, using the "total-etch-technique", and **Admira Bond**. Follow the instructions for use of **Vococid** and **Admira Bond**. Avoid any contamination of the etched enamel or the inhibition layer of **Admira Bond** as it is important for the adhesion to the filling material.

Application of Admira:

Let the material reach room temperature before application.

Apply selected shade of **Admira** (see "shade selection") with a suitable instrument. Use of a translucent matrix is advantageous (smooth surfaces).

Admira Caps:

Insert **Admira Caps** into the opening of the dispenser (follow instructions for use). Turn the cap in the desired direction and remove the protective top. Apply **Admira** directly into the cavity by slowly and evenly pressing the levers of the dispenser together.

Further processing:

Apply and polymerise fillings of more than 2 mm in layers, cure each layer for approx. 40 - 60 s, the opaque shade OA2 for at least 60 s, with a conventional polymerisation device. The light output should be a minimum of 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices.

Transparent strips and crowns etc. should be used. Care should be taken to ensure a good marginal seal. Hold the light emission tip of the light-curing device as close as possible to the surface of the filling. If the distance is more than 5 mm, the curing depth may be reduced and limited to only the area of the light cone. Incomplete curing may lead to discolouration and pulpitis-like complaints.

The filling can be finished and polished immediately after removing the shaping aids (e.g. fine or extrafine diamond burs, polishing disks) with cooling provided. The margin of the filling, or preferably the entire tooth should be fluoridated as a final step.

Admira can be used for direct and indirect inlays, according to established methods. The physical stability values can be improved by the usual external (additional) curing.

Warnings and precautions:

- No known side effects. Hypersensitive persons may experience sensitivity.
- Admira** should not be used if the prescribed application technique cannot be used or if patients have known allergies to individual ingredients (Bis-GMA, diurethane dimethacrylate, BHT, TEGDMA).
- Materials containing thymol, eugenol or oil of cloves interfere with the curing of filling composites. Zinc oxide eugenol cements or other materials containing eugenol should therefore not be used in combination with composites.

Storage:

Store at temperatures between 4 °C and 23 °C. **Admira** may be refrigerated, but must be allowed to reach room temperature before use (Caps at least 15 min).

To avoid exposure to light, syringes should be closed immediately after dispensing.

Do not use after the expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

Ormocer[®] is a registered trademark of Fraunhofer Institut, Würzburg, Germany.

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