- GB** Instructions for use
Aesculap Bone Fixation
- D** Gebrauchsanweisung
Aesculap Bone Fixation
- F** Mode d'emploi
Aesculap Bone Fixation
- E** Instrucciones de manejo
Aesculap Bone Fixation
- I** Istruzioni per l'uso
Aesculap Bone Fixation
- P** Instruções de utilização
Aesculap Bone Fixation
- NL** Gebruiksaanwijzing
Aesculap Bone Fixation
- S** Bruksanvisning
Aesculap Bone Fixation
- RUS** Инструкция по применению
Aesculap Bone Fixation
- CZ** Návod k použití
Aesculap Bone Fixation
- PL** Instrukcja użytkowania
Aesculap Bone Fixation
- SK** Návod na použitie
Aesculap Bone Fixation
- TR** Kullanım Kılavuzu
Aesculap Bone Fixation
- KR** 사용 설명서
Aesculap Bone Fixation



B | BRAUN
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 012829 2018-11 V6 Änd.-Nr. 59339

CE 0123 - DIR 93/42/EEC

Aesculap®

Aesculap Bone Fixation

Legend

- 1 Screw magazine
 - 1a Locked position of the screw magazine (closed)
 - 1b Operation of screw magazine
 - 1c Unlocked position of the screw magazine (open)
- 2 Screwdriver handle
- 3 Silver-colored steel end
- 4 Screw holding sleeve
 - 4a Screw not secured
 - 4b Screw secured
- 5 Screwdriver bit
- 6 Screw

Symbols on product and packages

	Follow the instructions for use
	Not for reuse in intended applications as defined by the manufacturer
	Non-sterile medical product
REF	Manufacturer's article number
QTY	Delivery quantity
LOT	Manufacturer's batch designation
	Manufacturer
	Date of manufacture
	CE-approved in accordance with Directive 93/42/EWG

Intended use

The Aesculap Bone Fixation screws are used for fixing a bone fragment to the jaw bone in implantological dentistry. The Bone Fixation spiral drill and the Bone Fixation screwdriver with bit and screw securing sleeve are provided for implantation of the Bone Fixation screws.

Operating principle

Traction hole-slide hole principle: the functionality of Bone Fixation uses traction screw technology to fix bone blocks.

To this end, a slide hole with a large diameter is created in the bone piece, opposite a traction hole in the transplant mounting with a small diameter. The insertion of a screw ensures the secure adaptation of the transplant on the mounting required for osseointegration.

Available sizes

Set

Art. no.	Designation
DX715	Bone Fixation Set

Uncrew screws

Art. no.	Designation	Slide hole drill	Traction hole drill
DX608T	Magazine with 8 Bone Fixation screws/1.5 mm, L = 8 mm, Phillips head	DX631R	DX630R
DX610T	Magazine with 8 Bone Fixation screws/1.5 mm, L = 10 mm, Phillips head		
DX612T	Magazine with 8 Bone Fixation screws/1.5 mm, L = 12 mm, Phillips head		
DX614T	Magazine with 8 Bone Fixation screws/1.5 mm, L = 14 mm, Phillips head		

Art. no.	Designation	Slide hole drill	Traction hole drill
DX646T	Magazine with 8 Bone Fixation micro-screws/1.2 mm, L = 6 mm, Phillips head	DX630R	DX629R
DX648T	Magazine with 8 Bone Fixation micro-screws/1.2 mm, L = 8 mm, Phillips head		
DX650T	Magazine with 8 Bone Fixation micro-screws/1.2 mm, L = 10 mm, Phillips head		

Drill

Art. no.	Designation
DX629R	Spiral drill, Ø 1.0 mm, L= 31/12 mm, with tool holder diameter 2.35 mm Dental
DX630R	Spiral drill, Ø 1.25 mm, L= 31/16 mm, with tool holder diameter 2.35 mm Dental
DX631R	Spiral drill, Ø 1.5 mm, L= 27/12 mm, with tool holder diameter 2.35 mm Dental

Screwdriver and accessories

Art. no.	Designation
DX604R	Screw holding sleeve for DX604R
DX605R	Screwdriver handle
DX606R	Screwdriver blade

Materials

The materials used in the implant are listed on the label:

- Titanium forged alloy Ti6Al4V acc. to ISO5832-3

The titanium implants are coated with a colored oxide. Slight changes in coloration may occur, but do not affect the implant quality.

Indications

Use for:

- Bone augmentation in the jaw region

Absolute contraindications

- Foreign body sensitivity to the implant materials
- Acute or chronic infections in the region of the implant anchoring
- Severe damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation

Relative contraindications

Do not use if the patient has:

- Fever
- Systemic diseases and metabolic disorders
- Acute osteoporosis or osteomalacia
- Medical or surgical conditions that could negatively affect the outcome of the implantation
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Degenerative bone diseases
- Inadequate patient compliance
- Anticipated overload on the implant
- Cases not listed under indications

Side effects and interactions

- Changes in position, loosening, wear and tear on, or fracture of implant components
- Primary and secondary infections
- Tissue reaction to implant materials
- Injury to blood vessels and nerves
- Haematomas and wound healing disorders

The application or improper use of this system may result in:

- Implant failure resulting from excessive load
 - Warping or bending
 - Loosening
 - Breakage
- Inadequate fixation
- Failed or delayed fusion

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- Aesculap is not responsible for any complications arising from, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- Aesculap Bone Fixation screws may only be combined with spiral drills DX629R to DX631R and screwdriver handle DX605R incl. screwdriver bit DX606R.
- Do not, under any circumstances, use damaged or surgically removed components.
- Implants already used must not be reused.
- The forces to be transferred by the implant must be kept to the achievable minimum.

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

Tissue damage due to heatup of Aesculap Bone Fixation components during MRT examinations!

- ▶ Preferably use low field strengths and sequence techniques with low SAR exposure.
- ▶ Prior to the examination, warn the patient of possible heat sensations and specifically monitor such heat sensations.
- ▶ Abort the examination immediately as soon as heat development occurs.

- Delayed healing can cause implant breakage due to metal fatigue.
- Bone screws once removed from the screw magazine must not be placed back into the screw magazine.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- The surgeon in charge of the therapy shall make any decision with regard to the removal of implant components that have been used.
- Replace the screwdriver blade as soon as it shows functional wear (screws not held by the blade anymore).

Sterility

- The implant components are supplied in unsterile condition.
 - The implant components are in packages of 8 pieces.
 - ▶ Store the implant components in their original packaging and only remove them from their original and protective packaging immediately prior to application.
 - ▶ Use the implant system storage devices for processing, sterilization and sterile setup.
 - ▶ Ensure that the implant components do not come into contact with each other or with instruments in their implant system storage devices.
 - ▶ Ensure that the implant components are not damaged in any way.
- For implant components that are to be resterilized:

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- ▶ Handle the implants with new gloves only.
- ▶ Keep the implant system storage devices covered or closed.
- ▶ Process implant system storage devices separately from instrument trays.
- ▶ Reprocess the implant components individually and separately if no implant system storage devices are available, ensuring that the implant components are not damaged in the process.
- ▶ Carry out the final rinse with distilled, demineralized, or fully desalinated water.
- ▶ Always observe the applicable hospital guidelines concerning the supply of sterile materials.

Disassembling the screwdriver

- ▶ Pull screw securing sleeve 4 from screwdriver bit 5.
- ▶ Retract sleeve on silver-colored steel end 3 and pull screwdriver bit out of the screwdriver handle 2.

Assembling the screwdriver

- ▶ Retract sleeve on silver-colored steel end 3 and insert six-sided end of screwdriver bit 5 into screwdriver handle 2.
- ▶ Turn screwdriver bit 5 until it couples in the screwdriver handle 2.
- ▶ Guide screw securing sleeve 4 completely over screwdriver bit 5.

Validated reprocessing procedure

General safety instructions

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.
The recommended chemistry was used for validation.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for plastics, surgical steel varieties, titanium alloy and pure titanium,
 - do not contain solvents (e.g. acetone or benzene).
- ▶ Do not apply processing procedures involving the use of hydrogen peroxide.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 60 °C.

Note

The Bone Fixation screws are processed in their screw magazines.

- ▶ Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection <ul style="list-style-type: none"> ■ Screwdriver bit DX606R ■ Screw magazine storage DX620 ■ Spiral drills DX629R to DX631R 	<ul style="list-style-type: none"> ■ Cleaning brush, e.g., TA011204 ■ 20 ml disposable syringe ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air 	Chapter Manual cleaning/disinfection and sub-chapter: <ul style="list-style-type: none"> ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection <ul style="list-style-type: none"> ■ Screw holding sleeve DX604R ■ Screwdriver handle DX605R 	<ul style="list-style-type: none"> ■ Cleaning brush, e.g., TA009800 or TA007747 ■ 20 ml disposable syringe ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air 	Chapter Manual cleaning/disinfection and sub-chapter: <ul style="list-style-type: none"> ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection <ul style="list-style-type: none"> ■ Screwdriver bit DX606R ■ Screw magazines DX608T, DX610T, DX612T, DX614T, DX646T, DX648T and DX650T ■ Screw magazine storage DX620 ■ Spiral drills DX629R to DX631R 	<ul style="list-style-type: none"> ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). 	Chapter Mechanical cleaning/disinfecting and sub-chapter:
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection <ul style="list-style-type: none"> ■ Screwdriver handle DX605R 	<ul style="list-style-type: none"> ■ Cleaning brush, e.g., TE654202 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). 	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: <ul style="list-style-type: none"> ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection <ul style="list-style-type: none"> ■ Screw holding sleeve DX604R 	<ul style="list-style-type: none"> ■ Cleaning brush, e.g. TE654202 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). 	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: <ul style="list-style-type: none"> ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
 FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
 RT: Room temperature
 *Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
 FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
 RT: Room temperature
 *Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfecting

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution <ul style="list-style-type: none"> - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water
 FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
 RT: Room temperature
 *Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
 RT: Room temperature
 *Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution <ul style="list-style-type: none"> - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

Note

Not applicable for implants.

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Check for compatibility with associated products.

Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min or as per locally approved policies and procedures
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Application

WARNING

Risk of implant malfunction due to incorrect application!

- ▶ Pre-drill the screw holes, observing the correct burr hole dimensions.
- ▶ Take care that the screwdriver movement causes as little libration (wobble) of the screw as possible when inserting the bone screws.
- ▶ Drive in the screw on the screwdriver axis. Do not tilt the screw.
- ▶ Remove the screws as soon as the healing process is complete.

WARNING

Injury to the digestive tract/dyspnea caused by swallowing/aspiration of bone screws!

- ▶ Use the screw holding sleeve.
- ▶ Take care to ensure that the bone screw does not drop from the screwdriver blade while guiding it into the operating field.

Note

The following warning notices apply for spiral drills DXG29R to DXG31R.

WARNING

Risk of injury and material damage due to inappropriate use of tools!

- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ When coupling/uncoupling, handle tools with cutting edges with care.

WARNING

Risk of burns to skin and tissue caused by blunt tools!

- ▶ Use tools only if they are in perfect condition.
- ▶ Replace blunt tools.

WARNING

Risk of injury to patients and users caused by tools bent or broken during operation!

- ▶ Avoid overstraining (e.g. bending) of tools during operation.

Burr hole dimensions

Screw system	∅ 1.2 mm (gray magazine)	∅ 1.5 mm (yellow magazine)
Pilot hole in the bone fragment/jaw bone	∅ 1.0 mm (blue marking)	∅ 1.25 mm (yellow marking)
Slide hole in the bone fragment	∅ 1.25 mm (yellow marking)	∅ 1.5 mm (red marking)

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- Appropriate clinical training as well as a theoretical and practical command of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.
- Information has been obtained from the manufacturer in the event of an ambiguous preoperative situation and if implants are present in the area to be treated

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The implant components must not be overloaded by extreme strains, hard physical labor, or sports.
- The patient must be informed about the limits of allowable strain on the implant, and be given appropriate behavior guidelines. The risks of transgressing these rules must be explained to the patient. The patient must undergo a regular medical follow-up examination.
- ▶ Pre-drill the bone fragment and the jaw bone with a pilot burr (∅ 1.0 mm or ∅ 1.25 mm).
- ▶ Drill out the pre-drilled bone fragment with a slide hole burr (∅ 1.25 mm or ∅ 1.5 mm).
- ▶ Fixate the bone fragment with bone screws.

Further information on implant systems is available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Technical Service

WARNING

Risk of injury and/or malfunction!

- ▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

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