algistar chromatic FAST SET / algistar chromatic REGULAR SET / algistar classic

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Alginate for dental impression - medical device

Intended use

Alginate for dental impression. Gluten free. Comply with ISO 21563:2013 DENTISTRY – Hydrocolloid impression materials.

Packaging: bag 453 g

Product description:
algistar chromatic *FAST SET*

Fast setting dust free chromatic alginate with chromatic phase indicator. Color variations indicate the respective work step for simple and precise use: violet during mixing, pink during working phase and positioning on the tray, light blue during the time in mouth while setting→ final color: light blue
Indication: creation of study models, situation models, models for orthodontic application.
Aroma: Vanilla-Mint.

algistar chromatic *REGULAR SET*

Dust free normal setting chromatic alginate with chromatic phase indicator. Color variations indicate the respective work step for simple and precise use: dark-orange during mixing, light-orange during the working phase and positioning on the tray, yellow during the time in mouth while setting→ final color: yellow.
Indication: creation of study models, situation models, models for orthodontic application.
Aroma: Tropical fruits

algistar classic

Dust free fast setting alginate.
Indication: creation of study models, situation models, models for orthodontic application.
Aroma: tutti frutti; Colour: orange

Area of application:

For any patient. The device is intended for use by dental professionals such as dentists and dental technicians.

Dosing accessories:

measuring spoons for water and powder.

Mixing ratio:

9,5 g powder : 20 ml water
19 g powder : 40 ml water
28,5 g powder : 60 ml water
38 g powder : 80 ml water

Procedure:

Avoid moisture contamination: open the refill package and transfer its contents into an air-tight container immediately after opening. Before any use, shake appropriately the container to loosen the alginate powder. Moisture and other contaminations lead to incorrect chemical reactions (working times, viscosity) and eventually to unacceptable impressions results. Use only well dry and clean tools (mixing bowl, spatula, mixer, etc.).

Directions:

Dispense the alginate powder from the container using the spoon provided and level the powder in the spoon with a spatula. Place the powder into a mixing bowl and add the corresponding water. Two spoons of alginate and two spoons of water are generally required for a complete impression. Mix the powder and water to an homogeneous consistency using a mixing spatula. Refer to the manufacturer's instructions when using automatic mixers.

Model fabrication:

Best results are achieved by rinsing the impression thoroughly after removal from the mouth, pouring it immediately after cleaning.

Compatible plasters

(ISO 6873:2013): it is recommended the use of the following gypsum products: Whip Mix Microstone -Type 3, Whip Mix Prima-Rock- Type 4, Whip Mix Hard-Rock -Type 5

Impression storage

If longer storage is required, the impression should be kept in a firmly sealed plastic bag (relative humidity of 100 %). This is essential for optimum results. The impression has a dimensional stability (without significant shrinking) of 5 days or 120 hours.

	Mixing time*	Working time* incl. mixing time	Minimum time in the mouth	Initial Setting time*	Total Setting time*
algistar chromatic <i>FAST SET</i>	35 sec.	80 sec.	40 sec.	80 sec.	120 sec.
algistar chromatic <i>REGULAR SET</i>	35 sec.	120 sec.	80 sec.	120 sec.	200 sec.
algistar classic	35 sec.	80 sec.	40 sec.	80 sec.	120 sec.

Times with *) apply from the start of mixing at a temperature of 23°C/73°F for both powder and water with deionized water. These times are reduced by higher temperatures and lengthened by lower temperature. Variations in water hardness may also result in deviations from these times. If tap water is preferred, some practical tests will give the needed information on times.

formation and dispersion in the air. Do not inhale any dust that may be released. Close the containers immediately after use. Ensure that the measuring and mixing area is well-ventilated. Collect any material that may have been spilled with mechanical or suction devices. If mixed with water, the material may form slippery films.

Disinfection:

Impressions may be treated with standard disinfectant solutions. After removal from the mouth and before disinfection thoroughly clean the impression with water, to remove biological residuals. Use a disinfectant intended for alginate impressions disinfecting with short-term contact time (max 10 min). No specific products or brands are recommended.

Procedure: Proceed by immersing in the disinfecting solution for the prescribed time. Alternatively, spray accurately the impression with disinfectant solution and then put the impression in a sealed plastic bag for the prescribed time, then clean with water. As a general recommendation, it is approved the use of water solutions releasing chlorine (for instance: sodium hypochlorites), with a Chlorine concentration of 5000 ppm. The approved procedures are by immersion or by spraying; in both cases with a recommended contact time of 10 min. Please follow the respective manufacturer's instructions for disinfection and processing.

Product storage:

Alginates should be stored in a closed air-tight container in a cool and dry area, protected by excess heat.

Shelf-life:

Our alginates, in the closed original package, have a shelf life of six years from manufacturing date, as printed on the packaging.

Contraindication and directions for patient safety:

Before using, a possible existing allergy potential of the patient should be clarified as well as an allergic reaction to the current product components in the past. In those cases the product must not be used. After removal of the impression, ask the patient to rinse well his mouth and check for any possible residuals left. Taking multiple impressions of the same jaw (more than 3 or 4 in the same session) can cause mucosal irritation.

CLP Safety information for users

The product is not subject to the classification and labeling requirements as it is a medical device in its final state, destined for professional use (EU Regulation No. 1272/2008, Art. 1.5 d). The product is, therefore, also excluded from obligations of information (Material Safety Data Sheet, downstream safety information) as it is excluded from the application of Title IV EU Regulation No. 1907/2006 Art. 2.6.c.

This preparation contains dust that up to 5% may be composed of respirable crystalline silica (Cristobalite, CAS nr. 14464-46-1); prolonged and massive inhalation of crystalline silica respirable dust may cause pulmonary fibrosis. In applying the following instructions, the quantity of dust released during daily dental practice should be such as to ensure the exposure levels below the risk thresholds are maintained. The material is, nevertheless, pretreated with binding agents in order to avoid the release of dust during measuring out and mixing processes.

Precautions and hygiene advice:

When measuring out the material, avoid dust

damage. Wear protective eyewear or mask; pay special attention if you are wearing contact lenses. In case of contact, do not rub or scratch; thoroughly flush with water or eye solutions for 10-15 minutes. Do not swallow. If swallowed, drink a lot of water. The ingestion of limited quantities is not dangerous. The inhalation of alginate powders can be dangerous. Powders can cause irritation or increase susceptibility to respiratory illnesses. Do not breathe the dusts. Avoid generating airborne powders and dusts while dosing and mixing the material. Close containers immediately after use. In case of massive inhalation move the victim to fresh air and facilitate the respiration. Blow the nose, rinse with water the mouth and the throat to eliminate the resting material.

Occupational exposure limits:

The preparation does not contain powdered ingredients for which European Union exposure limits in the workplace exist. However, such limits do exist in some countries in Europe and around the world. For information, the following are some of the occupational exposure limits highlighted:

Exposure limits

Ingredients	Quantity	In some EU countries	In the USA (ACGIH 2013)
Powder		TLV 10 mg/m ³ (inhalable) TLV 3 mg/m ³ (respirable)	
Crystalline silica – Cristobalite	< 5.0% (respirable fraction)	0.05 mg/m ³ or above	TWA 0.025 mg/m ³
Zinc oxide	< 2.0%	2 mg/m ³ (Zn, inhalable) 0.1 mg/m ³ (Zn, respirable)	TWA 2 mg/m ³ (respirable) STEL 10 mg/m ³ (respirable)
Calcium sulfate dihydrate	< 20.0%	4 mg/m ³ (inhalable) 1.5 mg/m ³ (respirable)	TWA 10.0 mg/m ³ (inhalable)
Talc	< 0.5%	0.5 mg/m ³ or above	2.0 mg/m ³ (inhalable)

Disposal information:

Discard damaged containers. Do not mix or use in combination with different products. Avoid contamination and cross-use. Dispose following local laws. Patient-contacted materials are at infection risk (CER cod 18.01.03*).

Waste classification (2014/955/CEE)

Item	Material	CER-Code	Namen
Primary bags and container	PET – ALU - PE	150106	Packaging for waste-mixed packaging
Packages	Paper/cardboard	150101	Packaging for waste – paper/cardboard
Alginate container and contaminated material	Alginates and Fillers	180107	Waste of chemical products and medicines
Dental impressions	Alginates and Fillers	180103*	Waste whose collection and disposal are subject to special requirements to prevent infections

NOTE:

The General Safety and Performance Requirements suggest that any serious incident that led or might led to the death of a patient or other person, to the temporary or permanent serious deterioration of a patient's or other person's health, or to a public health threat, should be reported to the National Competent Authority and to the manufacturer.

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Müller-Omicron GmbH & Co. KG