CLEARsplint® CAD/CAM Pro

astron®

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Advantages of astron® CLEARsplint® CAD/CAM Pro that speak for themselves!

- extremely break resistant, hence optimal protection of the teeth
- thermoactive flexibility which provides greater, comfortable and tension-free wearing comfort
- hypoallergenic, free of MMA (residual monomer content of MMA = 0%)
- self-adjusting. This superb accuracy permits the closest fit
- always transparent – will not yellow with age, since amine-free

astron® CLEARsplint® CAD/CAM Pro

The best solution for indications as easing the mandibular joint (bruxism), for the protection of implants, the adjustment of the occlusion etc.

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
PMMA blanks included in the Starter Kits serve as base blank and, for the first time, blanks have to be milled out and refilled with CLEARsplint® CAD/CAM Pro. Then the blank can be refilled repeatedly with CLEARsplint® CAD/CAM Pro until it wears out.

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Anatomical model with exactly fitting bite protection

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thermoactive bite protection self-adjusting and very flexible

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A deformed CLEARsplint® comfort splint adapts easily and precisely to the teeth thanks to the thermoactive property of the material and its excellent flexibility. Even under heat, the perfect surface hardness is maintained.

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Store powder and monomer in the refrigerator at 4–8 °C. (Practical tips: prior to use put mixing bowl into the refrigerator for 30 minutes to extend processing time – we recommend the use of ceramic- or glass bowls to avoid any chemical bond with the material). Thus working time will be extended of this heat-sensitive material.

In order to avoid uncontrollable chemical reactions please use only mixing bowls and -spatulas which do not get in contact with other resins. Processing repairs as usual. Prior to processing repairs, put some CLEARsplint® CAD/CAM Pro liquid on area to be repaired and let it dry for some seconds.

ATTENTION:

Please respect mixing ratio, polymerisation temperature, -pressure and -time thoroughly! Otherwise material risks not to cure properly or product-specific advantages may be influenced negatively!

- choice of the appropriate position and size of the area to be milled necessary for the filling with CLEARsplint® CAD/CAM Pro in the device-specific blank (base blank)

ATTENTION:

Only use blanks which are included in the Starter Kit or device-specific appropriate transparent blanks recommended for the manufacture of splints! Please make sure that the CLEARsplint® CAD/CAM Pro resin is well bonded to the base blank! A blank can be refilled repeatedly since material adheres to the residual material of the last filling (Reparability).

Mixing ratio: 1 g powder to 0.7 g monomer

CREATE THE BLANK:

For one splint you need according to size of the milled-out space approx. 30 to 40 g powder. This means you need for 30 g (0.07 lb) powder 21 g (0.05 lb) monomer, resp. for 40 g (0.09 lb) powder 28 g (0.06 lb) monomer.

- Isolate base plate with a thin layer of vaseline
- fix milled blank with the clamps onto the CLEARsplint® CAD/CAM Pro base (base plate)
- take powder and monomer out of the refrigerator just immediately prior to use. Mix thoroughly the appropriate amount for approx. 5–10 seconds as quick as possible and bubble-free
- fill the milled-out space as quickly as possible with the CLEARsplint® CAD/CAM Pro resin and wait until resin forms a membrane (start of the surface polymerisation). This is absolutely necessary to avoid water inclusions during polymerisation!
- Immerse the filled and fixed blank in the polymerisation pot. Make sure that the complete blank is covered with water
- polymerisation in 50–53 °C (122–127 °F) warm water in a pressure pot at 2 bar (29 PSI) for 45 minutes
- remove CLEARsplint® CAD/CAM Pro base from pressure pot and let it cool down under cold running tap water
- then clamp blank into the milling unit and start milling process of your CLEARsplint® CAD/CAM Pro splint
- splints made of CLEARsplint® CAD/CAM Pro resin can be trimmed and polished as usual

HINT:

We recommend to water the finished splint for approx. 3–4 hours, then to rinse and brush it under running water with a soft brush and to disinfect it afterwards with an alcohol-free disinfectant.

ATTENTION:

Please write down the lot number of the used resin, as well as the date of the manufacturing and the name of the patient. For the care of the splint we recommend to use DuraClean², a cleaner for splints especially developed from the astron company. Do not use denture cleaners, as well as abrasive tooth pastes, since they may damage the splints even up to complete destruction of the CLEARsplint® CAD/CAM Pro „comfort splint“.

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KENTZLER-KASCHNER DENTAL GmbH
Astron Dental Corp.