Celtra® Duo

CAD Blocks for CEREC® and INLAB®

DIRECTIONS FOR USE – ENGLISH

For dental use only.

EN

1. PRODUCT DESCRIPTION

Celtra® Duo Zirconia-Reinforced Lithium Silicate (ZLS) block for CEREC® and InLab® is a high-strength glass-ceramic material that, thanks to its translucent properties, can be used for the fabrication of highly esthetic all-ceramic restorations using a CAD/CAM procedure. Indirect restorations are fabricated by milling Celtra Duo using a Dentsply Sirona CAD/CAM system. Celtra Duo is a zirconia-reinforced lithium silicate (ZLS) block delivered in a fully crystallized, tooth colored state. It provides two different processing pathways where the material can be milled, polished and seated for a flexural strength of 210 MPa using a three-point bending test method; bi-axial test method is 413 MPa.

Alternatively, a heat firing cycle will increase the material's flexural strength to 370 MPa using a three-point bending test method; bi-axial test method is 560MPa. The use of stains and glaze is optional. In clinical situations where maximum strength is needed, it is recommended to fire the restoration. If the mill and polish option is chosen, all polishing steps must be followed to achieve the increase in strength; see section 3.3.1

1.1 Composition

Zirconia-reinforced lithium silicate (ZLS), type II dental ceramic material, Class 2a pursuant to DIN EN ISO 6872 (CTE 25–500°C: $11.6 \cdot 10^{-6} \text{ K}^{-1}$).

1.2 Indications

Celtra Duo is a fully crystalized zirconia-reinforced lithium silicate glass ceramic processed through CAD/CAM technique for the fabrication of:

- Single unit dental restorations, for example all-ceramic crowns, inlays, onlays and veneers.

Celtra Duo can be used as a material for fully contoured applications.

1.3 Contraindications

- Full veneers on molar crowns
- Very deep subgingival preparations
- Patients with a greatly reduced dentition
- Patients with bruxism

1.4 Compatible porcelains, stains and glazes

The use of stains and glaze is optional with Celtra Duo. Celtra Duo restorations are compatible with Dentsply Sirona's Universal Stains and Glaze System or Celtra Stains and Glaze System (available separately, see complete Directions for Use). The use of other stain and glaze systems is not recommended (see Warnings). Celtra Duo is compatible with and thermally matched to Celtra Correction Porcelain for desired contour adjustments (if necessary). Use of other correction porcelain with Celtra Duo is at the discretion and sole responsibility of the user.

1.5 Compatible luting cements

Celtra Duo restorations are compatible with adhesive bonding resin cement systems, including all Dentsply Sirona adhesive resin cement systems (available Separately (see complete Directions for Use)).

It is our recommendation that Celtra Duo be adhesively bonded regardless of indication. Full coverage crown restorations that have a retentive prep design and have been fired are compatible with self-adhesive resin cements, including all Dentsply Sirona self-adhesive resin cements (available separately, see complete Directions for Use). Inlays, onlays, non-retentive crowns should be adhesively bonded. Use of other cements or cement systems with Celtra Duo is at the discretion and sole responsibility of the dental practitioner.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other sections of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety message that follow this symbol to avoid possible injury.

2.1 Warnings

- If properly processed and used, adverse effects of these medical products are highly unlikely. However, reactions of the immune system (such as allergies) or localized paresthesia (such as an irritating taste or irritation of the oral mucosa) cannot be completely ruled out as a matter of principle. In case of skin sensitization

or rash, discontinue use and seek medical attention.

- Celtra Duo restorations are not suitable for patients with clinical symptoms of parafunctional habits or bruxism (see Contra-indications).
- Do not inhale dust particles during grinding. Wear suitable protective mask.
- Do not use ceramic pins, firing pastes, or any third-party porcelains or stains and glazes. Such use may result in cracks in Celtra Duo restorations.
- In patients with hypersensitivity to any of the ingredients, this medical device may not be used at all or only under the particular scrutiny of the dentist or physician in charge.

2.2 Precautions

- This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and is the sole responsibility of the practitioner.
- Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
- Contact with saliva, blood and/or some astringent solutions during adhesive procedures may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Celtra Duo restorations require adequate preparation reduction and restoration thickness. Insufficient wall thickness may lead to premature failure.
- Celtra Duo restorations are designed to be completely polished and/or fired before insertion. Direct insertion without following all of the polishing steps and/or firing may lead to failure.
- Use only in well ventilated areas.

2.3 Interactions

- Celtra Duo blocks are designed to be fabricated using a Dentsply Sirona CAD/CAM system. Cutting blocks using non-compatible CAD/CAM system may lead to inadequate or unacceptable restorations.

2.4 Adverse reactions

No adverse reactions have been reported for Celtra Duo. Should you hear or receive information about any adverse effects – even when doubtful – please notify Dentsply Sirona.

2.5. Storage conditions


Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Store in a dry place. Protect from moisture. Do not use after the expiration date.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Preparation

Proper reduction of the hard tissue of the tooth during preparation is essential for maximizing the strength, shade and retention of the finished restoration. When preparing anterior or posterior teeth, the anatomical form has to be reduced as shown.

Minimum wall thickness: The following diagram shows the specified minimum wall thickness for each indication. The minimum wall thickness must still be ensured after all manual adjustments have been made.

* 1.5-2mm for frameworks that are exclusively polished.

3.1.1 Inlays and onlays

- A conventional inlay/onlay design is recommended. Do not prepare undercuts. Ensure that the cavity walls form an angle of 5–6 degrees with the long axis of the tooth. Ensure that all sharp edges and angles are rounded. In centric and dynamic occlusion, reduce incisally/occlusally by 1.5 to 2 mm.
- Celtra Duo inlays and onlays are ideally delivered by fully adhesive cementing. Alternatively, highly retentive inlay or onlay restorations may be cemented with self-adhesive resin cements. See complete Directions for Use of the adhesive and cement products.

3.1.2 Full crowns

- Ensure that there is an axial reduction of 1.0 to 1.5 mm with the walls form an angle of 5–6 degrees with the long axis of the tooth. In centric and dynamic occlusion, reduce incisally/occlusally by 1.5 mm. The lingual shoulders must be extended at least 1.0 mm into the proximal contacts surfaces. It is recommended to use a shoulder preparation without a bevel: All angles must be rounded, and the preparation surfaces must be smooth.
- Celtra Duo full crowns can be delivered by fully adhesive or self-adhesive cementing. See complete Directions for Use of the adhesive and cement products.

3.1.3 Veneers

- The standard reduction is 0.6 mm for the labial surface and 0.4 mm in the gingival area (since the enamel is thinner in this region). Reduce the labiolingual incisal angle by 1.0 to 1.5 mm. The preparation margins should be located in enamel. A chamfer or rounded-shoulder preparation is recommended for all veneer margins. Proximal extensions must be located far enough proximally to conceal preparation margins from site and to avoid proximal gingival undercuts.
- Celtra Duo veneers are delivered by fully adhesive cementing. See complete Directions for Use of the adhesive and cement products. Self-adhesive cementing is not recommended for veneer restorations.

3.2 CAD/CAM processing

Celtra Duo restorations are produced with CEREC or inLab CAD/CAM systems by Dentsply Sirona. If you have any questions about these systems, please contact Dentsply Sirona Dental Systems GmbH.

3.2.1 Software requirements

- Celtra Duo is supported by CEREC version 4.2 and above. When using previous software versions, use of the IPS e.max® system standard grinding program is recommended (do not use the e.max® system speed program).
- In this case, the minimum wall thickness will have to be adapted to the Celtra Duo recommendation; see Preparation.

3.2.2 Grinder and block sizes

- In the CEREC software, select the Celtra Duo material. The Dentsply Sirona Celtra CAD/CAM grinding unit (milling unit) will prompt you to insert a Celtra Duo Block. For detailed processing steps, please consult the Instructions for Use and technical manuals of the appropriate CAD/CAM systems. Make sure to follow the manufacturers' recommendations.

3.3 Finishing the restoration

Celtra Duo provides two processing pathways where the material can be milled and polished or milled, glazed and fired for additional strength. Alternatively, the restoration can also be polished and fired for added strength. For polishing instructions go to section 3/3/1, and stain and glaze instructions go to section 3.3.3.

3.3.1 Polishing without the optional staining and/or glaze firing

- Clean the restoration in an ultrasonic water bath, or with a steam cleaner or clean the restoration using soap and water. Carefully blow-dry the restoration in an air stream.
- Try in the final restoration to examine its fit. If necessary, adjust the proximal and occlusal contacts.
- When making adjustments to the ceramic material, **make sure to avoid overheating**.
- Sprue removal should be made using a suitable grinding instrument, between 8,000-12,000 rpm's and using light pressure with a laboratory motor and hand piece. After sprue removal, use a coarse rubber wheel to remove any bur marks using a lab motor at a speed of between 8,000-12,000 rpm's with light to medium pressure.
- Diamond polishing bodies (< 60 µm) are recommended for polishing the occlusal surfaces using a lab motor at a speed between 8,000-12,000 rpm's with light to medium pressure. Do Not Overheat. Follow the steps below for complete polishing:
 - In clinical situations where maximum strength is needed, we recommend firing the restoration
 - The following steps must be followed in order to achieve the healing effect and increase in strength to 210MPa (three point bending test method) for the Mill and Polish option
 - Use a coarse grit to prepare the surface for polishing and to create a more uniform surface
 - Followed by a medium and then a fine wheel (8,000-12,000 rpm) using light to medium pressure
 - An optional final step is to go over the entire restoration with fine diamond paste using a soft-medium Robinson brush (4,000-8,000 rpm) with light pressure
- Celtra Duo quickly becomes shiny, there is a risk to under polish the material and not achieve the strength of 210 MPa. It is not the shiny surface that gives the material its strength, it is the effect of the smoothing and polishing that strengthens the material.

3.3.2 Staining, Glazing and Applying Corrections to milled Celtra Duo restorations

- If milled restorations require no corrections, proceed to step 3.3.3 for optional staining and glazing instructions.
- If milled restorations require corrections, proceed to step 3.3.4 for Application of Celtra Correction Porcelain instructions.

3.3.3 Instructions for optional staining and glaze-firing of Celtra Duo CAD/CAM restorations

- Custom staining of Celtra Duo restorations can be performed using Celtra Universal Stains and Glaze or Dentsply Sirona Universal Stains and Glaze instructions below.
- For best results, any sanding marks on the surface should be removed with a suitable coarse rubber wheel
- Before applying the stains and glaze, the restoration must be clean and free of grease. Clean the surface of the restoration with a steam cleaner or submerge for 10 minutes in an ultrasonic cleaner with distilled water. If a steam cleaner or ultrasonic

cleaner isn't available clean thoroughly with soap and water. Any contamination after cleaning must be prevented.

- Mix the mass thoroughly with a spatula. Do not use a metal spatula. Withdraw desired amount of Celtra stain or glaze and place it on the mixing palette. For a thinner consistency dilute the material with Dentsply Sirona Stain and Glaze Liquid.
- Use Dentsply Sirona stain and glaze liquid to clean the brush, do not use water to clean the brush as water can create an "orange peel" appearance.
- Apply a sufficient amount of glaze to the entire crown surface using a brush. Apply a thin layer of stain in the gingival area and spread with a brush toward the incisal region (or as needed). Check the shade using the shade tab. Use the brush to increase or decrease the stain amount to achieve the desired shade.
- Again using the brush, apply a thin layer of incisal stain to the incisal edge. Check the shade using the shade ring. Use the brush to increase or decrease the stain amount to achieve the desired esthetic incisal effects.

Notes:

- A more intensive shade effect can be achieved by repeating cycles of applying and firing the material. However, the use of too many layers of stain may result in an unnatural appearance.
- The cusps and fissures can be individually characterized with stains.
- The basic shade is determined based on the Vita shade groups (A, B, C, D) (see Table 1). Please use the Vita classical shade guide.

Table 1: Assignment of Shade Groups

Vita Shade Group	Bleach	A	B	C	D
Universal Stain Shade	Stain 0	Stain 1	Stain 2	Stain 3	Stain 4

Table 2: Use of Incisal Stains

Stain	Tooth Shade
Incisal Stain i1	A1, A2, A3, B1, B2, B3, B4, D2, D3, Bleach
Incisal Stain i2	A3.5, A4, C1, C2, C3, C4, D4

Use of glazes and stains by third-party manufacturers

Only the native range of glaze and stains (Dentsply Sirona Universal Stain & Glaze) in conjunction with the appropriate liquid is recommended for customizing and glazing Celtra Duo, as its compatibility with Celtra Duo has been tested and is guaranteed. Glazes and stains by third-party manufacturers may require other firing temperature and cycles than those recommended for Celtra Duo.

3.3.4 Instructions for Application of Celtra Correction Porcelain to Celtra Duo CAD/CAM restorations

- Mix correction porcelain with distilled water or Dentsply Sirona Modeling Liquid U to desired consistency (thick consistency). Apply porcelain where applicable. Fire to complete the restoration following the procedures and firing programs listed below:
- First Glaze and Correction Porcelain: Apply Dentsply Sirona Universal Glaze over the entire Celtra crown then apply CELTRA Correction porcelain to the required areas and fire at 820 °C according to the first correction firing program in Table 3. Polish to a high gloss, if needed.
- Second Glaze: If after first firing the required glaze appearance is not achieved, re-apply glaze, and fire according to 2nd and subsequent glaze firing program in Table 3.

3.4 Firing the restoration

3.4.1 Firing Universal Stains and Glaze Only

- Place the restoration on a firing pad and then on a honeycomb tray. Always use a firing pad for all restorations, placed onto a honeycomb sagger tray.
- If firing an anterior, place the restoration with the lingual facing down and with premolars, place the restoration with the interproximal facing down. If there are any imperfections left of the glaze from the firing pad, it can be easily corrected by polishing that area.

3.4.2 General firing recommendations

- The optional firing will increase the flexural strength of Celtra Duo to 370 MPa (three point bending test method). To increase the flexural strength it is only necessary to fire the restoration up to 820 °C. The application of glaze is an option but not necessary. In case of firing without glaze skip the Pre-drying and drying steps of the firing program. Make sure to completely polish the restoration first before firing.
- Starting temperature 500 °C
- Place the object on a firing pad then on honeycomb tray, then place on the firing tray/firing table of the furnace.
- Additional glaze firings may be performed at 770 °C in order to accentuate the shade or to correct it with glaze, or increase the gloss. It is necessary to coat the entire surface with Dentsply Sirona Universal Glaze (available separately, see complete Directions for Use) to obtain a uniform glossy finish.
- Check the shade against a shade tab and adjust if necessary.

- Note: After the first firing, the restoration may not appear sufficiently glazed. In this case, apply a thin layer of Dentsply Sirona Universal glaze to the entire surface with a brush. After applying the glaze, perform a second stain/glaze firing as defined in the firing recommendations (see table 3). Additional stain/glaze firing cycles can be performed using the same parameters.

Table 3: Firing Charts

Celtra Duo Firing Recommendations for the Programat CS/CS2, Programat EP 3000/5000, P 500								
	Standby Temperature B	Closing time S	Heating rate t ↑	Firing temperature T	Holding time H	Vacuum on/off	Long-term Cooling L	Cooling temperature tL
	°C	min	°C/min	°C	min	Vac 1 (°C) Vac 2 (°C)	°C	°C
Celtra Paint-On Glaze	500	3:30	60	820	1:00	off	750	50
2nd & Subsequent Glaze Firing - if needed	500	3:30	60	770	1:00	off	750	50
Spray Glaze (Indenco)	500	2:00	60	820	1:00	off	750	50
Polish and Fire	500	1:00	60	820	1:00	off	750	50
Correction Porcelain	500	3:30	60	820	1:00	V1= 500°C V2= 819°C	750	50

Celtra Duo Firing Recommendations for VITA ovens							
	Start Temperature	Pre-drying	Heating rate	End temperature	Holding time	Long-term Cooling	VAC
	°C	min	°C/min	°C	min	min	min
Celtra Paint-On Glaze	500	4:00	60	820	1:00	3:00	-
2nd & Subsequent Glaze Firing - if needed	500	4:00	60	770	1:00	3:00	-
Spray Glaze (Indenco)	500	3:00	60	820	1:00	3:00	-
Polish and Fire	500	2:00	60	820	1:00	3:00	-
Correction Porcelain	500	4:00	60	820	1:30	3:00	1:30

Note: Please visit Celtra-DentsplySirona.com for other furnace firing programs

3.4.3 Processing Celtra Duo in the Dentsply Sirona SpeedFire Oven

- Make sure to select the Celtra Duo material under Dentsply in the CEREC Software prior to scanning as this will ensure the correct firing program will be transferred.
- Place Celtra Duo directly on the round firing pad; if using a square firing pad, it will need to be trimmed to fit on the upper door insulation
- Select the Celtra Duo job press and start button on the display screen
- Please not that due to the hardware, the correct firing temperature of the CEREC SpeedFire is below the temperature listed in this document for conventional furnaces.
- If firing an anterior, place the restoration with the lingual facing down and with premolars, place the restoration with the interproximal facing down. If there are any imperfections left of the glaze from the firing pad, it can be easily corrected by polishing that area.

3.5 Cementing

3.5.1 Preparation of the Celtra restoration

- Clean the restoration with an ultrasonic or steam cleaner or with alcohol.
- Apply 5%-9% hydrofluoric acid etching gel (Available separately, see manufacturer's complete Directions for Use) to the interior of the restoration only. Etching time 30 seconds.
- **CAUTION:** Follow manufacturer's precautions. Do not allow tissue or eyes to come into contact with the acid! Remove the hydrofluoric acid as per the manufacturer's instructions.
- Dry the restoration in an air stream.
- It is recommended to silanate the etched surfaces immediately.
- At chairside, apply silane only to those surfaces required for adhesive cementing.
- Allow to soak for 60 seconds. If the silanize layer is no longer liquid, add more silane. Blow-dry in a powerful air stream. (Recommended material: Calibra® Silane Coupling Agent, Available separately, see complete Directions for Use).

3.5.2 Cementation

- Depending on the indication, Celtra Duo restorations can be self-adhesively or fully adhesively cemented. Compatible time-proven adhesive cementing materials are available as part of the Dentsply Sirona range of products. (Recommended material: Calibra® Ceram for fully adhesive and Calibra® Universal for self-adhesive Cements are available separately. See complete Directions for Use).

	Self-adhesive	Fully adhesive
Inlays	R	HR
Onlays	R	HR
Crowns	R	HR
Veneers		HR

R = recommended

HR = highly recommended

4. HYGIENE

CAUTION

Cross-contamination

Do not reuse single use products. Dispose of in accordance with local regulations.

5. LOT NUMBER AND EXPIRATION DATE

5.1 Do not use after expiration date. ISO standard uses: "YYYY-MM" or "YYYY-MM-DD"

5.2 The following numbers should be quoted in all correspondences:

- Reorder number
- Lot number
- Expiration date

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