Durafil® VS

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GB Instructions for use

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Manufacturer:
Kulzer GmbH
Luisenpark Straße 2
63450 Hanau (Germany)

Made in Germany

Distributed in USA / Canada exclusively by:
Kulzer LLC
4215 South Lafayette Blvd.
South Bend, IN 46614-2517
1-800-431-1785

Faqs:
Durafil® VS restricts this device to sale by or on the order of a dentist.

Importado e Distribuído no Brasil por:
Kulzer South America Ltd.
Rua Curitiba, nº 1000-02
R.uno Sena Strip, 27 - cj. 42
São Paulo - SP - CEP 05036-010
Braziliense Distribuidora S/A
Rua: Teófilo Otonari, 100
Res. Técnica: Dra. Regiane Marton -
CRP 70.705
N°ANVISA- vide embalagem

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Durafil® VS Instructions for use (GB)

Light curing microfilled composite for anterior restorations
Durafil® VS is a microfilled lightcured composite filling material with prepolymers in form of splinters. Durafil® VS may be finely polished and in connection with an adequate bonding system it is meant for adhesive filling therapy.

Composition: BisGMA, TEGDMA, UDMA, Triandiol, fluorescent pigments, metallic oxide pigments, organic pigments, ammoniumzinc sulfate, camphorquinone, BHT. Durafil® VS contains approximately 66% filler by volume (with 30% inorganic filler by volume), its filler particle size is 0.02 – 0.07 µm.

Indications: For fillings of cavities class II, IV in connection with the adhesive technique. Not intended for indications or uses that are not cited in these instructions for use.

Contraindications: Application to dentine close to the pulp (less than 1 mm). Use of this product is contraindicated in case of known or suspected allergies against components of this product. Do not use this product on known or suspected allergies against methacrylate compounds.

Warnings and safety instructions: This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer. May be irritating to eyes and skin. Avoid contact with eyes. May cause sensitization by skin contact. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with clothes. In case of contact with wear, rinse with water for several minutes. If symptoms persist, consult an eye specialist.

Application hints: Visually inspect for damage before use. Damaged products must not be used. Composites should only be used at room temperature. In the case of deep preparations in which only a minimal dentin exposure cover the dentine or direct pulp exposure cover the dentine or direct pulp exposure with a glass ionomer cement before applying the adhesive. Do not use eugenol containing materials in or around the cavity preparation being restored, because these may interfere with the proper curing of the composite. Use adequate protection for splinters to prevent from any contamination. Do not re-use after contamination. The product is not intended for single use only and reprocessing or re-use are not possible without damaging and contaminating the restorative material by those phenolic disinfectants solutions. To avoid cross contamination PLS must be disposed of immediately after each use. When using

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Indications: AT, AZ, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, C4, D3, D4, D2, O42, O43, O43.5, Incisal, SL, SSL.

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Durafil® VS PLS, wie eine gemille, even motion when exerting pressure on the applicator gun. Use of excessive force or sudden movement could result in unanticipated extrusion of the material or ejection of the PLS from the applicator gun, creating a potential hazard. We recommend using the Kulzer PLS Applikator STEEL. When using other PLS applicator guns see full instructions and precautions of these devices. Replace the cap of the syringe after every use.

Clinical treatment protocol: Select an appropriate shade or shade combination from the Durafil® VS Shade guide prior to placing the burr. This is important because teeth will significantly change color when desiccated. If the tooth surface being treated is not to be instrumented by a diamond stone, carborundum, drill or other abrasive instrument, the surface should be cleaned with pumice or a fluoride free polishing paste before restoration is placed. A non-sticking plastic, Teflon coated or gold plated plastic instrument (such as the Kulzer Pflmasarc Instrument) is recommended for use in clinically placing, sculpting and shaping the material.

1. For ideal results, the use of an ivory rubber dam (or similar) is recommended to control contamination from moisture, blood and saliva.

2. Prepare the cavity according to generally accepted restorative techniques. Then clean the cavity with the adhesive and prime the surface of the cavity.

3. Prior to application of Durafil® VS restorative material to the cavity preparation, pre-treatment of the tooth structure with an adequate enamel/dentine adhesive is mandatory (e.g. BOND® Universal, BOND® Total Etch by Kulzer). Once treated with the adhesive, avoid contamination with saliva or blood.

4. Apply Durafil® VS in thin layers (max. 2 mm) according to the following table with a fine-grained diamond stone or similar, or a light-curing LED light-curing unit (wavelength maximum at 460-470 nm, light output > 400 mW/cm²). The surface inhibition layer that forms after curing of Durafil® VS should not be removed or touched during the build-up phase of the restoration. It serves as chemical protection for any subsequent composite layers that are added.

Recommended Durafil® VS curing times with Translux® Curing Unit or other OHT light curing units:
Durafil® VS Shade Guide

Colour	Time of irradiation (s)	Layer thickness (mm) Translux®
A1, AZ, A3, B1, B2, B3, C1, C2, C3, D3, Incisal, SL, SSL	20	2.0
A3.5, A4, C4, O42, O43, O43.5	40	2.0

Hard Dentist Substance (mm)

Time of irradiation (s)	Layer thickness (mm) Translux®
20	2.0
40	2.0

• Curing below hard dent substance (all colours)

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