AESCULAP®

en	Instructions for use/Technical description Reusable filter
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusaifus.com . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung Dauerfilter
fr	Mode d'emploi/Description technique Filtre permanent
es	Instrucciones de manejo/Descripción técnica Filtro reutilizable
it	Istruzioni per l'uso/Descrizione tecnica Filtri permanenti
pt	Instruções de utilização/Descrição técnica Filtros permanentes
nl	Gebruiksaanwijzing/Technische beschrijving Permanente filters
da	Brugsanvisning/Teknisk beskrivelse Permanentfiltere
sv	Bruksanvisning/Teknisk beskrivning Permanentfilter
fi	Käyttöohje/Tekninen kuvaus Kestosuodatin
lv	Lietošanas instrukcijas/tehniskais apraksts Ilgstoši izmantojams filtrs
lt	Naudojimo instrukcija/techninis aprašas Ilgalaičio naudojimo filtras
ru	Инструкция по применению/Техническое описание Постоянный фильтр
cs	Návod k použití/Technický popis Permanentní filtr
pl	Instrukcja użytkowania/Opis techniczny Filtr trwały
sk	Návod na použitie/Technický opis Trvanlivé filtre
hu	Használati útmutató/Műszaki leírás Tartós szűrő
sl	Navodila za uporabo/Tehnični opis Trajni filter
hr	Upute za uporabu/Tehnički opis Trajni filter
ro	Manual de utilizare/Descriere tehnică Filtru permanent
bg	Упътване за употреба/Техническо описание Постоянен филтър
tr	Kullanım Kılavuzu/Teknik açıklama Sürekli filtre
el	Οδηγίες χρήσης/Τεχνική περιγραφή Φίλτρο διαρκείας
zh	使用说明 / 技术说明 长效过滤器

B | BRAUN

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.bbraun.com

AESCULAP® – a B. Braun brand

TA013138 2021-01 Change No. 63662

1. About this document

1.1 Scope

These instructions for use apply to the following reusable filters: JK089 JK090 JK091 and JK092.

- ▶ For article specific instructions for use and material compatibility and lifetime information, see B. Braun eFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

The permanent filters are used as germ barriers in Aesculap sterile containers. The hydrophobic PTFE filters can be used for steam sterilization in a sterilization procedure as per DIN EN 285, in a fractionated vacuum procedure, validated according to DIN EN ISO 17665. The filters can remain in the lid during the mechanical cleaning process and are suitable for up to 1 000 sterilization and wash cycles.

2.1.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.1.3 Contraindications

No known contraindications.

2.2 Safety information

2.2.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

The medical professional will make decisions on concrete applicability based on the warranted properties and technical data.

2.2.2 Sterility

The product is delivered in an unsterile condition.

- ▶ Clean the new product after removing its transport packaging and prior to its initial sterilization.

2.3 Application

⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Prior to each use, inspect the reusable filter for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Always carry out a function test prior to each use of the product.
- ▶ Secure the permanent filter in the sterile container using the appropriate filter holder.
- ▶ For round reusable filters (JK090) check for the correct filter holder (TE674/JK100).
- ▶ Ensure that the filter is not kinked when it is inserted.
- ▶ On first use of the filter, write the start and end of the usage period (equivalent to 1 000 cycles) with a suitable marker pen on the label tab.

Fig. 1 Write the start of use

Fig. 2 Write the end of use

- ▶ At the end of the usage period, change the filter.

3. Validated reprocessing procedure

3.1 General safety information

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeld-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessor is responsible for this.

Note

For up-to-date information about reprocessing and material compatibility, see also B. Braun eFU at eifu.bbraun.com. The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

3.2 General notes

Encrusted or fixated surgical residues can make cleaning difficult or ineffective. Therefore, no more than 6 hours should pass between use and preparation, pre-clean temperatures >45 °C liable to fusing should not be employed and disinfectants liable to fusing (aldehyde or alcohol-based) should not be used.

The cleaning agents must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for reprocessing the product. All of the chemical manufacturer's application specifications must be strictly observed.

- ▶ Do not clean with metal brushes or other abrasives that would damage the product surface.

3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

3.4 Product-specific safety information on the reprocessing method

Note

Aesculap recommends mechanical cleaning in an assembled condition.

When cleaning in a disassembled condition:

- ▶ Ensure that when removing the round filter JK090, the perforation in the center of the filter is not damaged.
- ▶ Ensure that the filter is not damaged during cleaning, such as by instruments that are cleaned along with it.
- ▶ Wash the filter in a disassembled condition separately.

3.5 Manual cleaning/disinfection

⚠ CAUTION

Risk of damage from manual cleaning!

- ▶ Clean the product as far as possible mechanically and in an assembled condition.
- ▶ Use only the cleaning agents for instruments and sterile containers used in the hospital.
- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process or, if necessary, machine the reusable filter.

3.5.1 Manual cleaning and wipe disinfection

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemistry
I	Cleaning	RT (cold)	-	-	D-W	-
II	Drying	RT	-	-	-	-
III	Wipe disinfection	-	>1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol
IV	Final rinse	RT (cold)	0.5	-	FD-W	-
V	Drying	RT	-	-	-	-

DW: Drinking water

FD-W: Fully desalinated water (demineralized, microbiological, at least of drinking water quality)

RT: Room temperature

Phase I

- ▶ Clean the product under running tap water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.

Phase II

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air).

Phase III

- ▶ Wipe all surfaces of the product with a single-use disinfecting wipe.

Phase IV

- ▶ Rinse disinfected surfaces under running demineralized water after the specified contact time has elapsed (at least 1 min).
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air).

3.6 Mechanical cleaning/disinfection

3.6.1 Mechanical neutral or mild-alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	D [°C/°F]	t [min]	Water quality	Chemistry
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Neutral: <ul style="list-style-type: none">■ B. Braun Helimatic Cleaner neutral<ul style="list-style-type: none">- pH neutral- Working solution 0.5 ‰ Mildly alkaline: <ul style="list-style-type: none">■ Concentrate:<ul style="list-style-type: none">- pH = 9.5- <5 ‰ anionic surfactant- Working solution 0.5 ‰ Alkaline: <p>Reprocessing is possible up to pH 10.5, provided the cleaning agent is manufacturer-approved for cleaning aluminum or plastic sterile containers.</p>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	Other process parameters may be feasible with agreement by the hospital hygienist.
V	Drying	-	-	-	According to mechanical cleaning program

DW: Drinking water
FD-W: Fully desalinated (demineralized) water

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.
- ▶ Repeat the cleaning/disinfection process if necessary.

Note

Temperatures of up to 120 °C are permitted for machine drying with hot air.

3.7 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

3.7.1 Visual inspection

- ▶ Remove the reusable filter.
- ▶ Ensure that all contamination has been removed.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. cracked, worn and severely scratched components.
- ▶ Check the product for missing or faded labels.
- ▶ Set aside the product if it is damaged.

3.7.2 Functional test

- ▶ Check for compatibility with associated products.
- ▶ Set aside the product if it is damaged.

3.8 Steam sterilization

- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

4. Technical service

⚠ CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ Do not modify the product.
- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1601
Fax: +49 7461 16-2887
E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

5. Disposal

⚠ WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

Note

The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

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