Enhance[®] & Enhance[®] mini Finishing Systems

Sistemas de acabamento

Putssystem

Fineersystemen

Pudse- og poleringssystem

Pusse- og poleringssystem

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Apdailos sistemos

Enhance[®] PoGo[®] Polishing System

Sistema de polimento

Poleringssystem

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Poleringssystem

Poleringssystem

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Poliravimo sistema

Instruções de utilização _____	Português	3
Bruksanvisning _____	Svenska	5
Gebruiksaanwijzing _____	Nederlands	7
Brugsanvisning _____	Dansk	9
Bruksanvisning _____	Norsk	11
Käyttöopas _____	Suomi	13
Lietošanas instrukcija _____	Latviešu	15
Naudojimo instrukcijos _____	Lietuvių k.	17

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For dental use only. USA: Rx only.

1. PRODUCT DESCRIPTION

The Enhance & Enhance mini Finishers are pre-mounted, single use, aluminum oxide impregnated, cured urethane dimethacrylate resin finishers designed for preparing composite surfaces for their final polish. Enhance PoGo Polishers are pre-mounted, single-use diamond impregnated cured urethane dimethacrylate resin polishing devices designed for use in the final polishing of all accessible composite resin and compomer restorations. When used with a conventional (low-speed) contra-angle latch handpiece, Enhance & Enhance mini Finishers will yield a smooth, contoured surface with pre-polished surface luster on composites and compomers. For optimal surface polish, follow the use of Enhance & Enhance mini Finishers with Enhance PoGo Polishers as outlined below, or Prisma™ Gloss™ Composite Polishing Pastes.

1.1 Delivery forms (Some delivery forms may not be available in all countries)

- Enhance Finishers are available in: Discs, Cups, Points
- Enhance mini Finishers are available in: mini Discs, mini Cups, mini Points
- Enhance PoGo Polishers are available in: Discs, Cups, Points

1.2 Composition

Enhance & Enhance mini Finishers: Polymerized Urethane Dimethacrylate Resin; Aluminum Oxide; Silicon Dioxide
Enhance PoGo Polishers: Polymerized Urethane Dimethacrylate Resin; Fine Diamond Powder; Silicon Dioxide; Plastic Latch-type mandrel

1.3 Indications

- The Enhance & Enhance mini Finishers are indicated for contouring and finalizing the form of composite and compomer restorations.
- The Enhance PoGo Polishers are indicated as the final step in polishing all properly contoured and finished composite resin and compomer restorations.

- Contraindications**
 - Do not reuse. Enhance & Enhance mini Finishers and Enhance PoGo Polishers components are single use, disposable items. Attempts to sterilize these items will adversely affect their performance.
 - Enhance & Enhance mini Finishers and Enhance PoGo Polishers **must not be used in High Speed (>20,000 RPM) handpieces.**

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

	Safety alert symbol This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.
---	---

2.1 Warnings

- Enhance & Enhance mini Finishers and Enhance PoGo Polishers are designed for use without water. As with any rotary instrument, heat will build up with prolonged contact. Use at conventional speed (10,000-15,000 RPM) with light intermittent pressure. Excessive speed and/or pressure may result in heat buildup and may endanger the stability of the instruments. Heat buildup may cause pulp inflammation. Application of an air-cooling stream during finishing and polishing assists with residue removal.
- Use of Enhance & Enhance mini Finishers and Enhance PoGo Polishers components with high speeds and/or excessive lateral pressure may dislodge finisher from the mandrel which may be aspirated, ingested or cause soft tissue injury.
- Aspiration: Seek appropriate medical attention.
- Ingestion: Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. If nausea or illness develop, seek medical attention immediately.
- Eyes, skin, oral mucosa: Thoroughly wash the affected area with water and seek appropriate medical attention.
- Care should be taken to have good control on any rotary instrument to protect patients from injury. Contact with soft tissues (skin, gingiva, mucosa) may result in injury to the tissue. If contact occurs, thoroughly wash the affected area with water and seek appropriate medical attention.
- Dental professionals and patients should wear protective eyewear while using the Enhance & Enhance mini Finishers and Enhance PoGo Polishers.
- Enhance & Enhance mini Finishers and Enhance PoGo Polishers should not be used with patients who have a history of severe allergic reaction to any of the components. If skin rash, irritation, sensitization or other allergic reaction occurs, seek medical attention immediately.

2.2 Precautions

- This product is intended to be used as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
- Use with excess lateral pressure may dislodge Finisher or Polisher from mandrel (see Warnings) or cause quick deterioration and decreased effectiveness. This type of excess pressure may cause gouging of the restoration surfaces.
- Use in a well-ventilated area. Avoid inhaling airborne particulates.
- Protect adjacent restorations from accidental contact with rotating instruments.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Insert the instrument as deeply as possible without excessive use of force into the chuck and check to ensure that it is secure.

- The hand piece must be in proper working condition.
- Immediately remove any instruments that are damaged, bent or no longer run concentrically.
- **Interactions:** None Known.

2.3 Adverse reactions

- Product misuse may cause failure resulting in respiratory, digestive, eye, skin, or mucosa injury. (See Contraindications, Warnings, and Precautions)
- Product misuse may cause excessive heat. (See Contraindications and Warnings)
- Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. (See Warnings)

2.4 Storage

Dry storage is recommended. Store the Enhance & Enhance mini Finishers and Enhance PoGo Polishers between 50°F-80°F (10°C-27°C) Avoid direct sunlight. Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS FOR USE

- Place restoration or composite veneer as per standard techniques. Complete gross reduction of excess and general outline form of the restoration using Prisma™ Finishing Burs, or diamond finishing instrument.

- Excessive heat and/or instrument failure**
 - Use only conventional speed (10,000-15,000 RPM) handpieces
 - Use moderate to light, intermittent pressure

- Insert an Enhance or Enhance mini Finisher into a conventional speed contra-angled handpiece and continue finishing.
NOTE: Enhance & Enhance mini Finishers MUST NOT BE USED IN HIGH SPEED HANDPIECES. The aggressiveness of the Enhance & Enhance mini Finishers is controlled by the pressure applied to the surface of the composite. Use moderate to light, intermittent pressure in a buffing motion. Frequently raise the finisher from surface contact.
- Once final contours are achieved, proceed to polishing with Enhance PoGo Polishers as outlined below, or Prisma™ Gloss™ Composite Polishing Pastes per manufacturer's directions for final high surface luster.

TECHNIQUE HINTS

- The aggressiveness of the Enhance & Enhance mini Finishers is controlled by the pressure applied to the surface of the composite. The greater the pressure, the more material is removed; lighter pressure leaves a smooth surface without removing bulk.
- Prior to proceeding to polishing with Enhance PoGo Polishers or Prisma Gloss Composite Polishing Pastes, the restoration should be finished to a final contour. Surface should be smooth and defect free.

For best results in larger restorations, the use of an intermediate finishing instrument such as an Enhance or Enhance mini Finisher in a conventional contra-angled handpiece will aid in removing scratches and establishing the desired smooth, defect free surface.

- Insert an Enhance PoGo Polisher into a conventional speed contra-angled handpiece.

- Excessive heat and/or instrument failure**
 - Use only conventional speed (10,000-15,000 RPM) handpieces
 - Use moderate to light, intermittent pressure

NOTE: Enhance PoGo Polishers MUST NOT BE USED IN HIGH SPEED HANDPIECES.

- Begin polishing by applying light intermittent pressure at moderate speed to the restoration surface. When using the disc, the flat broad surface of the Enhance PoGo Polisher should be used in contact with the restoration to optimize its effectiveness. To increase surface luster, decrease the pressure applied to the restoration surface using a light buffing motion.
Technique Tip: The flexibility of Enhance PoGo Polishers, when used lightly, allows polishing of established line angles without alteration of form.

Shape	Recommended Application
Enhance Disc	Facial/Smooth surfaces
Enhance mini Disc	Buccal and Lingual interproximal surfaces and embrasures
Enhance PoGo Disc	Facial surfaces and established line angles
Enhance & Enhance PoGo Cup	Cervical and proximal line angle areas
Enhance mini Cup	Subgingival Class V and indirect restoration margin finishing; establishing cusp anatomy
Enhance & Enhance PoGo Point	Posterior occlusal and concave lingual surfaces
Enhance mini Point	Anterior lingual fossa, fine posterior occlusal anatomy

4. HYGIENE

- Cross-contamination**
 - Do not reuse single use products. Dispose of in accordance with local regulations.

5. LOT NUMBER AND EXPIRATION DATE

- Do not use after expiration date. ISO standard used: "YYYY/MM"
- The following numbers should be quoted in all correspondences:
 - Lot number
 - Lot number on package
 - Expiration date on package

	Manufactured by Dentsply LLC 38 West Clarke Avenue Milford, DE 19963 USA Tel.: 1-302-422-4511 www.dentsplysirona.com	EC REP Dentsply DeTrey GmbH De-Trey-Strasse 1 78467 Konstanz, Germany Tel.: 49-7531-583-0 www.dentsply.eu	Distributed by DENTSPLY Canada Ltd. 161 Vinyl Court Woodbridge, Ontario Canada, L4L 4A3
---	---	---	---

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