## Glossar / Glossary

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Dreve Dentamid GmbH is not liable for any damages caused by improper application of the material.

**Dreve Dentamid GmbH**  
Max-Planck-Str. 31 • 59423 Unna/Germany  
www.dentamid.dreve.de

- 2.3 Printing Parameters  
 a) Slice thickness: 50 µm  
 b) Optimal Orientation: 0–15 degree tilted orientation  
 c) Support point size: varies based on support type chosen  
 d) Support density: perimeter of denture base and labial region.  
 The labial side faces the build platform.

- 2.4 Environmental Conditions  
 a) MAX UV print room temperature: 35°C ± 3°C  
 b) Humidity: 30–70%

### 3. Rapid Shape printer (D20II; D30II; D90II) and software

3.1 Hardware  
 Please refer to the printer manufacturer's manual for this information.

3.2 Rapid Shape Printer Software (Autodesk netfabb®)  
 Please refer to the printer manufacturer's manual for this information.

- 3.3 Printing Parameters  
 b) Optimal Orientation: 0–15 degree tilted orientation  
 c) Support point size: varies based on support type chosen  
 d) Support density: perimeter of denture base and labial region.  
 The labial side faces the build platform.

- 3.4 Environmental conditions  
 a) Print temperature: Please refer to the printer manufacturer's manual for this information.  
 b) Humidity: 30–70%

### 4. Recommended curing light equipment (Post curing units)

4.1 UV curing Flood Type System

Manufacturer / Model	PCU LED N <sub>2</sub> (Dreve)	FotoDent® flash (Dreve)	Otoflash G171 (NK- Optics)	RS cure (Rapid Shape)
Supply Voltage	100–240 V	115–230 V	100, 117, 230 V	24 V
Lamp Power	N/A	N/A	2 x 100 W	220 W
Light Intensity	80%	Not adjustable	Not adjustable	According to respective parameter set for FotoDent® denture
Wavelength	410 nm	300–1000 nm	280–700 nm	According to respective parameter set for FotoDent® denture
Curing time	12 min in total	12 min + 1 min gas pre flow	2 x 4000 flashes, turnaround after 4000 flashes, UV-B block bowl and nitrogen	According to respective parameter set for FotoDent® denture

### 5. Notification

The specifications and biocompatibility have been validated using printer software (Autodesk Netfabb®, Composer) and aforementioned printer and process parameters indicated.

## Processing

### Preparation

- Optimal orientation of denture base for printing: 0–15 degree tilted orientation. The labial side faces the build platform, for example as shown in the picture below.


figure 1


figure 2

- During the application of FotoDent® denture we recommend wearing personal protective equipment i.e. suitable gloves, safety goggles etc.
- FotoDent® denture should be well shaken for a minimum of 30 seconds, 1 hour before use.
- Carefully pour FotoDent® denture into the pre-determined container of the production unit.
- Please see further steps within printer manual.
- Remove all bubbles with a cleaned object.

Processing temperature	
Rapid Shape DII series	23 ± 3°C
Asiga MAX UV	35 ± 3°C

- Select the buildstyle/INI file for FotoDent® denture in the printer software. Make sure that all software is up to date.

Rapid Shape: FotoDent® denture pink opaque #D35501  
 FotoDent® denture pink transparent #D35500

Asiga: FotoDent® denture pink opaque\_D35501.ini  
 FotoDent® denture pink transparent\_D35500.ini

### Further processing

- After completion of the building process, direct further processing (see following steps) is recommended.
- After the platform has been raised, a dripping time of approx. 10 minutes is recommended.
- Prewash the printed parts in 97 % Isopropanol (IPA) in an ultrasonic bath for 2 minutes in a container, afterwards wash for additional 2 minutes in a container with fresh Isopropanol. Dry the printed part with compressed air. For a particular effective cleaning with isopropanol, the use of explosion-protected ultrasonic units is recommended.

### Post Curing

PCU LED N <sub>2</sub>	12 min, 80% light output (incl. nitrogen)
FotoDent® flash	12 min + 1 min gas pre flow
Otoflash G171	2 x 4000 flashes, turnaround after 4000 flashes, UV-B block bowl and nitrogen
RS cure	According to respective parameter set for FotoDent® denture.

- Afterwards temper the denture base in an oven at 90 °C for 10 minutes.
- Recommended bonding adhesive for denture teeth: VITA VIONIC® BOND or IvoBase CAD Bond.
- Remove the container from the printer and filter the resin through a fine 190 Micron paint strainer, if:
  - print has failed partially or completely, or
  - particulates of polymerized residues are visible in the container or stick to the bottom
- To avoid bubbles, let the material rest for 1 to 2 hours prior to use, after filtration.
- Discard and replace the FotoDent® denture material with a new batch, if contamination, evident gelation or polymerization is observed after filtering.
- Do not mix different batches of FotoDent® denture material.

**To avoid compromising material quality, do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may impair biocompatibility, lead to modified mechanical characteristics and/or color variations of the FotoDent® denture material.**

## Caution

The specifications and biocompatibility have been validated using printer software (Autodesk Netfabb®, Composer) and aforementioned printer and indicated process parameters, herein. Any unauthorized changes to the process equipment, parameters or software may result in a device that is out of specifications.

Please contact Dreve Dentamid GmbH for a list of validated software and process hardware.

## Disposal

Dispose FotoDent® denture in accordance with local/regional/national/international regulations.

See also information in safety data sheet.

## Hazards identification

Before use observe the Material Safety Data Sheet of FotoDent® denture.

**Hazard statements: H317** May cause an allergic skin reaction. **H411** Toxic to aquatic life with long lasting effects.

**Precautionary statements: P261** Avoid breathing dust/fume/gas/mist/vapors/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/eye protection/face protection. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P321** Specific treatment (see on this label). **P501** Dispose of contents / container in accordance with local/regional/national/international regulations.

**Contains:** Neopentylglycol(PO)2 Diacrylate, Aliphatic Urethantrimethylacrylate, pigment paste

## Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® denture packaging. In case of complaints please always indicate the lot number of the product. Do not use the product after the expiry date.

## Cleaning

For cleaning dentures made of FotoDent® denture the use of MaxiFresh™ cleaning tablets is recommended. Please be aware of the MaxiFresh™ instruction for use. Long term contact of dentures with coloring beverages (e.g. tea or coffee) may discolor the denture.

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## Intended use

EN

**Indication:** FotoDent® denture is a methacrylate-based material for DLP systems with 385nm-LED for production of denture base.

**Indication for use:** FotoDent® denture is a light curing resin intended for manufacturing of full and partial removable dentures.

**Contraindication:** FotoDent® denture is contraindicated for individuals with known allergy to methacrylates and acrylates.

**Patient target group:** Patients with a partially toothed or completely edentulous upper and/or lower jaw.

**Intended users:** The dental technicians or specialist staff in the dental laboratory is responsible for the production of denture bases in 3D printing. The dentist uses the finished denture base on the patient.

## Requirements

- Printer:
  - Rapid Shape D20II, D30II, D90II, D20+, D30+ (385nm)
  - Asiga MAX UV (385nm)
- Software: Autodesk Netfabb®, Composer
- Post curing unit: PCU LED N<sub>2</sub>, FotoDent® flash, Otoflash G171, RS cure

## Basic Material Properties

Characteristics/Properties	Units	Specifications
<b>Before curing (liquid state)</b>		
Viscosity	cps	>1000 and <2000
Density	g/cm <sup>3</sup>	>1.05 and <1.2
<b>After curing</b>		
Density	g/cm <sup>3</sup>	>1.15 and <1.25
Flexural strength	MPa	>80
Flexural Modulus	MPa	>2000

## Specific Manufacturing Considerations

### 1. Digital denture base file

- 1.1 File format: STL file
- 1.2 Digital design: denture base or baseplate

### 2. Asiga printer (MAX UV) and software

- 2.1 Hardware
  - Please refer to the printer manufacturer’s manual for this information.
- 2.2 Asiga Printer Software (Composer)
  - Please refer to the printer manufacturer’s manual for this information.