- For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website www.dentamid.dreve.de.

Post curing FotoDent® model2 beige-opaque	
PCU LED N ₂	
Rapid Shape DII / Asiga MAX UV	8 minutes, 80 % light power (vacuum)
FotoDent® flash	
Rapid Shape DII / Asiga MAX UV	8 minutes + 1 minute gas pre-flow / 12 minutes + 1 minute gas pre-flow
PCU 90	8 minutes (incl. 4 minutes gas pre-flow)
RS Cure	According to respective parameter set for FotoDent® model2
Post curing FotoDent® model2 black	
PCU LED N ₂	
Rapid Shape DII / Asiga MAX UV	12 minutes, 80 % light power (vacuum)
RS Cure	According to respective parameter set for FotoDent® model2
Post curing FotoDent® model2 white	
PCU LED N ₂	
Rapid Shape DII / Asiga MAX UV	6 minutes, 40 % light power (vacuum)
RS Cure	According to respective parameter set for FotoDent® model2

- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® model2. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with FotoDent® model2.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the FotoDent® model2 material.

Possible dangers

Please note the safety instructions for processing FotoDent® model2.

Hazard statements: **H315** Causes skin irritation. **H319** Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H413** May cause long lasting harmful effects to aquatic life.

Precautionary statements: **P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: 2-hydroxyethyl methacrylate, Diurethandimethacrylat, hydroxypropylmethacrylate, monoester with propane-1,2-diol, pigment paste

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® model2 packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Intended use

FotoDent® model2 is a material based on (meth-)acrylate resin for DLP systems with 385 nm LED for the production of dental models. For the post curing of FotoDent® model2 built parts, the light curing units PCU LED N₂, FotoDent® flash and RS cure are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® model2 we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® model2 should be well shaken approx. 1 hour before use.
- Carefully pour FotoDent® model2 into the pre-determined container of the production unit.
- Please see further steps within printer manual.
- Remove all bubbles with a cleaned object.

Processing temperature	
Rapid Shape DII printer	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C

- Select the product specific parameter settings for FotoDent® model2 for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2x6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended. Alternatively, cleaning with isopropanol can be carried out in a separate beaker under pouring conditions for approx. 2x6 minutes. Recommendation: Beaker and magnetic stirrer, laboratory shaker or ultrasonic cleaner. In case of insufficient cleaning (such as shiny areas on the object after drying) repeat cleaning process.

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