FREEPRINT®

ortho

3D PRINTING MATERIAL

MED RESIN

385 nm

drilling templates

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Resin for dental 3D printing

INDICATION

Orthodontic bases, surgical guides and X-ray templates for implant dentistry, occlusal orthos, fixation and transfer keys

PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

INTENDED USERS

Dentist, dental technician

SUITABLE FOR THE FOLLOWING DLP-PRINTER

Asiga 385 nm 405 nm		MiiCraft 385 nm	Rapidshape 385 nm 405 nm		W2P 385 nm
MAX UV Pro2 UV Pico2 UV	MAX 405 Pro2 405 Pico2 405	MiiCraft y-Serie	DII-Serie	D-Serie	Solflex-Serie

Printers may only be operated using material parameters authorised by DETAX!

PROCESSING

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 10).
- ▶ After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- ▶ Minimum material thickness for the design is 1.5 mm
- ▶ Maximum curing depth* at direct post-exposure: clear: 6 mm
* In case of large objects and exposure on both sides, the material thickness can be up to 12 mm (Example FREEPRINT® ortho –with a curing depth of 6 mm).
- ▶ Polish surface mechanically.
- ▶ Processing temperature 23 °C ± 2 °C.

FREEPRINT® ORTHO 385 - STERILIZATION OF DRILLING TEMPLATES

- ▶ Locally applicable regulations and hygiene guidelines are to be observed.
- ▶ Important: after main cleaning and prior to post-curing, the components to be sterilized are to be dried, e.g. to be heated in an oven to approx. 40 °C for 30 minutes to remove remaining solvent from the cleaning process
- ▶ FREEPRINT® ortho 385 is validated for the „W&H Lisa 522, Program UNIVERSAL 121“ sterilization process (technical data on the right). Steam sterilization may only be performed with equipment that complies with the standards EN 13060 and EN 285.

- ▶ Responsibility for sterility lies with the user. Safety Information

FREEPRINT® ORTHO (405 NM)

- ▶ Do not use heat-based methods for disinfection or sterilization. This could possibly deform the workpiece.

TECHNICAL DATA FOR THE STERILISATION PROCESS

- ▶ A fractionated pre-vacuum with a plateau phase at 121 °C/ 15 minutes is to be selected, whereby this S to the requirements for a B-cycle small steriliser in terms of its venting performance.

SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

CLEANING

By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

By the patient

For daily cleaning, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards. Cleaning tablets can be used following the manufacturer's instructions.

Note

- ▶ Do not use toothpaste; the abrasive particles are coarse and scratch the surface.
- ▶ Do not use any acids or solvents.

NOTE

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.
- ▶ Please follow the instructions on the safety data sheet!

For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.de and to the competent authority of the Member State in which the user and/or patient is established.

STORAGE

FREEPRINT® ortho is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

CONTRAINDICATION

Contains (meth)acrylics and phosphine oxide.

Some ingredients of **FREEPRINT® ortho** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **FREEPRINT® ortho** only insert intraorally in completely polymerised state.

ADVERSE EFFECTS

Product may cause allergic reactions.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Symbols:

MD Medical Device

Storage:

Application:

At 23 °C ± 2 °C

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath.

Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath. Prior to post-exposure, check the openings and additional bore holes for residues. Then blow off with compressed air.

Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

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GmbH & Co. KG