# ROEKO Gelatamp

Instructions for use	EN
----------------------	----

**Properties:**

Gelatamp (14 × 7 × 7 mm) is a haemostatically acting gelatine sponge to which 5% colloidal silver is added. It facilitates optimum wound treatment when applied to a surgical cavity and can be cut to the required size to fit smaller wound cavities. The evenly porous foam structure absorbs its own weight in blood several times over, promotes thrombocyte aggregation due to the large surface and fills the wound cavity. The plug thus formed has a constant volume, fits snugly and stabilises blood coagulum. This prevents the formation of fissures and secondary cavities which, without Gelatamp, could form by contraction of the blood coagulum and trigger infection due to the invasion of contaminated saliva. Callus formation is not hindered in this way. Impaired wound healing in larger surgical cavities is thus avoided. Gelatamp remains in the wound and is completely absorbed within four weeks. The addition of colloidal silver has an antimicrobial effect. Unlike other potential antimicrobial additives, colloidal silver cannot be washed away from the sponge so that its insolubility produces a long lasting depot effect. Gelatamp is sold gamma-sterilised.

**Composition:**

10 mg Gelatamp gelatine sponge contain:
Hardened gelatine Ph. Eur. 9,5 mg
Colloid silver Ph. Eur. 0,5 mg

**Indication:**

The treatment of wound cavities after tooth extractions.

**Instructions for use:**

After opening the aluminium seal, Gelatamp is removed under aseptic conditions and is immediately ready for use. The size of the small sponge can be adjusted to fit the wound cavity if need be, without squeezing. Two Gelatamp plugs can be used for larger wounds.

The sponge is immediately placed dry and unsqueezed in the fresh, blood-filled wound cavity and should be allowed to absorb as much blood as possible (visual inspection). The sponge must not protrude over the inner gingival margin so as not to adversely affect the desired epithelisation. This constitutes definitive treatment of the alveolus and there is no need for subsequent rinsing.

**Sterility:**

Gelatamp is supplied gamma-sterilised. The sterility of undamaged, unopened packs is guaranteed until the first application. If packages are damaged prior to use they have to be discarded.

**Contra-indications:**

Infected, secretory and purulent wounds; hypersensitivity to silver or gelatine.

**Side effects:**

Recovery may be delayed if used incorrectly, especially in the case of contaminated alveolar spaces.

**Interactions with other drugs:**

Concomitant application of alveolar pastes, wound plugs or other inserts and moistening with liquid substances or medication such as anion-active surfactants, antiseptics containing phenol, ethacrine lactate, dilute mineral acids and concentrated saline solutions can hinder wound recovery and should therefore be avoided. Please avoid contact with hydrogen peroxide as this will bleach the sponge.

**Storage instructions:**

Gelatamp possesses moisture-absorbing properties. The package must therefore be closed immediately after use and stored dry, protected from light. Keep medical devices out of the reach of children. Gelatamp has a shelf-life of 3 years. It can be used unopened up to the expiry date. Use within four weeks after opening. The product should no longer be used once the date on the packaging and label has elapsed.

Unused sponges can easily be disposed of with laboratory waste whereas, the plastic container, lids and aluminium foil, etc. can be recycled.

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