**GRADIA® DIRECT LoFlo**LIGHT-CURED FLOWABLE COMPOSITE  
light-cured composite, micro-filled hybrid type  
composite resin

Thanks to its HGF (High Density Radiopaque Pre-Polymerized Filler) technology, it contributes to achieve excellent abrasion resistance. The thicker consistency (lower flowability) is excellent for cavities of class III, IV and root surface caries.

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

- Restoration of Class I, II, III, IV cavities (particularly for small Class I cavities / shallow Class V cavities / other similar cavities).
- Restoration of root surface caries.
- Restorations in deciduous teeth.
- Filling tunnel shaped cavities.
- Sealing hypersensitive areas.
- Liner / base filling in cavity undercuts.
- Sealer.
- Splicing mobile teeth.
- Addition to composite restorations. (Fig. 1)

**CONTRAINDICATIONS**

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- CLEANING AND DISINFECTION:**  
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged or discolored.
- DO NOT IMMERSE:** Thoroughly clean device to prevent drying and accumulation of contaminants.
- Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

**CAUTION**

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Take special care that the patient should not swallow material.
- Do not mix with other similar products.
- Dosieringur nařízeník nemožné do báry sterilizaci v autoklavu ani chemikálií.
- Personal protective equipment (PPE) such as gloves, face masks and eye wear should always be worn.
- Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.europeo.com>. May also be obtained from your supplier.

They can also be obtained from your supplier.

**GRADIA® DIRECT LoFlo**  
LIGHT-CURED FLOWABLE COMPOSITE

GRADIA® DIRECT LoFlo is a multi-purpose light-cured, radiopaque, micro-filled hybrid type composite resin. Thanks to its new HDR (High Density Radiopaque Pre-Polymerized Filler) technology, it contributes to achieve excellent abrasion resistance. Also the thicker consistency (lower flowability) is excellent for cavities of class III, V and root surface caries.

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

- Restoration of Class I, II, III, IV cavities (particularly for small Class I cavities / shallow Class V cavities / other small cavities).
- Restoration of root surface caries.
- Restorations in deciduous teeth.
- Filling tunnel shaped cavities.
- Sealing hypersensitive areas.
- Liner / base / filling in cavity undercuts.
- Sealant.
- Splinting mobile teeth.
- Additions to composite restorations. (Fig. 1)

**CONTRAINDICATIONS**

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**DIRECTIONS FOR USE**

- Preparations
  - Hold the syringe upright and remove the wing cap by turning counter-clockwise. Take care not to expose material to direct light from the dental lamp or natural light (Fig. 2).
  - Promptly and securely attach the dispensing tip (plastic or needle type) to the syringe by turning clockwise (Fig. 2)(Fig. 3).
- Note : Take care not to attach the dispensing tip too tightly. This may damage its screw.
- After attaching the dispensing tip, protect it with the cover until ready for use in order to avoid exposure to light (Fig. 2).

## 2. Cavity Preparation

Prepare cavity using standard techniques. Note : For pulp capping, use calcium hydroxide.

3. Shade Selection  
Select shade from 7 shades of A1, A2, A3, A3.5, A3, CV and BW.

A shades are based on Vita® shade.

\*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

4. Use of Light-Curing Bond System

For bonding GRADIA® DIRECT LoFlo to tooth structure, use a light-curing bonding system such as G-ænial® Bond, G-Premio BOND or G-BOND.

Note : When using light-curing bonding system, follow manufacturer's instructions for use.

5. Placement of GRADIA® DIRECT LoFlo

- Remove cover from the dispensing tip on the syringe (Fig. 5).
- Prior to extrusion of material, gently trial push the syringe tip to decide the moment to make sure that the dispensing tip is securely attached to the syringe. To remove any air from the dispensing tip, with the tip pointing upwards gently push forward the syringe plunger until material reaches the mouth of the tip (Fig. 6).

Note : If there is air inside the dispensing tip, air bubbles may be formed at the time of injection.

3) Place the dispensing tip as close as possible to the cavity, and slowly push the plunger to inject material into it (Fig. 7). Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.

Note :

- When attaching the dispensing tip, make sure that no material is sticking to the joint between the tip and the syringe in order to ensure a tight connection.
- If the dispensing tip does not glide smoothly, remove the dispensing tip and extrude material directly from the syringe to make sure that material is coming out.

c. The material will start to harden if exposed to the dental light or ambient light. Be sure to protect it from light when working from a mixing pad.

d. After use, immediately remove the dispensing tip and tightly close the syringe with the wing cap.

**Clinical Hint 1**

In order to inject effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

**Clinical Hint 2**

Where higher abrasion resistance and thicker consistency are required, it is recommended to place GRADIA® DIRECT LoFlo instead of GRADIA® DIRECT LoFlo or low viscosity flowable composites.

**Clinical Hint 3**

When filling a large cavity, it is recommended to place material incrementally into the cavity. Another effective method is to use GRADIA® DIRECT LoFlo for filling in undercuts or as a liner/base, and then to place composite resin (GRADIA® DIRECT, G-ænial®, Essentia®, etc.) on top.

## 6. Light Curing

When curing the GRADIA® DIRECT LoFlo using a light curing unit (Fig. 8). Refer to the following chart for Irradiation Time and Effective Depth of Cure.

Note : When light curing material, wear protective glasses.

**Irradiation Time and Effective Depth of Cure**

Irradiation time	Shade	A1, A2, A3	A3.5, BW	A03, CV
20 sec. Halogen/LED (700mW/cm²)		2.0 mm	1.5 mm	1.5 mm
10 sec. High Power LED (more than 1200mW/cm², e.g. GC-D-Light Pro)		3.0 mm	2.5 mm	3.0 mm
40 sec. Halogen/LED (700mW/cm²)		3.0 mm	2.5 mm	3.0 mm

## 7. Shaping and Polishing

Shape and polish using standard techniques.

**SHADES**

A1, A2, A3, A3.5, A03, CV (Cervical color), BW (Bleaching White)

**STORAGE**

Recommended for optimal performance, store in a cool and dark place (4 - 25°C / 39.2 - 77°F) (Shelf life : 2 years from date of manufacture)

**PACKAGING**

- Syringe packages: syringe 1.3g (0.8ml) x 2 (in one shade), 4 dispensing tips (needle type), 2 light protective covers.
- Dispensing tip package:
  - 20 dispensing tips (needle type), 2 light protective covers
  - 20 dispensing tips (plastic type), 2 light protective covers

**CAUTION**

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.

3. Take special care that the patient should not swallow material.

4. Do not mix with other similar products.

5. The dispensing tip cannot be sterilized in an autoclave or chemiclaire.

6. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be more hazardous according to GHS. Always familiarize yourself with the safety data sheets available at <http://www.geceurope.com>

They can also be obtained from your supplier.

**CLEANING AND DISINFECTION:**

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

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