HeraCeram Zirkonia

| | | |
|------|----------------------------------|----|
| (DE) | Gebrauchsanweisung | 2 |
| (GB) | Instructions for use | 4 |
| (FR) | Mode d'emploi | 6 |
| (ES) | Instrucciones de uso | 8 |
| (IT) | Istruzioni per l'uso | 10 |
| (PT) | Instruções de uso | 12 |
| (NL) | Gebruiksaanwijzing | 14 |
| (SE) | Bruksanvisning | 16 |
| (DK) | Brugervejledning | 18 |
| (NO) | Bruksanvisning | 20 |
| (FI) | Käyttöohjeet | 22 |
| (CZ) | Návod k použití | 24 |
| (HU) | Használati utasítás | 26 |
| (LV) | Lietošanas instrukcija | 28 |
| (LT) | Naudojimo instrukcija | 30 |
| (PL) | Instrukcja obsługi | 32 |
| (HR) | Upute za uporabu | 34 |
| (RO) | Instrucțiuni de utilizare | 36 |
| (RU) | Инструкция по применению | 38 |
| (SK) | Návod na použitie | 40 |
| (BG) | Инструкции за употреба | 42 |
| (UA) | Інструкція по застосуванню | 44 |
| (EE) | Kasutusjuhend | 46 |

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Stand: 2017-07

HeraCeram Zirkonia is the veneering ceramic for frameworks made of zirconium dioxide, CTE ($\alpha_{25-500^{\circ}\text{C}}$) 10.5 $\mu\text{m}/\text{mK}$.

Framework preparation: A very thin and even layer of Zr-Adhesive paste is applied and fired in a vacuum (see firing table) to ensure optimal adhesion between the zirconium dioxide and HeraCeram Zirkonia.

Liner application: If necessary, a liner can also be applied in an evenly thin layer and fired to cover the framework and provide base colouration. This step can be repeated if the liner does not cover it completely. Alternatively, the HeraCeram Zirkonia increasers can also be used for this base layer.

Dentine/incisal layers: The ceramic veneer is designed in a two-layer technique or, depending on the desired aesthetic requirements, with additional individual materials.

Note: Wear a face mask and eye protection when grinding ceramic materials and work with a suction device. Avoid breathing in the ceramic dust.

Glaze firing: The degree of glazing of the ceramic surface can be defined by the temperature and the holding time at the final temperature during glaze firing. For this reason, the specifications for the glaze firing temperature and holding time are for orientation only and may need to be adjusted to achieve the respective required result. It is also possible to polish HeraCeram Zirkonia mechanically. Our Signum® HP pastes have proven themselves ideal for final polishing.

Colour mapping table

| | A1 | A2 | A3 | A3,5 | A4 | B1 | B2 | B3 | B4 | C1 | C2 | C3 | C4 | D2 | D3 | D4 | BL 1 | BL2 | BL3 | BL4 |
|---------------------------|-----------------------|-------|-------|---------|-------|-----------------------|-------|-------|-------|-----------------------|-------|-------|-------|----------|----------|-----------------------|----------|----------|----------|----------|
| Liner | LA1 | LA2 | LA3 | LA3,5 | LA4 | LB1 | LB2 | LB3 | LB4 | LC1 | LC2 | LC3 | LC4 | LD2 | LD3 | LD4 | L BL 1-2 | L BL 1-2 | L BL 3-4 | L BL 3-4 |
| Increaser | IN A1 | IN A2 | IN A3 | IN A3,5 | IN A4 | IN B1 | IN B2 | IN B3 | IN B4 | IN C1 | IN C2 | IN C3 | IN C4 | IN D2 | IN D3 | IN D4 | | | | |
| Increaser | | | In C | In C | | | | In S | In S | | IN T | IN T | IN T | IN T | IN T | IN T | | | | |
| Dentine | D A1 | D A2 | D A3 | D A3,5 | D A4 | D B1 | D B2 | D B3 | D B4 | D C1 | D C2 | D C3 | D C4 | D D2 | D D3 | D D4 | D BL1 | D BL2 | D BL3 | D BL4 |
| Incisal | S1 | S1 | S2 | S2 | S4 | S1 | S1 | S2 | S4 | S 1 | S3 | S3 | S3 | S1 | S2 | S2 | S BL | S BL | S BL | S BL |
| Mamelon secondary dentine | MD 1 | MD 1 | SD 2 | SD 3 | SD 4 | MD 2 | MD 2 | MD 3 | MD 3 | MD 2 | SD 1 | SD 2 | SD 2 | MD 1 | MD 3 | MD 1 | | | | |
| Value | VL 1 | VL 2 | VL 3 | VL 4 | VL 4 | VL 1 | VL 2 | VL 3 | VL 4 | VL 1 | VL 2 | VL 3 | VL 4 | VL 2 | VL 3 | VL 4 | | | | |
| Opal incisal | OS1 | OS1 | OS2 | OS2 | OS4 | OS1 | OS1 | OS2 | OS4 | OS 1 | OS3 | OS3 | OS3 | OS1 | OS2 | OS2 | | | | |
| Shoulder ceramics HM/LM | 1 | 1 | 2 | 2 | 6 | 3 | 3 | 4 | 4 | 5 | 5 | 6 | 6 | 1 | 2 | 4 | | | | |
| Mono Liner | | ML-A | ML-A | ML-A | ML-A | | ML-B | ML-B | ML-B | | ML-C | ML-C | ML-C | | ML-A | ML-B | | | | |
| Mono Body | MB Light | MB-A | MB-A | MB-A | MB-A | MB Light | MB-B | MB-B | MB-B | MB Light | MB-C | MB-C | MB-C | MB Light | MB Light | MB-C | | | | |
| Body Stains/Glaze | 50% BS A 50% Glaze | Glaze | BS A | BS A | BS A | 50% BS B 50% Glaze | Glaze | BS B | BS B | 50% BS C 50% Glaze | Glaze | BS C | BS C | BS A | BS C | 50% BS C 50% Olive | | | | |

General firing table

| | | Zr-Adhesive | Liner | Shoulder ceramic HM 1 | Shoulder ceramic HM 2 | Primary dentine | Secondary dentine | Glaze | Mono Dentin and glaze firing | Universal stains | Correction material | Shoulder ceramic LM |
|----------------------------|--------|-----------------|-------|-----------------------|-----------------------|-----------------|-------------------|-------|------------------------------|------------------|---------------------|---------------------|
| Start temperature | °C | 600 | 600 | 600 | 600 | 600 | 600 | 600 | 600 | 600 | 600 | 600 |
| Pre-drying and pre-heating | min | 6 | 6 | 4 | 3 | 5 | 5 | 4 | 5 | 6 | 4 | 4 |
| Temperature increase | °C/min | 100 | 100 | 100 | 100 | 100 | 100 | 100 | 100 | 100 | 100 | 100 |
| Final temperature | °C | 1050 | 880 | 870 | 860 | 860 | 850 | 850 | 850 | 850 | 810 | 790 |
| Holding Time | min | 10 ¹ | 1 | 1 | 1 | 1 | 1 | 1 | 0,5-1 | 0,5-1 | 1 | 1 |
| Vacuum start | °C | 600 | 600 | 600 | 600 | 600 | 600 | – | 600 | 600 ² | 600 | 600 |
| Vacuum stop | °C | 950 | 880 | 870 | 860 | 860 | 850 | – | 850 | 850 ² | 810 | 790 |

¹ = under vacuum

² = vacuum only for Mono

Safety instructions

This product or one of its components may, in particular cases, cause hypersensitive reactions. If suspected, information about the ingredients can be obtained from the manufacturer. To be used only by dental technicians for its intended purpose. The data given are furnished to the best of our knowledge and are based on in-house test results. Never swallow or ingest this product! Avoid contact with eyes. In case of contact with eyes, rinse immediately and thoroughly with plenty of water and contact an ophthalmologist. The use of this product is contraindicated in cases of known or suspected allergies to components of the ceramic. In the case of a suspected allergy, we recommend taking an allergy test prior to use.

Please observe the information in the corresponding safety data sheet.

Please quote the batch number in all correspondence about the product.

Batch number Lot: See note on packaging.

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Kulzer GmbH
Kulzer, LLC