HeraCeramPress

(DE) Gebrauchsanweisung	2
(GB) Instructions for use	5
(ES) Instrucciones de uso	8
(IT) Istruzioni per l'uso	11
(PT) Instruções de uso	14
(FI) Käyttöohjeet	17
(PL) Instrukcja obsługi	20

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Purpose

Manufacture of single and multi-surface pressable ceramic inlays, onlays and veneers as well as full and veneered ceramic crowns using the stain or layering technique.

- **Stain technique:** Fully anatomical full and partial crowns, inlays, onlays and veneers. The final effect is achieved by staining with HeraCeram Stains and Glaze universal.
- **Layering technique:** semi-anatomical pressable ceramic inlays, onlays and veneers and anterior and premolar crowns veneered with HeraCeram.

Contraindications

- Bridge restorations
- Very deep subgingival preparations
- Patients with parafunctions (e.g. bruxism)
- Severely reduced residual dentition

Preparation

1. Model preparation

- Definition of the preparation margin
- For the adhesive gap, apply the spacer
 - for crowns and veneers in two layers up to approx. 1 mm from the preparation margin
 - for inlays and onlays of up to three layers, right up to the preparation margin of the die
- Insulate die from wax
- In the case of a known or suspected allergy to the ceramic ingredients, the use of this product is contraindicated.

2. Modeling

- For modeling only use residue-free, full burn-out modeling materials and separating liquids intended for the pressing technique.
- Objects for the layering technique must be 0.8 mm thick minimum.
- Premolars: reduced anatomical shape (occlusal profile) for the layering technique.
- Framework thickness should be at least 50% of the total thickness.
- For molar crowns the stain technique is to be used.

3. Spruing

The objects must be sprued with approx. 5–6 mm long wax wires. For smaller inlays and crowns 2.5–3.0 mm thick in \emptyset , for more voluminous restorations 3.5 mm thick in \emptyset .

- Apply sprue at the point with the greatest volume, without tapering
Pressable object and sprue ideally form a line = direction of flow.

Before waxing on the channel former, determine the wax weight of the objects for the quantity of pellets to use:

- wax weight of up to 0.6 g → 1 pellet
- wax weight of up to 1.4 g → 2 pellets
- Position pressable objects all at the same level, therefore the object margins are all on one level – as a result, the sprues may vary in length.
- Position sprues on the base with a slight outward inclination (angle of approx. 60–70°C).
For inlays and onlays, the basal surface should face outwards.

4. Investment

- Please follow the manufacturer's instructions when processing the investment materials (e.g. Heravest Press).
- Lightly vibrate the investment material while pouring it bubble-free into the mould former and put on the ring gage. The ring gage ensures the correct height of the mould and the parallel orientation of the mould to the pressing direction of the pressing furnace.
- Allow the investment material to set free of vibration.

5. Pre-heating

- The ring and mould former should be removed after the setting period. For safe positioning of the mould, clean the top surface of the mould.
- Put mould and aluminium oxide press plunger in the preheating furnace.
HeraCeramPress pellets and Heravest disposable press plunger must not be preheated.
- Heat the mould (Linear or Speed heating – according to the instructions for use of the investment material) to 850°C.
- Holding times: 200 g mould → 60 min
100 g mould → 60 min

Note: With Speed heating, the holding time starts when the furnace has reached again the final temperature!

6. Pressing

- Activate the press programme.
- Take mould out of the preheating furnace when the press furnace has reached the start temperature of 700°C.
- Place 1–2 HeraCeramPress pellets in the mould channel.
- Place press plunger in the channel (disposable press plunger or aluminium oxide press plunger).
- Put immediately the mould so assembled in the press furnace and start the press programme.
- After pressing, remove the mould from the press furnace and allow to cool to room temperature.

Avoid cooling of the mould between removal from the preheating furnace and the start of pressing procedure

General Press Recommendation

Start temperature	Vacuum start	Vacuum stop	Vacuum	Heat rate	Final temperature	Hold	Pressing	Pressure force
°C	°C	°C	hPa	°C/min	°C	min	min	bar
700	700	1030	30	50	1030	15	8	depending on type of furnace

7. Deflasking

- Determine and mark the position of the pressed objects with the aid of a press plunger. Then make a deep cut into the investment material with a cut-off wheel at the end side of the press plunger. Separate the investment parts from each other through a rotating movement.
- Remove investment material by means of cautious sandblasting.

Note: Avoid localised sandblasting which will cause overheating of the surface.

Do not use aluminium oxide

- Rough investment removal with glass beads (4 bar, 50 µm, do not use a recycling blasting unit) - do not completely expose the objects yet.
- Complete investment removal: continue to blast extensively on the surfaces using reduced pressure (1.5 bar max., 50 µm, glass beads).

Finishing

Pressable ceramic must be processed in (water-)cooled conditions. Apply low pressure so that the pressable ceramic do not overheat and no cracks occur. Use diamond abrasives which are preferably in mint condition. Please pay attention to the minimum thickness of the inlays (2.0 mm minimum in the isthmus) and the veneer shells or crown copings (0.8 mm).

- Avoid localised overheating when separating and grinding the sprues.
- Fitting the objects on the model: fit carefully on the die. In case of precontacts, thinly stain the die with occlusal paste / spray. Remove precontacts carefully using fine diamonds.

Note: Wear a face mask and eye protection when grinding ceramic materials and work with an extraction unit. Avoid inhaling the ceramic dust.

- Clean object with a steam cleaner. Avoid overheating.

8. Veneering

- **Layering technique:** Apply a thin layer of adhesive paste to the surfaces to be veneered and fire. (See below for firing instructions)
- The HeraCeram materials are then layered and fired according to the HeraCeram user instructions.
- **Stain technique:** The restoration must be clean and grease-free. A roughened surface is of advantage.
- Personalize restorations with HeraCeram Stains universal and fire. For more intense chromatic results, the procedure can be repeated.
- Finally, glaze with HeraCeram Glaze universal and fire.
- **Firing procedure:** Veneers and inlays are positioned and fired on the smallest possible piece of firing wool placed of top firing pins.
- Crowns are fired on suitable firing pins, fan pins, special firing wires or pins with adapted platinum foil or firing wool.

Firing recommendations

 **The firing temperatures given are intended as guidelines. Deviations are possible due to varying furnace performance.**

General Firing Programme	Pre-heat or start temperature	Pre-dry and pre-heat time	Heat rate	Final temperature	Holding time	Vacuum start	Vacuum stop
	°C	min	°C/min	°C	min	°C	°C
Adhesive	600	5	100	860	1	600	860
Stain and glaze firing using Stains / Glaze universal	600	6	100	860	0,5–1	–	–

Please follow our detailed processing instructions for processing of HeraCeramPress and the instructions in our safety data sheet.

Safety instructions:

Hypersensitivity to the product or its components cannot be ruled out in individual cases – in a suspected case, please contact the manufacturer with regard to the ingredients. In the case of a suspected allergy, an allergy test is recommended before treatment. The data indicated corresponds to the best of our knowledge and is based on internal investigations. Do not swallow or consume this product. Avoid contact with the eyes. If the product does come into contact with the eyes, rinse them thoroughly with plenty of water immediately and consult an ophthalmologist.

If you wish to contact us with regard to the product, please always indicate the lot designation.

Designation lot LOT: See note on packaging.

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Kulzer GmbH