
HeraCeram NP-Primer

- DE** Gebrauchsanweisung
- NL** Gebruiksaanwijzing
- GB** Instructions for use
- SE** Bruksanvisning
- FR** Mode d'emploi
- FI** Käyttöohjeet
- ES** Instrucciones de uso
- LT** Naudojimo instrukcija
- IT** Istruzioni per l'uso
- HU** Használati utasítás
- PT** Instruções de uso
- RU** Инструкция по применению

99000712/11

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GB Instructions for use

Purpose

NP Primer serves to condition non-precious metal alloy surfaces for veneering with HeraCeram.

Its feature is the prevention of uncontrolled, excessive oxidation that may affect the bond between the alloy and the ceramic. The oxide formation in non-precious metal alloys is controlled such that the surface of the framework is perfectly covered by the opaquer.

NP Primer is available as a ready-to-use paste.

NP Primer instructions for use

Trim the non-precious metal surface to be veneered in the usual manner, then blast it with aluminium oxide (100-150µm, 3-4 bar), and clean it. Oxide firing is not required (Fig. 1).

Use the paste opaque brush to apply NP Primer to the surface to be veneered in a very thin layer (high brush pressure!) and fire under vacuum at 950°C for 1 minute (Fig. 2).

Notice:
NP-Primer paste is thixotrope. It has to be stirred before use!

After firing, the coated surface has a silky lustre (Fig. 3).

The surface is now prepared for ceramic veneering. Perform the next steps in accordance with the instructions for use for HeraCeram (Fig. 4-6).

	START	↑	→	V _{start}	V _{stop}	°C ↗ _{min}	High Temp	→	↘
	°C	min	min	°C	°C			min	min
NP-Primer	600	2	1	600	950	100	950	1:00	0

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