Made in Germany

Manufacturer:

Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)

**KULZER**  
MITSUI CHEMICALS GROUP

# Hera

## HERAVEST<sup>®</sup> Press

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DE

Instructions for use

GB

Mode d'emploi

FR

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ES

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IT

Kullanma talimatı

TR

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**HERAVEST Press is a phosphate bonded, graphite free, special investment material for speed and linear preheating for the whole range of dental press ceramics.**

**Due to a low degree of reaction layer HERAVEST Press is also suitable for lithium disilicate press ceramics, such as IPS e.max Press.**

**Properties**

- very good controllable expansion
- good mixing and flow properties
- smooth surface of the pressed objects
- simple de-vesting by sandblasting with beads
- very thin reaction layer (lithium disilicate)
- possibility of speed heating
- program-controlled preheating

**Storage**

Keep the powder in cool, dry place.  
Do not store the mixing liquid at less than 5°C (41°F) as it is sensitive to frost.

**GB**

**- 6 - Working temperature**

Approx. 22°C (72°F) (room temperature).

**Mixing ratio**

100 g Powder – 22 ml Liquid.

**Mixing ratio for sachets:**

Powder 100 g, Liquid 22 ml.

Powder 200 g, Liquid 44 ml.

For accurate dispensing please use our dispensing syringes.

**Liquid concentration (examples) for 1000 ml liquid**

Concentration (%)	Liquid (ml)	Distilled water (ml)
90	900	100
80	800	200
70	700	300
60	600	400
50	500	500
40	400	600

**Concentrations**

Concentration (%)	Sachet 1 x 100 g		Sachet 2 x 100 g (200 g)	
	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)
100	22	–	44	–
95	21	1	42	2
90	20	2	40	4
85	19	3	38	6
80	18	4	36	8
75	17	5	34	10
70	15	7	30	14
65	14	8	28	16
60	13	9	26	18
55	12	10	24	20
50	11	11	22	22
45	10	12	20	24
40	9	13	18	26

### Expansion control (%)

	Liquid (ml)	Dist. water (ml)
Molar Crowns	75 %	25 %
Pre-Molar Crowns	75 – 80 %	25 – 20 %
Front Crowns	75 – 80 %	25 – 20 %
Veneers	70 %	30 %
Inlays 1-layer and 2 layers	40 – 50 %	60 – 50 %
MOD – Inlays	70 %	30 %
Onlays	85 – max. 90 %	15 – 10 %
Up to 3-unit bridges (IPS e.max Press)	60 – 70 %	40 – 30 %

For Pre-Molar Crowns and Front Crowns attention has to be paid to the fact that better fitting can be obtained with slim and lean preparations at the higher expansion rate of 80 %.

Expansion figures given above are guidelines and have to be adjusted accordingly i.e. due to different preparation models, press ceramic blanks and press temperature.

The figures given by our research and development department are guidelines and can be individually adjusted.

GB

- 7 -

### Stirring

Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.

### Mixing in vacuum

60 seconds. Place in vacuum for 15 seconds but without switching the paddles on.

Mix for 60 seconds. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the casting.

### Working time

Approx. 6 minutes at approx. 22°C (72°F) (Room temperature).

The working time depends on the room temperature. Warmth curtails the working time.

### Investing

**Fill the ring with investment:** only use a vibrator if it is required to improve the flow.

**Avoid high vibration!** This will lead to the formation of bubbles and breakdown of the mixture.

### Pre-drying time

30 minutes after start of mix.

### Placing the ring

1. Remove casting ring former and sprue former base.
2. Straighten bottom of ring (plaster knife / sand paper).
3. After 30 minutes place the ring in the 850°C (1562°F) preheated furnace (speed heating) or in cold furnace (program controlled preheating)

### Preheating

850°C (1562°F).

### Hold-times speed heating

Hold-time of ring once preheating temperature 850°C (1562°F) is attained again.

100 g ring: 60 minutes.

200 g ring: 60 minutes.

## Hold-times and steps of program controlled preheating

Program level	Temperature	Heating rate (°C / min)	Holding time (min) and usable ring size	
			100 g ring	200 g ring
1	ca. 23°C * – 270°C (73°F * – 518°F)	8	30	40
2	270°C – 580°C (518°F – 1076°F)	8	30	40
3	580°C – 850°C (1076°F – 1562°F)	8	60	60

\* Room temperature

## Investing

### Initial De-vesting

Remove with glass beads 50 microns, 4 bar, the investment until short distance before the press objects.  
Do not uncover the objects.

### De-vesting

Before reaching the pressed object reduce the air pressure down to 1.5 bar and proceed with sandblasting.

**Note:** Avoid high speed punctual blasting to prevent overheating!

GB

### - 8 - Note lithium disilicate:

Pressing lithium disilicate material creates a reaction layer with the investment.

To eliminate this reaction layer please refer to the Instructions for use of IPS e.max Press from Ivoclar.

## Precautions

- Do not open the furnace during the wax burnout phase.
- Investments contain quartz. **Do not inhale the dust!** Risk of lung damage (silicosis, lung cancer).  
**Recommendation:** Wear a protective face mask, type FFP 2 EN 149:2001.  
Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- Use an extractor with a fine dust filter when sandblasting.

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials. The results may vary if other materials (e.g. casting ring liners, waxes, equipment etc.) or different techniques are used.

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