Hera

HERAVEST® Speed

| | | |
|----|----------------------------|----|
| DE | Gebrauchsanweisung | 3 |
| GB | Instructions for use | 6 |
| FR | Mode d'emploi | 9 |
| ES | Instrucciones de uso | 12 |
| IT | Istruzioni per l'uso | 15 |
| NL | Gebruiksaanwijzing | 18 |
| SE | Bruksanvisning | 21 |

KULZER
MITSUI CHEMICALS GROUP

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Graphite-free, phosphate-bonded precision investment material for fast preheating for crown and bridgework variably preheatable with and without steel ring.

Storage

Keep the powder in a cool, dry place

Do not store the mixing liquid at less than +5°C (41°F) as it is sensitive to low temperatures.

Working temperature

approx. 22°C (71.6°F) (Room temperature)

Mixing ratio

100 g Powder – 25 ml Liquid

60 g Powder – 15 ml Liquid

Mixing ratio for 160 g Sachet:

Powder 160 g, Liquid 40 ml

For accurate dispensing please use our dispensing syringes.

| Liquid concentration (examples) | | |
|---|------------------|----------------------|
| Liquid concentration for 1000 ml liquid | | |
| Concentration | Concentrate (ml) | Distilled water (ml) |
| 50 % | 500 | 500 |
| 65 % | 650 | 350 |
| 75 % | 750 | 250 |
| 95 % | 950 | 50 |

| Concentrations | | | | | | |
|-----------------|---------------------|------------------|-----------------------------|------------------|-----------------------------|------------------|
| Concentration % | Sachet 1 x 160 g | | Sachet 2 x 160 g (320 g) | | Sachet 3 x 160 g (480 g) | |
| | Liquid (ml) | Dist. water (ml) | Liquid (ml) | Dist. water (ml) | Liquid (ml) | Dist. water (ml) |
| 100 | 40 | – | 80 | – | 120 | – |
| 95 | 38 | 2 | 76 | 4 | 114 | 6 |
| 90 | 36 | 4 | 72 | 8 | 108 | 12 |
| 85 | 34 | 6 | 68 | 12 | 102 | 18 |
| 80 | 32 | 8 | 64 | 16 | 96 | 24 |
| 75 | 30 | 10 | 60 | 20 | 90 | 30 |
| 70 | 28 | 12 | 56 | 24 | 84 | 36 |
| 65 | 26 | 14 | 52 | 28 | 78 | 42 |
| 60 | 24 | 16 | 48 | 32 | 72 | 48 |
| 55 | 22 | 18 | 44 | 36 | 66 | 54 |
| 50 | 20 | 20 | 40 | 40 | 60 | 60 |
| 45 | 18 | 22 | 36 | 44 | 54 | 66 |
| 40 | 16 | 24 | 32 | 48 | 48 | 72 |

| ► Expansion control (%) | | | | | | | |
|---|-----------------|-------------------------------------|--|---|--------------------------------------|---|------------|
| Indication/ Pattern material | | alloy | Casting alloys, high gold content Universal bonding alloys | Casting alloys, reduced gold content Universal bonding alloys | Bonding alloys, high gold content | Bonding alloys, reduced gold content | NEM-alloys |
| | | 1. Inlays, Onlays Partial crowns | - Wax | 50 +/-5 | 50 +/-5 | - | - |
| 2. Single crowns and small-span bridges | - Wax | 65 +/-5 | 65 +/-5 | 75 +/-5 | 75 +/-5 | 90 +/-5 | |
| 3. Outer copings telescopic crowns | - Pattern resin | 85 +/-5 | 85 +/-5 | 90 +/-5 | 90 +/-5 | - | |
| 4. Outer copings tapered crowns | - Pattern-resin | 45 +/-5 | 45 +/-5 | - | - | - | |
| The figures given by our research and development department are guidelines and can be individually adjusted. ► To achieve higher levels of expansion we recommend using Investment Special BS Liquid 3. | | | | | | | |

Stirring

Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.

► Mixing in vacuum

60 seconds

Place in vacuum for 15 seconds but without switching the paddles on. Mix for 60 seconds. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the casting.

Working time

Approx. 5–6 min. at ap-prox. 22°C (Room temp.)

The working time depends on the room temperature. Warmth curtails the working time.

| ► Suitable ring systems | |
|---|---|
| Steel ring: X1 / X3 / X6 / X9 | Casting ring–X1 and X3 steel: We recommend using one dry cellulose ring liner with Vaseline. Casting ring–X6 and X9 steel: We recommend using two dry cellulose ring liners with Vaseline. |
| Ringless: X1 / X3 / X6 | Use Kulzer silicone rings for investment. |

Investing

Fill the ring with investment: only use a vibrator if it is required to improve the flow. Avoid high vibration! This will lead to the formation of bubbles and breakdown of the mixture.

Placing the rings in a preheated furnace

20 minutes after start of mix. Roughen the top of the mould before placing it in the furnace!

| ▶ Preheating Speed preheating | | | | | | |
|--|-------------------|------------------|---|----|----|-----|
| Type of alloy | Set-up temp. (°C) | Final temp. (°C) | Hold-times (min.) (commences once the final temp. is attained again) | | | |
| | | | X1 | X3 | X6 | X9* |
| Casting alloys, high gold content** | 700 | max. 700 | 30 | 30 | 45 | 60 |
| Universal bonding alloys, high gold content** | 750-800 | max. 800 | 30 | 30 | 45 | 60 |
| Casting alloys, reduced gold content** | 700 | max. 700 | 30 | 30 | 45 | 60 |
| Universal bonding alloys, reduced gold content** | 750-800 | max. 800 | 30 | 30 | 45 | 60 |
| Bonding alloys, high gold content | 850-900 | max. 950 | 30 | 30 | 45 | 60 |
| Bonding alloys, reduced gold content | 850-900 | max. 950 | 30 | 30 | 45 | 60 |
| Non-precious alloys | 900-950 | max. 950 | 30 | 30 | 45 | 60 |

*Fit casting ring at 850°C and heat up to the final temperature of 950°C at maximum heating rate.

Transfer technique:

To achieve better casting results and surface qualities with high-gold content alloys (e.g. Bio Maingold SG) and gold-reduced alloys (e.g. Hera SG) with a preheating temperature of 700°C, first preheat the investment compound to 900°C [depending on casting ring size maintain the temperature (see table)] and then transfer to the 700°C hot oven or reduce the temperature from 900°C to 700°C.

Please follow instructions for the respective alloys.

⚠ Safety Instructions

- Do not open the furnace during the wax burnout phase.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer).
Recommendation: Wear a protective face mask, type FFP 2 EN 149:2001.
Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when devesting.
- Use an extractor with a fine dust filter when sandblasting.

⚠ Changes

Revision marking November 2007 / Changes in relation to Issue 02/2005.
Expansion control, vacuum stirring time, suitable casting ring systems, preheating process

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials. The results may vary if other materials (e.g. casting ring liners, waxes, equipment etc.) or different techniques are used.

The figures given by our research and development department are guidelines and can be individually adjusted.

We appreciate your comments, feedback, and suggestions.
Contact in the countries and more information are available by the displayed QR code or at our website www.kulzer-dental.com

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Kulzer GmbH