
- GB** Instructions for use
- USA** Instruments with box-locks
- D** Gebrauchsanweisung  
Instrumente mit Durchsteckschluss
- F** Mode d'emploi  
Instruments à mécanisme de fermeture de sécurité
- E** Instrucciones de manejo  
Instrumental con tapa-pasador
- I** Istruzioni per l'uso  
Strumenti con meccanismo di incastro
- P** Instruções de utilização  
Instrumentos com fecho de inserção
- NL** Gebruiksaanwijzing  
Instrumenten met insteekopening
- S** Bruksanvisning  
Instrument med skänkelgenomföring
- RUS** Инструкция по применению  
Инструменты с соединением бранш типа "box lock"
- CZ** Návod k použití  
Nástroje s uzávěrem prostrčení
- PL** Instrukcja użytkowania  
Instrumenty ze złączem przetykanym
- SK** Návod na používanie  
Nástroje s prestrkávacim uzáverom
- TR** Kullanım Kılavuzu  
Anahtar kilitli aletler

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

TA-Nr. 013729 11/12 V6 Änd.-Nr. 46258

Technical alterations reserved

## Aesculap® Instruments with box-locks

### Symbols on product and packages

Caution, general warning symbol  
Caution, see documentation supplied with the product

### Scope

Multipiece instruments with one or more box-locks with and without industry springs, which, except for the box-lock, have accessible and visible surfaces.

### Intended use

These instructions for use describe an instrument for universal surgical use in all specialist areas.

### Safe handling and preparation

#### CAUTION

##### Federal law restricts this device to sale by, or on order of a physician!

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.

### Validated reprocessing procedure

#### General safety instructions

##### Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

##### Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

##### Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

##### Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

##### Note

If there is no final sterilization, then a virucidal disinfectant must be used.

##### Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at [www.aesculap-extra.net](http://www.aesculap-extra.net)

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

#### General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to Publications, Red Brochure – Proper maintenance of instruments.

#### Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

## Cleaning/disinfection

### Product-specific safety notes on the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning agents and disinfectants according to the manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 55 °C.

### Validated cleaning and disinfection procedure

| Validated procedure                                   | Specific requirements   | Reference  |
|---|---|--|
| Manual cleaning with immersion disinfection           | <ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>          | Chapter Manual cleaning/disinfection and sub-chapter:<br><ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with immersion disinfection</li> </ul>               |
| Mechanical alkaline cleaning and thermal disinfection | <ul style="list-style-type: none"> <li>■ Mount jaws protection on the product.</li> <li>■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> </ul> | Chapter Mechanical cleaning/disinfecting and sub-chapter:<br><ul style="list-style-type: none"> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul> |

### Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

### Manual cleaning with immersion disinfection

| Phase | Step                     | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical   |
|-------|--------------------------|--------------|------------|--------------|------------------|--|
| I     | Disinfecting<br>cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| II    | Intermediate<br>rinse    | RT (cold)    | 1          | -            | D-W              | -  |
| III   | Disinfection             | RT (cold)    | 15         | 2            | D-W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| IV    | Final rinse              | RT (cold)    | 1          | -            | FD-W             | -  |
| V     | Drying                   | RT           | -          | -            | -                | -  |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

#### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

## Mechanical cleaning/disinfecting

### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

## Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                      | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical/Note   |
|-------|---------------------------|--------------|------------|------------------|---|
| I     | Prerinse                  | <25/77       | 3          | D-W              | -   |
| II    | Cleaning                  | 55/131       | 10         | FD-W             | <ul style="list-style-type: none"><li>Concentrate, alkaline:</li><li>- pH = 13</li><li>- &lt;5 % anionic surfactant</li><li>0.5 % working solution</li><li>- pH = 11*</li></ul> |
| III   | Intermediate<br>rinse     | >10/50       | 1          | FD-W             | -   |
| IV    | Thermal disin-<br>fecting | 90/194       | 5          | FD-W             | -   |
| V     | Drying                    | -            | -          | -                | According to the program for cleaning and disinfection device   |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

## Inspection, maintenance and checks

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Check for compatibility with associated products.

## Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

## Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
  - Steam sterilization through fractionated vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

## Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10<sup>-6</sup>, Aesculap recommends the following parameters:

| Aesculap Orga Tray/Sterile container (perforated bottom) |               |       |                     |
|--|---------------|-------|---------------------|
| Minimum cycle parameters*                                |               |       |                     |
| Sterilization method                                     | Temp.         | Time  | Minimum drying time |
| Prevacuum  | 270 °F/275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

## Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

## Technical Service

**Risk of injury and/or malfunction!**

- ▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

### Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 (7461) 95-1602  
Fax: +49 (7461) 16-5621  
E-Mail: ats@aesculap.de  
Or in the US:  
Aesculap Inc.  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042  
Aesculap Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420  
Other service addresses can be obtained from the address indicated above.

## Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

## Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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