
Указания за работа с материала	1	BG
Pokyny pro zpracování	5	CZ
Verarbeitungshinweise	9	D
Bearbejdelsesindikationer	13	DK
Consejos para su utilización	17	E
Indications de mise en œuvre	21	F
Instructions for use	25	GB
Υποδείξεις επεξεργασίας	29	GR
Feldolgozási útmutató	33	H
Napomene za obradu	37	HR
Note per la lavorazione	41	I
Apdorojimo nuorodos	45	LT
Norādes sakarā ar pārstrādi	49	LV
Instructies voor de bewerking	53	NL
Instruções de utilização	57	P
Informacje dotyczące obróbki	61	PL
Instrucțiuni de prelucrare	65	RO
Рекомендации по работе с материалом	69	RUS
Beredningsanvisningar	73	S
İşlem talimatları	77	TR
تعليمات بخصوص التحضير	81	UAE
加工说明	85	RC
加工について	89	J
가공처리 주지사항	93	ROK

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LFC

### Indications for use

- LFC is a hydrothermal dental ceramic specially formulated for low fusion ceramic repairs and corrections to crowns and bridges veneered with Duceram® Plus or Duceram® Kiss, as well as for inserting inlays, onlays, and facets on the master cast.
- For dental use only

### Contraindications

- Suitable only for the indications listed above

### Warnings

Adverse effects from these medical products are very rarely observed if they are correctly processed and used. However, immune reactions (e.g. allergies) and/or localised abnormal sensations (e.g. taste irritation or mucosal stimulation) cannot be a priori excluded absolutely. We ask you to contact us in the event adverse effects occur, even those you may be doubtful about. This medical product must either not be used or used only under the close supervision of the physician/dentist in the case of patient hypersensitivity to LFC for veneers or to any of its components.

Known cross-reactions or interactions between this medical product and other medical products or agents currently being used orally must be taken into consideration by the physician/dentist when this medical product is being used. Please provide all of the information listed above to the clinician involved if you are processing this medical product for a tailored restoration.

- Do not inhale grinding dust

### Precautions

Follow the Directions for Use and the Safety Data Sheet when using the product.

### Adverse reactions

We are not aware of any risks and/or side effects related to LFC for veneers.

### Technical data

- CTE – dentin: 11.1  $\mu\text{m}/\text{m}\cdot\text{K}$  (25–400 °C)
- Dental ceramics, type 1, Class 2–8 according to DIN EN ISO 6872
- Flexural strength and chemical solubility according to DIN EN ISO 6872

### Transport and storage conditions

- Protect liquids from freezing. Store containers tightly sealed at a temperature above 10 °C. 10 °C ↓ Lower temperature limit
- Protect powder from light and moisture and store vibration-free.
  - Keep dry
  - Keep away from sunlight

### Please pay attention to the following symbols on the product labels

- REF Catalogue number
- LOT Batch code
- Use by
- Pay attention to instructions for use

### Combinable liquids

- Bonder/powder opaques:
  - Ducera® Liquid SD
  - Ducera® Liquid Form
- Shoulder porcelain:
  - Ducera® Liquid Quick
- Dentins/Incisals etc.:
  - Ducera® Liquid SD
  - Ducera® Liquid Form
- Stains/Glaze:
  - Ducera® Liquid Stain improved
- Isolation:
  - Ducera® Sep Isolating Fluid

### Ceramic furnace

To ensure superior results you must ensure that the required firing temperatures and times are reached. If necessary, you may have to adjust the furnace parameters to get the correct setting.

### 1 2 3 General ceramic repairs and corrections

Abrade the metal framework carefully with aluminium oxide 110–150 µm. Then clean the bridge with care. LFC Opaquer should never be mixed with opaquer fluid but only with Modelling Liquid. The exposed metal parts are coated with LFC Opaquer. Fire in accordance with the LFC wash firing directions (Fig. 1).

Next, layer with LFC materials (Fig. 2). LFC materials can be applied directly to the smooth surface and fired without prior surface treatment.

Lightly smooth and polish the transition zones with a rubber polisher and diamond paste. A bridge repaired with LFC after the firing procedure (Fig. 3).

### 4 5 6 Corrections and surface finishing of ceramic shoulders

Enhanced precision of fit can be obtained also through post-glazing application of LFC or SM-L materials. All shade corrections should be completed before using these materials. Use only SEP Isolating Liquid to separate the die. The die must not be sealed; this reduces ease of separation.

LFC materials selected according to the LFC indicator are dabbed with the specified liquids and applied to the shoulder in need of correction (Fig. 4). The workpiece is mounted on the master cast using gentle tapping and concomitant suctioning. Once the LFC material is well dried, the workpiece is removed carefully from the cast and fired according to the directions for firing (Fig. 5).

The transition zone between the LFC and the metal-veneer ceramic should be finished with a rubber polisher and then smoothed and polished with diamond paste.

Finished crown (Fig. 6).

### 7 8 9 10 11 12 13 14 15 16

#### Ceramic inlays, onlays, facets

To fabricate inlays, onlays and facets, model systems are available that use state-of-the-art methods to enable FFM dies to be mounted on the master cast. Using LFC materials, however, a method has been developed that permits one to skip the use of such specialized model fabrication systems. One can fabricate ceramic inlays and facets by a simple method which is both intuitive and economical.

Because of their brilliant light diffusion and extremely homogeneous surfaces, LFC materials are setting the standard in aesthetics, occlusal repairs, and periodontic acceptability.

Once the duplication model has been fabricated with Ducera® Lay Superfit die material, a ceramic veneer of Duceram® Kiss is fired on the die, and this is then transferred to the master cast. When using LFCs, layer the workpiece to the desired form and shade, remove it from the master cast, and fire it on furnace wadding.

Based on the wide differences in firing patterns between Duceram® Plus or Duceram® Kiss and LFC, the close fit of the ceramic veneer is not affected at all, even with geometrically awkward inlays.

The work sequence in detail:

Prepare the master die.

Develop an appropriate cementation gap on the master die with Ducera® Lay Spacer. The spacer should come no closer than approx. 1 mm to the preparation margin. Block out any undercuts (Fig. 7).

## Instructions for use

Duplicate the die with high quality duplication material, e.g. from DeguDent (Fig. 8). Once you have carefully devested the duplication model from the master die, pour it with Ducera® Lay Superfit die material. The investment compound die should be kept as small as possible. Make sure to observe the mixing ratios specified in the Work Instructions for the Ducera® Lay Die Material.

Carefully devest the Ducera® Lay Superfit die. One will find the hardening time and a description of the sinter firing procedure in the Instructions for Use for Ducera® Lay Superfit.

Once the sinter firing is complete and the Ducera® Lay Superfit die has cooled, liberally apply Connector Paste except on the preparation margins, where a light application is sufficient to prevent possible lifting of the margins during subsequent firings (Fig. 9).

The Connector firing that follows is intended to seal the Ducera® Lay Superfit die and provide a good bond between the die and the ceramic veneer to follow.

The Connector firing is usually performed twice using the same firing protocol and can be repeated until the Connector surface is uniformly sealed and glossy (Fig. 10). When performing the Connector firing, make sure to follow the Instructions for Use for the Ducera® Lay Superfit or Connector Paste products.

Following the Connector firing, cover the cooled die with a base layer of the appropriate ceramic (Fig. 11). Be aware that the refractory die absorbs a very high amount of heat. Because of this, elevate the **firing temperature by about 30 °C** to ensure that the veneer ceramic is thoroughly sintered. **The temperature elevation rate should be at least 99 °C/min.**

On completion of firing, compensate for shrinkage loss until the ceramic layer is flawless. Observe the firing protocol. Make sure that the veneer ceramic foundation for the high fusion ceramic is properly dimensioned. Next, carefully abrade the Ducera® Lay Superfit with 25–50 µm glass beads at approx. 1.5 bar and check for fit on the master cast (Fig. 12).

Separate the model, if necessary, using SEP Isolation Liquid only. Then carefully mount it again on the master cast and complete the contouring with LFC materials directly applied to the master cast (Fig. 13). The ceramic layer should be bonded with LFC die adhesive since this is the only way to ensure residue-free firing of the adhesive (a smidgen of adhesive is fully adequate).

Next, carefully detach the restoration from the master cast and set it on furnace wadding for firing (Fig. 14).

Corrections may be performed in the manner described above (Fig. 15).

On completion of the finishing procedure, make the ceramic workpiece ready for glaze firing (Fig. 16). Additional surface sealing will be provided by the glazing materials.

**General Firing Protocol for Ceramic Corrections (core firing or base layer)**

	Pre-heat temp. (°C)	Drying time (min.)	Heat-up rate (°C/min.)	Firing temp. (°C)	Hold time (min.)	Vacuum (hPa)
Wash firing	450	1	55	700	1	50
Dentine firing 1	450	6	55	680	1	50
Dentine firing 2	450	3	55	670	1	50
Glaze firing	450	2	55	650	1-2	–

**General Firing Protocol for Inlays and Veneers (core firing or base layer)**

	Pre-heat temp. (°C)	Drying time (min.)	Heat-up rate (°C/min.)	Firing temp. (°C)	Hold time (min.)	Vacuum (hPa)
Sinter firing	575	5	55	1100	5	–
Connector firing	620	4	55	970	1	50
Base layer	575	6	99	940	1	50
LFC dentine 1*	450	6	55	680	1	50
LFC dentine 2*	450	3	55	670	1	50
LFC glaze*	450	2	55	650	1-2	–

\* These firing data apply after removal from the refractory die.

In general, one must use or adapt both the tempering and stress relief cooling times for the particular alloy. For further information on firing parameters, see the Abbreviated Instructions or the Instructions for Use for the particular Kiss ceramic.

The settings given here are guidelines and pointers only. Variations in firing results may occur. Firing results depend on the performance of a specific furnace and can vary by manufacturer and furnace age. Therefore, guidelines must be adapted each time to the particular firing. We recommend a test firing to check the furnace. We carefully prepare and check all data but we cannot issue guarantees.

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DeguDent GmbH