## OPTIGLAZE™ color

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Use only by a dental professional in the recommended indications.

### RECOMMENDED INDICATIONS

1. Characterization of direct & indirect composite restorations, acrylic denture base and acrylic denture relines. 2. For obtaining surface smoothness and wear resistance of restorations made of composite resin, acrylic denture base and artificial acrylic teeth. 3. For polishing on glass ionomer cements. Not intended for characterization of natural teeth.

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### DIRECTIONS FOR USE

Surface characterization: 1) After contouring the resin surface, finish with a carbide bur or coarse silicone point. The oxygen inhibition layer must be removed as it inhibits the polymerization of OPTIGLAZE color. Laboratory use: Sandblast with 25-50µm alumina (0.15MPa/1.5bar), clean and dry. NOTE: 1) If surface still contains oil residue, the surface can be cleaned further with ethanol (alcohol). The use of ethanol (alcohol) on acrylic resin may cause micro cracks. A mild detergent can be used to clean the acrylic resin. 2) When applying OPTIGLAZE color on hybrid ceramic blocks such as Clear Clear HV, use CERAMIC PRIMER II. 2) Shake the bottle of OPTIGLAZE color well. Dispense a few drops in a dispensing dish. Dilute by adding color to the desired shade. Apply a thin layer of dentin primer or primer. Do not air blow. 3) Light cure with a suitable light curing device. (See table 1) Place the light curing device as close as possible to the surface. Laboratory use: Set a 2cm platform in the curing device and place the prosthesis on the surface facing the light. Reverse the prosthesis and repeat the curing. Repeat the light curing to fully cure the undercut or area in shadow.

Internal characterization: 1) Build up the dentin layer with composite and light cure. Follow the instructions of your composite system.

Do not remove the oxygen inhibition layer. Apply OPTIGLAZE color on the oxygen inhibition layer. If the oxygen inhibition layer is removed, treat the surface with CERAMIC PRIMER II (blow dry) and COMPOSITE PRIMER (light cure) before applying OPTIGLAZE color. 1) Shake the bottle of OPTIGLAZE color well. Dispense a few drops in a disposable dispensing dish. Apply a thin layer of OPTIGLAZE color of choice to the dentin surface with an instrument or brush (e.g. the convex and concave of dentin margin or incisal area). Mix multiple colors to obtain the desired shade if needed. Do not air blow. 2) Light cure with a suitable light curing device. (See table 1) 4) Continue building up with enamel shade of the oxygen inhibition layer. Follow the instructions of your composite system. 5) Apply the finish with CERAMIC CLEAR or CLEAR HV shade to the surface.

Repair of restoration (if the gloss is lost): 1) Roughen the surface of restoration using a carbide bur and paper cone 2) Cleanse and finish with a steam cleaner or ultrasonic cleaner and then dry with free air. NOTE: 1) Resin is contaminated by oil residue, the surface can be cleaned further with ethanol (alcohol). The use of ethanol (alcohol) on acrylic resin may cause micro cracks. A mild detergent can be used to clean the acrylic resin. 2) Apply the finish with CERAMIC CLEAR or CLEAR HV shade to the surface. Do not air blow. 3) Light cure with a suitable light curing device. (See table 1)

Table 1: Light Cure Time	90 sec.
GC LABLIGHT DUO	90 sec.
GC LABLIGHT LV-III	5 min.
GC STEPLIGHT SL-I	Distance from light source <3cm 10 sec.
	Distance from light source >3cm 20 sec.
Halogen Light	
Plasma Arc	4 sec.
LED Light (with light source having wavelength lower than 430nm)	8 sec.
Table 2: Depth of Cure	
Orange, Yellow, Red/brown, Red	0.07 mm
A Plus, B Plus, C Plus, White, Ivory White, Pink Orange, Pink, Clear, Lavender, Grey, Blue	0.15 mm
Clear, Clear HV	0.5 mm

**STORAGE**  
Recommended for optimal performance, store at room temperature of 4-25°C (39-77°F).

**PACKAGE:**  
Complete set:  
A Plus, B Plus, C Plus, White, Ivory White, Pink Orange, Pink, Clear, Lavender, Grey, Blue  
Clear HV, Clear HV 5ml, jednorazová dávková miska, plochá štětka, křída

**CAUTION**  
1. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. 2. Particular avoid contact of uncured material with skin. 3. Avoid getting material on clothing. 4. OPTIGLAZE color is highly flammable liquid and vapor. DO not use near open flame. Keep away from source of ignition. 5. OPTIGLAZE color is volatile and may cause respiratory irritation. Use in a well-ventilated space. Replace cap immediately after dispensing. 6. It is not recommended to use OPTIGLAZE color on high pressure bearing area. 7. Avoid contact with eyes, flush immediately with water and seek medical attention. 8. In case of contact with skin, immediately remove with a sponge or cotton ball. 9. The oxygen inhibited layer of the resin surface to be coated should be removed before applying OPTIGLAZE color as it may hinder the proper curing of the product. 10. Dispense a few drops of OPTIGLAZE color into a disposable dispensing dish. Do not dispense directly from the syringe. 11. Do not use OPTIGLAZE color on volatile solvent. 12. Store locked up. 13. Dispose all waste according to local regulation. 14. Do not use the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**CLEANING AND DISINFECTING**  
1. Use disinfectant to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. 2. Rinse essential to avoid cross-contamination. 3. Do not IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level disinfectant. 4. Rinse essential to avoid cross-contamination. 5. Follow regional/national guidelines.

Some products referred to in the present IVD may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.geplore.com>

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They can also be obtained from your supplier.

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events during the use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country according to the following link.

EN  
Prior to use, carefully read the instructions for use.

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## OPTIGLAZE™ color

LIGHT-CURED CHARACTERIZATION MATERIAL

For use only by a dental professional in the recommended indications.

### RECOMMENDED INDICATIONS

- Characterization of direct composite restorations, acrylic denture base and artificial acrylic teeth.
- For obtaining surface smoothness and wear resistance of restorations made of composite resin, acrylic denture base and artificial acrylic teeth.

Note:  
Not intended for characterization of natural teeth.

### CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to acrylate / methacrylate monomers or acrylate / methacrylate polymers.

### DIRECTIONS FOR USE

Surface characterization  
1. After contouring the resin surface, finish with a carbide bur or coarse silicone point. The oxygen inhibition layer must be removed as it inhibits the polymerization of OPTIGLAZE color.  
Laboratory use:  
Sandblast with 25-50 µm alumina (0.15 MPa/1.5 bar), clean and dry.

- If surface still contains oil residue, the surface can be cleaned further with ethanol (alcohol). The use of ethanol (alcohol) on acrylic resin may cause micro cracks. A mild detergent can be used to clean the acrylic resin.
- When applying OPTIGLAZE color on hybrid ceramic blocks such as CERASMART, treat the surface using CERAMIC PRIMER II.
- Shake the bottle of OPTIGLAZE color well. Dispense a few drops in a dispensing dish. Dilute by adding CLEAR to the desired shade. Apply a thin layer to the resin surface using the brush provided. Do not air blow.
- Light cure with a suitable light curing device. (See table 1) Place the light tip as close as possible to the surface.  
Laboratory use:  
Set a 2 cm platform in the curing device and place the prosthesis on it with the coated surface facing the light. Reverse the prosthesis and repeat the light curing to fully cure the undercut or area in shadow.

### Internal characterization

- Build up the dentin layer with composite and light cure. Follow the instructions of your composite system.  
Note:  
Do not remove the oxygen inhibition layer. Apply OPTIGLAZE color on the oxygen inhibition layer. If the oxygen inhibition layer is removed, treat the surface with CERAMIC PRIMER II (flow dry) and COMPOSITE PRIMER (light cure) before applying OPTIGLAZE color.
- Shake the bottle of OPTIGLAZE color well. Dispense a few drops in a disposable dispensing dish. Apply a thin layer of OPTIGLAZE color of choice to the dentin surface with an instrument or brush (e.g. the convex and concave of dentin mamon or incisal areas). Mix multiple colors to obtain the desired shade if needed. Do not air blow.
- Light cure with a suitable light curing device. (See table 1)
- Continue building up with enamel shade on the oxygen inhibition layer. Follow the instructions of your composite system.
- Finish and polish using standard technique or apply OPTIGLAZE color CLEAR or CLEAR HV on the entire surface.

### Repair of restoration (e.g. when the glass is lost)

- Roughen the surface of restoration using a carbide bur and paper the face.
- Clean the surface with a steam cleaner or ultrasonic cleaner and then dry with oil free air.

Note:  
If surface is contaminated by oil residue, the surface can be cleaned further with ethanol (alcohol). The use of ethanol (alcohol) on acrylic resin may cause micro cracks. A mild detergent can be used to clean the acrylic resin.

- Apply OPTIGLAZE color CLEAR or CLEAR HV shade to the surface. Do not air blow.
- Light cure with a suitable light curing device. (See table 1)

Table 1: Light Curing Time	
GC LABOLIGHT DUO	90 sec.
GC LABOLIGHT LV-III II	5 min.
GC STEPLIGHT SL-I	Distance from light source <3 cm 10 sec. Distance from light source >3 cm 20 sec.
Halogen Light	40 sec.
Plasma Arc	8 sec.
LED Light (with light source having wavelength lower than 430 nm)	40 sec.
Table 2: Depth of Cure	
Orange, Yellow, Red/Brown, Red	0.07 mm
A Plus, B Plus, C Plus, White, Ivory White, Pink Orange, Pink, Olive, Lavender, Grey, Blue	0.15 mm
Clear, Clear HV	0.5 mm

**STORAGE**  
Recommended for optimal performance, store at room temperature of (4-25°C) (39-77.0°F).

**PACKAGE:**  
Complete set:  
A Plus 2.6 mL, B Plus 2.6 mL, C Plus 2.6 mL, White 2.6 mL, Ivory White 2.6 mL, Yellow 2.6 mL, Orange 2.6 mL, Pink Orange 2.6 mL, Pink 2.6 mL, Red Brown 2.6 mL, Olive 2.6 mL, Lavender 2.6 mL, Grey 2.6 mL, Blue 2.6 mL, Red 2.6 mL, Clear 5 mL, Clear HV 5 mL, Disposable Dispensing Dish, Flat Brush, Round Brush.

- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In particular avoid contact of material with skin.
- Avoid getting material on clothing.
- OPTIGLAZE color is highly flammable liquid and vapor. Do not use near open flame. Keep away from sources of ignition.
- OPTIGLAZE color is volatile and may cause respiratory irritation. Use in a well ventilated space. Replace cap immediately after dispensing.
- It is not recommended to use OPTIGLAZE color on high pressure bearing area.
- In case of contact with skin, flush immediately with water and seek medical attention.
- In case of contact with eyes, immediately remove with a sponge or cotton pellet. Flush with water.
- The oxygen inhibited layer of the resin surface to be coated should not be removed before applying OPTIGLAZE color as it may hinder the proper curing of OPTIGLAZE color.
- Dispense a few drops of OPTIGLAZE color into a disposable dispensing dish immediately before applying as the material contains a volatile solvent.
- Store locked up.
- Dispose all waste according to local regulation.
- In rare cases the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**CLEANING AND DISINFECTING:**  
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection.  
Immediately after use inspect device and label for deterioration. Discard device if damaged.  
DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:  
<http://www.gceurope.com> or  
<http://www.gcamerica.com>  
They can also be obtained from your supplier.

Undesired effects- Reporting:  
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:  
[https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en) as well as to our internal vigilance system: [vigilance@gc.dental](mailto:vigilance@gc.dental)  
In this way you will contribute to improve the safety of this product.

Last revised: 07/2019

Rx Only

MANUFACTURED BY  
GC DENTAL PRODUCTS CORP.  
2-263 Toritama-cho, Kasugai, Aichi 486-0841, Japan  
RESPONSIBLE MANUFACTURER IN CANADA:  
GC AMERICA INC.  
3737 West 127th Street, Atlap, IL 60603 U.S.A.

DISTRIBUITO DA  
GC CORPORATION  
76-1 Hasunuma-cho, Tabashi-ku, Tokyo 174-8566, Japan  
EU: GC EUROPE NV  
Renaatweg 9, Leuven-Louvain 1200, Belgium  
TEL: +32 16 70 10 90

GC AMERICA INC.  
3737 West 127th Street, Atlap, IL 60603 U.S.A.  
TEL: +1 708 599 1000  
GC SOUTH AMERICA  
Av. Santa Rosa 1500, P.O. Box 57, BRASIL  
CEP: 02022-051 - TEL: +55 11 2925-0905  
FAX: +55 11 2925-0906  
RESP: TEL: Mayra de Santos Ribeiro - CROSP 105.1682

GC ASIA DENTAL PTE. LTD.  
5 Tampines Central 1, #06-01 Tampines Plaza, Singapore 529541  
TEL: +65 6546 7588  
GC AUSTRALASIA DENTAL PTY. LTD.  
173 Boscawen Road, Balaclava, NSW 2016, Australia  
TEL: +61 2 9391 6200

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