Oscillating Instruments

DE

Oszillierende Instrumente

Intensiv Bevelshape
Intensiv Cavishape
Intensiv Proxoshape
Intensiv Proxoshape Coarse
Intensiv Proxoshape Flexible
Intensiv Rootshape

Zu verwenden in Kombination mit:

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To be used in combination with:

- Intensiv Swingle, WG-69 LT Cons, with light
- Intensiv Swingle, WG-69 A Cons, without light

Alternatively:

- KaVo 61 LRG or 61 LG
- KaVo contra-angle E61C

Files

Intensiv Bevelshape:

Margin beveling in cavity, crown and chamfer preparations
Intensiv Swingle WG-69 LT / KaVo 61 LRG / KaVo 61 LG / KaVo E61C

Intensiv Cavishape:

Contouring and finishing of cavity margins
Intensiv Swingle WG-69 LT / KaVo E61C / KaVo 61LG

Intensiv Proxoshape:

Contouring, finishing and pre-polishing of proximal filling surfaces
Intensiv Swingle WG-69 LT / KaVo 61 LRG / KaVo 61 LG / KaVo E61C

Intensiv Proxoshape Coarse:

Modelling of interproximal aesthetic restorations
Intensiv Swingle WG-69 LT / KaVo E61C / KaVo 61LG

Intensiv Proxoshape Flexible:

Finishing and polishing of aesthetic restorations with difficult to access proximal areas. Disposal of excess and over hangings filling
Intensiv Swingle WG-69 LT / KaVo 61 LRG / KaVo 61 LG / KaVo E61C

Intensiv Rootshape:

Rootplaning and odontoplasty
Intensiv Swingle WG-69 LT / KaVo 61 LRG / KaVo 61 LG / KaVo E61C

Recommendations for use

- Set the instrument to the recommended speed before applying it to the tooth.
- Do not stop the instrument in the interdental space.
- In case of very tight contact points: separate the interdental space by placing a wooden wedge cervically.
- Insert the instrument as far as possible.
- Set the micro motor at 40'000 RPM/min max. With Intensiv Swingle this corresponds to 20'000 oscillations per minute operation. (2:1).
- Provide an adequate water spray cooling (minimum 50ml/min.).
- Regulate water flow in case of water excess or shortage.
- Pay attention to applied load (around 50 grams): the instrument must not be bent.
- Remove the strip from the interdental space whilst the contra-angle is moving.

Maintenance and Sterilization

- Clean the instruments and remove debris after each use, so as to maintain their abrasive properties.
- Instruments are delivered non sterile. They must be sterilized prior to first use and disinfected, cleaned, and sterilized after each use.
- Use a single-pulsed or fractionated vacuum autoclave and sub-vacuum drying. Chemiclave sterilizers may also be used. Hot air sterilizers are not suitable for diamond instruments.
- Sterilization must be carried out according to validated procedures.
- Concerning the sterilization process we refer to the ISO standard 17664.
We suggest to follow the indications below:
 - Cycles at 134°C:
Tmin = 134°C – Tmax = 138°C
Pressure = 3.15 bar abs
Time = 4 min (raisable)
 - Cycles at 121°C:
Tmin = 121°C – Tmax = 125°C
Pressure = 2.10 bar abs
Time = 16 min (raisable)

Risk Warnings

- Avoid bending and levering over the limit: use of excessive force (strong levering effect with the contra-angle) could lead to the breakage of the instrument.
- Avoid jamming or levering actions when oscillating, as this increases the risk of instrument breakage.
- Oscillating over 20'000 oscillations per minute could cause breakage of the instrument.
- Avoid temperatures above 180°C which may affect the durability of the instrument.
- Avoid applying excessive loads that would lead to the bending of the instrument, as this could cause damage to the instrument itself.
- To ensure traceability of the instruments during their entire application, we recommend keeping the packaging.

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