Palodent® V3

Sectional Matrix System

DIRECTIONS FOR USE ENGLISH 2

Système de Matrices Sectorielles

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Sistema de matrizes seccionais

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Sectional Matrix System

DIRECTIONS FOR USE - ENGLISH

Caution: For dental use only. USA: Rx only.

1. PRODUCT DESCRIPTION

The Palodent® V3 Sectional Matrix System employs a nickel-titanium retaining ring, interproximal wedge, WedgeGuard and sectional matrices for placement of restorations in the posterior region. Sectional matrices are suitable for use with all posterior restorative materials. Palodent® V3 EZ Coat matrices have a micro-thin, non-stick finish that minimizes the chance that bonding agent or composite sticks to the band, making matrix removal easier. The retaining rings are available in narrow, for most premolar and small molars or universal, for larger molars. Rings may be used either singly or in tandem, for simultaneous restoration of multiple surfaces. The advantages of the sectional matrices and rings include natural contours for better control of contact areas and embrasures, ease of placement, better visualization of the operative field, and more comfort for the practitioner and patient. Eliminating the leverage of a conventional matrix retainer allows the practitioner to use less tension, avoiding excessive compression of tooth and easier formation of anatomical proximal contact areas. When placed prior to preparation, the innovative WedgeGuard combination wedge and shield protects approximating surfaces from inadvertent damage during preparation. Using the pin tweezers, the shield easily slides away, leaving the wedge securely in place.

1.1 Delivery forms (Some delivery forms may not be available in all countries)

- Sectional matrices are available in regular polished stainless steel or color-coded, non-stick Palodent® V3 EZ Coat in 5 sizes: 3.5mm, 4.5mm, 5.5mm, 6.5mm and 7.5mm. These matrices are suitable for use with all posterior restorative materials.
- Retaining rings are available in narrow (dark blue) and universal (light blue) sizes suitable for use in the posterior region either singly or in tandem.
- Wedges and WedgeGuards are available in small (dark blue), medium (medium blue) and large (light blue) sizes.

1.2 Composition

- Sectional Matrices and WedgeGuard shield: Polished stainless steel;
Palodent® V3 EZ Coat matrices: Whitford 8820HR-coated polished stainless steel
- Retaining rings: Nickel-Titanium with Glass Fiber Reinforced Tines
- Wedge and WedgeGuard: Polyacetal Thermoplastic
- Placement forceps and pin tweezers: 420 Stainless steel

1.3 Indications

- The Palodent® V3 System is indicated for Class II restorations to facilitate the contouring of a direct placement restorative materials.
- The Palodent® V3 WedgeGuard is indicated for pre-operative placement to protect surfaces adjacent to Class II, Class III, veneer and crown preparations.

1.4 Contraindications

None known.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

- The Palodent® V3 System consists of small, sharp parts. As with all intraoral instruments, care must be exercised to ensure stable handling and placement. Placement of rubber dam or other suitable isolation is advised prior to operative procedures. Wear suitable protective eyewear and gloves.

Aspiration: If aspiration of matrix, wedge, WedgeGuard, or retaining ring occurs seek immediate medical attention.

- Ingestion:** If accidental swallowing occurs, encourage patient to seek medical assistance to locate swallowed part. If nausea or illness develop, seek medical attention immediately.
- Laceration:** The Palodent® V3 matrices and WedgeGuard shields are sharp, and may cut unprotected skin. Handling matrices and WedgeGuard shields with gloved hands and/or pin tweezers or similar instrument is recommended. In the event a clean matrix or WedgeGuard shield causes a skin laceration, employ usual first aid measures. Discard and properly dispose matrix or WedgeGuard shield.
- Used matrices and accessories are contaminated. Follow proper exposure control plans. In the event of an accidental exposure, follow recommended post-exposure procedures. Following use of the Palodent® V3 System, proper handling and disposal techniques are required. Used matrices are contaminated and sharp. When handling or disposing, use precautions as when handling or disposing other contaminated sharps.
- Care should be taken to have good control on any instrument to protect patients from injury.
- When properly placed, retaining rings exert a clamping force for stabilization. Use of narrow rings on large prepared molars, or use of any ring on teeth compromised by extensive preparation and/or with unsupported tooth structure may lead to tooth crazing or fracture. For compromised teeth, consider alternative matrix retainer.

2.2 Precautions

- This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
- Care must be taken when placing the retaining ring, due to the spring nature of the ring. Ensure ring is fully seated in forceps before placement. Before withdrawing forceps, check that the ring is seated securely apically (see Step-by-step Instructions). Tie dental floss to the Palodent® V3 ring for added security.
- Excessive pressure on matrix, wedge, WedgeGuard, or retaining ring in an apical direction may result in injury to periodontal tissues.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Do not use accessory tools other than those designed, supplied, and/or outlined for use with the Palodent® V3 System. Injury may result. (See Step-by-Step Instructions)

2.3 Adverse reactions

- Product misuse may cause failure resulting in respiratory, digestive, eye, skin, or mucosa injury. (See Contraindications, Warnings and Precautions)
- Failure to properly employ accepted handling and disposal techniques may present biohazard. (See Warnings and Precautions)
- Use of narrow retaining rings on large teeth, or use of any retaining ring on weakened or unsupported tooth structure may lead to tooth fracture (See Warnings)

2.4 Storage conditions

Store at temperatures between 10-24°C/50-75°F. away from moisture or excessive humidity. Keep out of sunlight. Keep away from flame or source of ignition.

3. STEP-BY-STEP INSTRUCTIONS

1. Pre-operative Sterilization

Matrices, wedges and WedgeGuards are intended for single use. If desired, these components may be steam autoclaved one time (134°C, 3 minutes minimum holding time) prior to clinical use. Retaining rings, forceps and pin tweezers are not supplied sterile, and should be steam autoclaved prior to first and each subsequent use (see Hygiene section).

2. Place rubber dam or other suitable isolation technique.

3. Pre-operative Placement (optional)

1. WedgeGuard placement

Prior to preparation, select the appropriate sized WedgeGuard, and place interproximally. Palodent® V3 WedgeGuards contain a placement tab with hole, which is securely grasped by the pin tweezers. Do not use excessive force.

OR

2. Ring Placement

To effectively use a contoured sectional matrix, there must be enough interproximal space to position the matrix without force or distortion. The retaining ring, when placed prior to and left on during cavity preparation, provides the necessary separation with no additional chair time. **NOTE:** Do not place retaining ring pre-operatively if WedgeGuard is used. Follow these ring placement instructions after preparation completion, shield removal and matrix placement, (step 7).

To place the ring, grasp it within at its widest diameter with the provided forceps, being sure the ring rests within the channel on the forceps beaks. Spread the forceps to expand ring only enough to clear tooth height of contour. Do not over-expand. Over-expansion may result in permanent ring deformation. The locking slide on the forceps may be engaged to facilitate carrying and placement of the ring.

Place the ring with one tine in each interproximal space from buccal to lingual adjacent to the surface to be restored. The direction the ring points, mesial or distal, does not matter, though placing it so that the ring points mesially is usually more easily achieved due to easier access.

With the tines properly positioned interproximally, release the locking slide, and slowly release the forceps tension. When ring is securely seated, release the forceps completely and withdraw forceps.

Small parts under tension

- Be sure ring is securely seated in the forceps before introducing into the patient field.
- Be sure ring is securely seated apically before withdrawing forceps.

4. Cavity Preparation

Complete cavity preparation as usual. Once the preparation is complete, remove the ring or the metal shield of the WedgeGuard (if applicable). Use the pin tweezers to grasp the hole on the WedgeGuard shield, and slide shield laterally to remove. The wedge may remain in place if apical to the gingival cavosurface margin. Properly dispose of the shield with contaminated sharps.

5. Matrix Selection

The 5.5mm Palodent® V3 contoured sectional matrix is recommended for most applications. The 4.5mm and 3.5mm are designed for use with adolescents and patients with poorly erupted posterior teeth. They are also ideal for the direct-indirect technique. The 6.5mm matrix is designed for larger cavities occlusogingivally. It augments a contoured matrix with a gingival apron and an adjustable continuation of the contour at the marginal ridge. The 7.5mm matrix is useful for patients with periodontal issues that require more coverage toward the gingival margin.

Technique Tip: Rings, Matrices, Wedges and WedgeGuards may be used with most popular matrix systems, including other sectionals. The Rings, Wedges and WedgeGuards may be used with most circumferential matrices. The benefits of the Ring, Wedge and WedgeGuard contours assist proper proximal contour.

6. Application of Matrix and Wedge

- Grasp the matrix with the pin tweezers using the hole located on the occlusal tab. If the ring was placed properly before tooth preparation, then the matrix should slide easily into place adjacent to the space to be restored from the occlusal approach. Ensure matrix is apical to preparation margin with no gap.

Technique Tip: If WedgeGuard wedge prevents matrix from fully seating gingivally, slide wedge out slightly, place matrix, then reposition wedge.

- Close the gingival margin by placing a contoured wedge or anatomic wooden wedge. Palodent® V3 wedges contain a placement tab with hole, which is securely grasped by the pin tweezers. This assures secure placement and minimizes the likelihood of dropping small accessories. **Technique Tip:** Palodent® V3 wedges' hollow design allows for easy placement of a second wedge from the opposing side (if needed) without dislodging the first. **Technique Tip:** Alternatively, the fork design allows placement of the Palodent® V3 retaining ring first if desired. The wedge can then be easily placed through the fork tine without dislodging the ring.

- Replace the Palodent® V3 retaining ring in the proximal plane, engaging the matrix and the wedge, thus stabilizing it against tooth structure. The forks of the tines should be sitting over the wedge.
- Lightly burnish the matrix against the adjacent tooth both buccally and lingually to further adapt it to form proper contours. **NOTE:** Over burnishing may distort matrix and compromise contours. The pin tweezer has a built-in ball burnisher for this purpose. The matrix placement tab may be bent over the adjacent marginal ridge to facilitate restoration marginal ridge contour. **Technique Tip:** When restoring an MOD preparation, two rings may be placed facing opposite directions, one mesially and one distally or both pointing in the same direction, usually mesially, one over the other placing the mesial matrix and ring first.
- Placement of Restorative
Place restorative material and establish proximal contact as per usual technique. Light cure and/or allow material to reach initial set. Gross contouring may be accomplished with matrix in place.
- Matrix Removal
Remove ring(s) using forceps, remove wedge(s) and sectional matrices using pin tweezers. Finish contouring and verify appropriate proximal contact. **Technique Tip:** It is advisable to "post-cure" light cured restorative materials following matrix removal from all surfaces, including proximo-facial and lingual line angles.

4. HYGIENE

Cross-contamination

- Do not reuse single use products. Dispose of in accordance with local regulations. During incineration, Wedges and WedgeGuards may produce irritating and highly toxic gasses.
- Reprocess reusable products as described below.

Instructions for Disinfecting/Sterilizing the Palodent® V3 placement forceps, pin tweezers and retaining rings

Warnings	<ul style="list-style-type: none"> High level disinfection by itself is not appropriate for the Palodent® V3 placement forceps, pin tweezers and retaining rings. It is not recommended to submerge Palodent® V3 placement forceps, pin tweezers and retaining rings into disinfection solutions. Do not allow the device to exceed 134°C. The Palodent® V3 sectional matrix band, Wedges and WedgeGuards are intended for single use only. Do not attempt to reprocess and/or re-use. Properly dispose used Wedges, WedgeGuards and matrices with other contaminated sharps. Properly dispose of the used and/or contaminated matrices Wedges or WedgeGuards in accordance with local regulations. During incineration, Wedges and WedgeGuards may produce irritating and highly toxic gasses. Ensure that phosphoric acid etchants are removed from the retaining ring tines prior to steam sterilization.
Limitations on	<ul style="list-style-type: none"> The Palodent® V3 placement forceps, pin tweezers and retaining rings can be subjected to the cleaning and Reprocessing automated disinfection up to 500 times, and sterilization process up to 1000 times. Repeated processing has minimal effect on the Palodent® V3 placement forceps, pin tweezers and retaining rings. Cold liquid disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.
Point of Use	<ul style="list-style-type: none"> Remove excess soil with disposable cloth/paper wipe. Most ultrasonic cleaning solutions will not affect the Palodent® V3 rings, although chemicals in some cleaners may cause tine discoloration. The discoloration is only cosmetic, but if that is a concern put the rings in a beaker with a water/detergent solution and place that in the ultrasonic cleaner to prevent contact with chemicals. It is recommended that the Palodent® V3 placement forceps, pin tweezers and retaining rings are reprocessed as soon as is reasonably practical following use.
Containment & Transportation	<ul style="list-style-type: none"> No particular requirements.
Cleaning & Disinfection: Automated:	<ul style="list-style-type: none"> Use only a properly maintained, inspected, calibrated, and approved washer/disinfector (according to ISO 15883). Perform a cleaning program cycle (the following cycle is recommended): <ul style="list-style-type: none"> 1 minute pre-wash/rinse in cold tap water. 5 minutes wash with enzymatic detergent at 55°C. 6 minute enzymatic neutralization. 1 minute rinse in deionized water. Thermal disinfection with minimum A0 of 3000 (e.g. 5 minutes at 90°C). 30 minutes drying at 60°C. Follow manufacturer's recommendation for use of detergent and neutralizer, observing concentrations and contact times.
Cleaning & Disinfection: Manual: detergent.	<ul style="list-style-type: none"> As an alternative to automated cleaning and disinfection, the Palodent® V3 placement forceps, pin tweezers and retaining rings have to be cleaned by scrubbing with hot water and soap or <ul style="list-style-type: none"> Use a new clean pair of examination gloves After cleaning thoroughly wipe all device surfaces with a single-use cloth in combination with an alcohol-based, bactericidal, virucidal and fungicidal disinfection solution approved according to local regulations and use according to disinfectant solution manufacturer's Instructions for Use. Make sure disinfectant solution is compatible with cleaning solution. Pay special attention to device seams and insertions. Remove disinfectant solution residue with a cloth soaked with water. Discard gloves according to local regulations and disinfect hands with an appropriate bactericidal, virucidal, and fungicidal hand disinfectant solution according to local regulations. Use according to disinfectant solution manufacturer's Instructions for Use. Dry device with a lint-free single-use cloth.
Packaging Sterilization	<ul style="list-style-type: none"> Paper/plastic steam sterilization pouches may be used, but are not required. Use a steam autoclave. Do not exceed 134°C. Steam autoclaving – Full cycle: wrapped or unwrapped at 134°C, 2 bar/273°F, 29 psi for 3 minutes.
Drying dry before storage.	<ul style="list-style-type: none"> Use the drying cycle of the autoclave. Do not exceed 134°C. Allow the Palodent® V3 placement forceps, pin tweezers and retaining rings to fully <ul style="list-style-type: none"> Avoid leaving the Palodent® V3 placement forceps, pin tweezers and retaining rings in an autoclave overnight or for long periods, as this may cause staining from the autoclave itself.
Maintenance	<ul style="list-style-type: none"> If the Palodent® V3 placement forceps, pin tweezers or retaining rings are damaged, worn, or distorted they should be discarded. Before use, visually inspect the tines for damage, such as cracking or pitting. The Palodent® V3 forceps and pin tweezers require no additional maintenance and should not be oiled or greased. If a Palodent V3 retaining ring is deformed by over stretching during placement or removal, the ring should be discarded.
Storage	<ul style="list-style-type: none"> Store the Palodent® V3 placement forceps, pin tweezers and retaining rings at room temperature, away from moisture or excessive humidity. Palodent® V3 placement forceps, pin tweezers and retaining rings should be disinfected/sterilized and dried prior to storage. Sterilization of new matrices, Wedges and WedgeGuards pre-use (if desired) may be accomplished (one time) by steam autoclaving, 134°C, 3 minutes minimum prior to storage and/or use.
Manufacturer	<ul style="list-style-type: none"> Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact Contact your local Dentsply Sirona representative.

5. LOT NUMBER AND EXPECTED USEFUL LIFE

- The following numbers should be quoted in all correspondence:
 - Reorder Number
 - Lot number on the package
 - Copy of proof of purchase

Manufactured by
Dentsply Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
1-302-422-4511
www.dentsplysirona.com

Dentsply DeTrey GmbH
De-Trey-Strasse 1
78467 Konstanz
Germany
49-7531-583-0
www.dentsplysirona.com

Distributed by
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada

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© 2017 Dentsply Sirona Inc. Todos os direitos reservados. **559702v** (R 10/5/17)