EN PermaShade™ LC Light Cure Veneer Cement

Description: PermaShade™ LC Veneer Cement is a light-cured, moderately viscous, permanent resin cement. It is 66.5% filled by mass, is available in four shades, and has very low shrinkage and color change properties.

Indications for Use: PermaShade™ LC Veneer Cement can be used for permanent cementation of porcelain, zirconia, composite and other indirect veneers. Directions for Use:

1. Remove PermaShade™ LC from refrigerator and bring to room temperature before use.
2. Clean preparation, rinse and dry (Consepsis® Scrub with a STARBrush™ to clean, disinfect, and remove residual cement).
3. Verify prosthetic fit.
4. Prepare inside surface of prosthesis.

- a. Ceramic / porcelain: WARNING: Hydrofluoric acid gives off a vapor that irritates respiratory passageways and on a chronic basis can damage them. Hold vacuum tip adjacent to surface while etchants are dispensed. WARNING: Proper isolation technique must be used when using hydrofluoric acid. WARNING: Patient, clinician and assistant must wear protective eyewear when using this product.

i. Apply hydrofluoric acid (UltraDent Porcelain Etch with Inspira® Brush Tip) to inside surface of prosthesis for 90 seconds, rinse and dry.

ii. Verify the color and consistency of the product flow before applying by expressing a small amount onto a mixing pad or gauze. Discard any clear liquid that is expressed.

iii. Apply phosphoric acid (Ultra-Etch™ with Blue Micro® Tip) for 5 seconds to remove porcelain salts and debris formed by hydrofluoric acid etching. Rinse and dry.

iv. Apply silane (UltraDent Silane with Black Mini® Brush Tip) to inside surface of prosthesis for 60 seconds, dry and set prosthesis aside. Do not rinse.

b. Zirconia veneers:

i. Do not use hydrofluoric acid or silane.

ii. Air abrade internal surface and apply a zirconia primer as per manufacturer's instructions.

c. Composite veneers:

i. No preparation to the inside surface of composite veneer is needed.

5. Prepare tooth surface

a. Self-Etch Technique

i. Rinse prepared tooth and leave moist.

ii. OPTIONAL - Apply Consepsis® (where available) with Black Mini® Brush Tip to preparation with a scrubbing motion and air dry.

iii. Apply Peak® SE with Black Mini® Brush Tip - 20 seconds to tooth preparation using continuous agitation.

iv. Thin/dry for 3 seconds using 3-way syringe or by placing the high volume suction directly over preparation.

v. Proceed to "Bond" Step

b. Total Etch Technique

i. Apply phosphoric acid (Ultra-Etch™ with the Blue Micro® Tip) to all surfaces of tooth preparation for 20 seconds.

ii. Rinse thoroughly for 5 seconds and leave damp.

iii. OPTIONAL - Apply Consepsis® (where available) with Black Mini® Brush Tip to preparation with a scrubbing motion and air dry.

iv. Proceed to "Bond" Step

6. Bond:

a. Apply a puddle coat of Peak® Universal Bond with Inspira® Brush Tip - gently agitate for 10 seconds.

b. Thin/dry 10 seconds using high air pressure.

Preparation should appear shiny.

c. Light cure with VAL® LED Curing Light on Standard Power for 10 seconds. (For lights with output <600mW/cm² cure 20 seconds or for lights with output >600mW/cm² cure 10 seconds.)

7. Belestion

a. Bring PermaShade™ LC into the zone of contact with Inspira® Brush Tip - gently agitate for 10 seconds.

b. Thin/dry 10 seconds using high air pressure.

Preparation should appear shiny.

c. Light cure with VAL® LED Curing Light on Standard Power for 10 seconds. (For lights with output <600mW/cm² cure 20 seconds or for lights with output >600mW/cm² cure 10 seconds.)

8. Cement

a. Express a thin layer of the PermaShade™ LC onto the inside bonding surface of the veneer only, not on the prepared tooth.

b. Carefully position and seat the veneer using gentle pressure. Tack with curing light (VALO with Point Cure lens) for 2 seconds to avoid shifting.

c. Clean the excess cement from margins. Optional: Apply DeOx around margins to prevent oxygen inhibition layer. Light cure with VAL® LED Curing Light on Standard Power for 10 seconds. (For lights with output <600mW/cm² cure 20 seconds or for lights with output >600mW/cm² cure 10 seconds.)

9. Warnings and Precautions:

1. WARNING: Hydrofluoric acid gives off a vapor that irritates respiratory passageways and can cause damage if exposed on a chronic basis.

2. WARNING: Proper isolation technique must be used when using hydrofluoric acid.

3. WARNING: Patient, clinician and assistant must wear protective eyewear when using this product.

4. Do not avoid inhalation of irritating Porcelain Etch vapors, hold vacuum in close proximity of area being prepared.

5. Verify the color and consistency of the product flow before applying by expressing a small amount onto a mixing pad or gauze. Discard any clear liquid that is expressed.

6. Carefully read and understand all instructions before using.

7. Use only recommended tips.

8. If etchants come in contact with any soft tissue, immediately rinse area with any amount of water.

9. Verify flow prior to applying any product intraorally. If resistance is felt, replace tip and re-check. Use only recommended tips.

10. Do not allow Porcelain Etch or Silane to be ingested.

11. Re-cap syringes with the Luer Lock Cap. To avoid cross contamination, disinfect and wipe syringe with an intermediate level disinfectant between uses or use a disposable syringe cover.

12. Dispose of used tips and empty syringes properly. To avoid cross-contamination, tips should not be re-used.

13. Resins can be sensitizing. Wear gloves and avoid repeated contact of uncured dental resin with skin. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis or rash develops, consult a physician.

14. Do not pull back on syringe stem after expressing product; doing so produces an air bubble within the syringe which can cause leakage and deterioration of the product.

15. Keep product out of direct sun and heat.

16. Re-cap product immediately after use.

17. Keep out of reach of children.

18. Do not use after expiration date noted on syringe.

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