DE	Gebrauchsanweisung	2
GB	Instructions for use	14
FR	Mode d'emploi.....	26
ES	Instrucciones de uso	38
IT	Istruzioni per l'uso	50
PT	Instruções de uso	62
NL	Gebruiksaanwijzing	74
SE	Bruksanvisning	86
DK	Brugervejledning.....	98
NO	Bruksanvisning	110
FI	Käyttöohjeet	122

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Stand: 2020-08

Medical device – only for use by dental professionals.

Purpose

RetraXil is an astringent retraction paste for the temporary displacement of the marginal gingiva and drying of the gingival sulcus. The paste, which contains aluminium chloride, is inserted into the gingival sulcus. In addition to the displacement effect of the paste, the astringent aluminium chloride ensures a haemostatic effect. Depending on the clinical situation and the practitioner's individual preference, the paste can be used as an alternative or in combination with retraction cords.

Indications

For the temporary displacement of the marginal gingiva and drying of the gingival sulcus, e.g. for

- conventional or digital impressions,
- cementation of temporary and permanent restorations
- and the creation of Class II and V fillings.

Contraindications

The use of this product is contraindicated in cases of known or suspected allergies to components of the product.

Do not use on patients with periodontal disease, furcation exposure or exposed bone.

Potential side effects

Hypersensitivity reactions to the product or one of its components cannot be excluded. If suspected, information on the ingredients can be obtained from the manufacturer. Local tissue damage if the application time is too long.


Product features**Accessories to the packaging**

- A. Twisting aid (can be used multiple times after sterilisation)
- B. Application cannulas (for single use only – please discard after use)

Composition

Phyllosilicates, water, aluminium chloride hexahydrate, mineral-based oleogel, vinyl siloxane, Pigment Blue 15:3

Product application

Preparation

- Thoroughly rinse the sulcus with water and gently air dry before using the paste. Remove and save the cap (Fig. 1).
- Depending on the clinical situation, a retraction cord can be inserted into the sulcus before the application of the retraction paste.
- Screw on one of the application cannulas supplied (Fig. 2), tighten it using the twisting aid provided and check that it is securely attached (Fig. 3).
- Apply a small amount of the retraction paste to a mixing pad and discard.

Note: Do not bend the cannula tip.

Application

- Insert the application cannula into the sulcus and apply the retraction paste circularly around the tooth so that excess material wells out of the sulcus (Fig. 4).

Note: Do not use on prepared tooth surfaces – RetraXil can impair the adhesive bond of dentine-bonding agents.

- Leave the retraction paste in the sulcus for two minutes (Fig. 5).
 - No moisture should be allowed in during this time.

Note: A longer application time can lead to tissue damage.

Rinse / remove

- Once the application time has ended, remove the retraction paste from the sulcus using an air-water mixture and a suction device; take care not to leave any residue (Fig. 6).
- Check the retraction result and remove any residue from the patient's sulcus/mouth.

Note: If RetraXil is not rinsed out of the mouth sufficiently, it can negatively affect the quality of the impression.

Closure

- Remove the application cannula after application using the twisting aid provided (Fig. 7)

- Discard the used application cannula and reclose the syringe with the supplied cap (Fig. 8)
- Remove the detachable finger support from the syringe (Fig. 9)
- Clean the syringe with a suitable, commercially-available disinfectant (Fig. 10)
- The twisting aid supplied and the detachable finger support must be sterilised after use (Fig. 11). **For more information, see "Cleaning, disinfection and sterilisation".**
- After the sterilisation is complete, reattach the finger support to the syringe (Fig. 12).

Note: Close RetraXil carefully after each use.

General notes

Please carefully observe the following safety instructions as well as all processing instructions included in other sections. This product is only to be processed in accordance with the instructions for use. Any alteration of these instructions during use is at the discretion of and under the sole responsibility of the dental professional.

Visually inspect for damage before use. Damaged products must not be used. Do not reuse the syringes if there is evidence that the paste has been contaminated with body fluids.

Warning and safety instructions

- The product may irritate eyes. Avoid contact with eyes. Immediately and thoroughly rinse with plenty of water in case of contact with eyes and contact an ophthalmologist.
- The product may irritate skin. Avoid contact with skin. If the product comes into contact with the skin, thoroughly wash with soap and water immediately.
- Wear appropriate protective clothing, eye protection and gloves.
The patient should use eye protection.
- The retraction paste must not be swallowed by the patient; this may result in an intolerance.
- Immediately clean paste residues from metal surfaces that have come into contact with the retraction paste, as they can corrode.

Storage conditions

Use within ten weeks of opening the product. Do not store above 25°C (77°F). Do not use the material after the expiry date.

Disposal instructions

Recommendation: Dispose of the content(s) in accordance with official regulations. Do not dispose of contents or packaging that have not been emptied in the household waste or allow them to enter the sewage system.

European Waste Catalogue: 180106 Chemicals consisting of or containing dangerous substances.

Cleaning and disinfection

WARNING: Please observe the country-specific regulations on hygiene standards and on the cleaning/disinfection of medical devices in dental practices. Wear suitable gloves. The twisting aids are supplied unsterilised. The twisting aid must be sterilised prior to first use and cleaned, disinfected and sterilised after each use.

Processing restrictions

Oxidising disinfectant solutions are not suitable. Strong acids and strong bases can damage the plastic.

SYRINGE instructions

In treatment room

Unscrew the cannula with the help of the twisting aid and discard it. Put on the grey cap. Remove the finger support.

Retention and transport

For transport outside of treatment room use a closed container.

It is recommended that the medical devices be reprocessed as soon as possible, but no more than two hours after use.

Preparing for decontamination

Wear personal protective equipment (sturdy gloves, water-repellent protective gown, face mask or safety goggles and mask).

Cleaning and disinfection: by hand

Equipment: non-protein-binding VAH-listed, virucidal disinfectant with cleaning action based on alkylamine derivatives and quaternary compounds (e.g. FD 300 top wipes, Dürr Dental).

- 1) Take out a disinfectant wipe and wipe the surfaces of the medical device thoroughly.
- 2) Make sure that all surfaces and existing grooves and notches of the medical device can be reached during the cleaning process and are completely wetted.
- 3) Check for cleanliness; if contamination is still visible, repeat the above steps.
- 4) Use another disinfectant wipe to wipe the surfaces of the medical device thoroughly.
- 5) Make sure that all surfaces and existing grooves and notches of the medical device can be reached during the disinfection process and are completely wetted.
- 6) Allow the disinfectant to take effect for the time specified by the disinfectant manufacturer (e.g. FD 300 top wipes, 5 minutes).
- 7) Allow the medical device time to dry.

Only validated processes may be used for the cleaning/disinfection of medical devices.

TWISTING AID instructions

In treatment room

Unscrew the twisting aid from the syringe and discard the cannula. Remove the finger support.

Retention and transport

For transport outside of the treatment room use a closed container.

It is recommended that the medical devices be reprocessed as soon as possible, but no more than two hours after use.

Preparing for decontamination

Wear personal protective equipment (sturdy gloves, water-repellent protective gown, face mask or safety goggles and mask).

A) Cleaning and disinfection: by hand

Equipment: non-protein-binding VAH-listed, virucidal disinfectant with cleaning action based on quaternary compounds (e.g. ID 212, Dürr Dental), lint-free cloths (e.g. Braun Wipes Eco), plastic brush (e.g. Interlock), if possible fully demineralised water (fully demineralised water, according to KRINKO/BfArM-recommendation free from facultative pathogenic microorganisms), tub for cleaning agents and disinfectants.

- 1) Prepare the cleaning solution according to the manufacturer's instructions (e.g. ID 212, 2%).
- 2) Immerse the medical devices completely in the cleaning solution.
- 3) While immersed, brush the hard-to-reach areas of the immersed medical devices with a soft brush. Pay attention to the critical hard-to-reach areas and cavities where it is not possible to visually assess the cleaning effect.

- 4) Immerse in the cleaning solution according to the manufacturer's instructions.
- 5) Remove the medical devices from the cleaning solution and rinse each one thoroughly (e.g. 30 seconds) under running water.
- 6) Check for cleanliness; if contamination is still visible, repeat the above steps.
- 7) Prepare the disinfectant according to the manufacturer's instructions (e.g. ID 212, 2%).
- 8) Immerse the medical devices completely in the disinfectant.
- 9) Immerse in the disinfectant according to the manufacturer's instructions (e.g. 5 minutes).
- 10) Remove the medical devices from the disinfectant and rinse each one thoroughly (e.g. 30 seconds) with demineralised water.
- 11) Dry with a lint-free disposable cloth and/or with medical compressed air according to the European Pharmacopoeia.

Only validated processes may be used for the cleaning/disinfection of medical devices.

B) Cleaning and disinfection: Automated WD [washer disinfecter]

Equipment: cleaning and disinfection device according to DIN EN ISO 15883-1 and 2 with thermal program (temperature 90°C /194°F to 95°C /203°F),

Cleaning agent: mildly alkaline cleaning agent (e.g. neodisher® MediClean forte, Dr. Weigert).

- 1) Brush the hard-to-reach areas of the medical devices for one minute under running potable water with a soft brush (e.g. Interlock). Pay attention to the critical hard-to-reach areas and cavities where it is not possible to assess the cleaning effect visually.
- 2) Place all medical devices in a suitable mesh tray or place them on the load carrier so that all internal and external surfaces of the instruments are cleaned and disinfected.

Programme step	Water	Amounts required	Time	Temperature
Pre-wash	CW (cold water)		5 min	
Cleaning solution dosing		Acc. to the manufacturer's specifications		Acc. to the manufacturer's specifications
Cleaning	DW (fully demineralised water)		10 min	55°C (131°F)
Rinsing	DW (fully demineralised water)		2 min	
Disinfection	DW (fully demineralised water)		3 min	A0 value > 3000 ¹ (e.g. 90°C / 194°F, 5 min)
Drying			15 min	Up to 120°C (248°F)

¹ The authorities may issue other parameters for disinfection performance in their area of responsibility.

- 3) Remove all medical devices at the end of the programme.
- 4) Check that the load is dry and, if necessary, dry with medical compressed air in accordance with the European Pharmacopoeia or with a lint-free cloth.

- 5) After removal from the washer disinfectant, perform a visual check for cleanliness. If contamination is still visible, clean the medical devices by hand. Automated processing of the cleaned medical devices must then follow.

Maintenance, inspection and testing

All instruments are to be checked visually for cleanliness, integrity and functionality, if necessary, using an illuminated magnifying glass (3-6 dioptre).

All medical devices must be checked for damage and wear. Damaged products must not be reused and must be disposed of.

C) Sterilisation

Device: Steriliser according to DIN EN 285 or small steam steriliser according to DIN EN 13060, Type B process.

Packaging:

- Equipment: Foil-paper packaging (e.g. steriCLIN, Item no. 3FKFS230110 and 3FKFS230112), sealing device (e.g. HAWO, Type 880 DC-V)
A suitable procedure (sterile barrier system) must be used for packaging the medical devices. Packaging according to DIN EN ISO 11607 (individual) or DIN 58953-9 (sets).
- Individual: A sterile barrier system (e.g. foil-paper packaging) according to DIN EN ISO 11607 should be used, one that is intended by the manufacturer for steam sterilisation. The packaging must be large enough so that the sealing seam is not under tension. (During the validation, double packaging is carried out in foil-paper packaging).
- Note: After the heat sealing process, the sealing seam must be checked visually for possible defects. In the event of errors, the packaging must be opened and the instrument repackaged and sealed.

- Sets: Sort medical devices into the appropriate sieves so that they cannot be mechanically damaged during transport. Sterile barrier systems in accordance with EN ISO 11607 (e.g. sterilisation fleece) (the instruments can be packed individually or as a double pack) or sterile goods containers in accordance with DIN 58953-9 are suitable for packaging the sieves. The weight of the sterile goods container must not exceed 10 kg.

Procedure: Steam sterilisation with fractional pre-vacuum, 134°C (273°F), holding time at least 3 min or 132°C (270°F) at least 4 min. Longer holding times are possible. In Germany, a holding time based on the KRINKO/BfArM-recommendation of at least 5 minutes is recommended.

- 1) Place the packaged medical devices in the sterilisation chamber.
- 2) Start the programme.
- 3) Remove the products at the end of the programme and allow them to cool.

Check the packaging for any damage or moisture penetration. Compromised packaging is to be assessed as non-sterile. The instruments must be repackaged and sterilised.

Storage: Storage (protect from dust and moisture) and storage duration at the user's premises according to the specifications.

Only validated processes may be used for the sterilisation of medical devices.

Please quote the batch number and article number in all correspondence about the product.

Article number **REF** and batch number **LOT**: See note on syringe/packaging.

According to the Medical Device Regulations, users/patients are obliged to report serious events with a medical device to the manufacturer and the competent authority in the country in which they occurred.

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Kulzer GmbH
Kulzer, LLC