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
WARNING AND PRECAUTIONS
DO NOT USE IF PACKAGE IS OPENED OR DAMAGED OR IF EXPIRATION DATE HAS BEEN EXCEEDED.
Do not overfill defects.
The maximum amount per site of implantation is 10cm³.
The device should be sealed to prevent motion and migration, used in areas where the graft can be adequately contained.
Do not leave defect open (e.g. muscle well sutured and use of appropriate membrane).
As with any surgical procedure, care should be exercised in handling patients with preexisting conditions that may affect the success of the surgical procedure.
This device does not adhere after implantation until healing. This device is a void filler used in cavity with respected cortices or associated with an osteosynthesis device.
Do not compromise blood supply to the defect area.
Maximum contact between the product and the recipient bone must be established.
This device has to be in direct contact of at least 3 surfaces.
Adequate post-operative immobilization is necessary to prevent movement that may lead to soft tissue ingrowth.
The Synthetic Bone Substitute should be manipulated with precaution in order to preserve the porous structure. It must not be out of contact before implantation.
The delivery anology is made of glass. A care must therefore be taken: the syringe cannot be used for other ancillary purposes.
The Synthetic Bone Substitute is intended for use by physician familiar with bone void filling and rigid fixation techniques.
The Synthetic Bone Substitute is to be handled or implanted by trained qualified physicians having read this instruction for use.
Caution: Federal law in USA and Canada restricts this device to sale by or on the order of a physician.

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SYNTHETIC BONE SUBSTITUTE
Maxillofacial/Dental application - Instructions for use

DESCRIPTION
The Synthetic Bone Substitute is a bone void filler. This product is a microporous and macroporous biphasic calcium phosphate ceramic consisting of 60% hydroxyapatite (HA) and 40% beta-Tricalcium Phosphate (β-TCP).
The Synthetic Bone Substitute is presented in a porous form required for the biological exchanges particularly for bone ingrowth.
The Synthetic Bone Substitute may be used with physiological saline, patient's own serum, whole blood, or bone marrow aspirate (BMA).
The Synthetic Bone Substitute is available in various shapes packaged in vials or syringes to be easily adapted to several bone defect designs. This product is provided sterile for single patient use.
The Synthetic Bone Substitute is progressively replaced by new bone according to the remodeling process of osteoblasts.

INTENDED USE
The Synthetic Bone Substitute is intended for use to reconstruct bony voids or bone gaps in particular in maxillofacial applications:

- Comblement de déficits parodontaux/intra-osseux
- Periodontal/intrabony defects filling
- Ridge augmentation
- Elevation of the sinus
- Relief of the cavities clefts
- Cystic cavities filling

The surface of the Synthetic Bone Substitute are the filling of bone defects and the bony ingrowth from local osseous tissue onto the surface of the product (osteointegration process).
The total resorption of the product depends on numerous factors such as its size and volume, the location of the defect, the surgical technique, the health status of the patient. The bone reconstruction time is always inferior to the total resorption time.

SURGICAL INDICATIONS
The Synthetic Bone Substitute is indicated for bony voids surgically created or bone gaps osseous defects created from traumatic injury to the bone.
The Synthetic Bone Substitute can be used with autograft as a bone void filler extendender.
The target population is adults, excluding pregnant women. The Synthetic Bone Substitute has not been tested on pregnant women and neither for use in pediatric population.

CONTRAINDICATIONS
The Synthetic Bone Substitute is contraindicated where the device is intended as structural support in the skeletal system.
Conditions representing contraindications include:

- Osteomyelitis
- Implantation in necrotic surgical sites
- Degenerative bone disease
- Long-term steroidal therapy or treatment acting on the calcium or phosphorus metabolism.

ADVERSE EFFECTS

- Post-operative discomfort
- Hematoma
- Infiltration
- Risque d'allergie
- Infection
- Foreign body reaction
- Pain
- Fever
- Inflammatory reaction

- Osmotic shock
- Adverse effects on blood coagulation and blood components
- Thrombosis
- Necrosis
- Non-union/pseudarthrosis due to displacement of product or fibrous encapsulation

Some of adverse effects may be treated by drugs or require reintervention.


STERILISATION AND CONSERVATION
The Synthetic Bone Substitute is provided sterile and is sterile unless the inner packaging has been opened or damaged. The Synthetic Bone Substitute must not be resterilized. The sterilization process of this product using radiation was validated during its conception. We do not guarantee the sterility of the product if a new sterilization is performed.
The Synthetic Bone Substitute is for single patient use only and should never be reused. The non-respect of the single use can conduct to septic issues.
The Synthetic Bone Substitute must be stored in its original packaging and kept at a temperature between 15°C and 30°C (59°F and 86°F).
Any product or container which has been in contact with the patient or is contaminated, should be treated as clinical waste and disposed of in accordance with the local regulations for the disposal of clinical waste.

WARNING AND PRECAUTIONS
DO NOT USE IF PACKAGE IS OPENED OR DAMAGED OR IF EXPIRATION DATE HAS BEEN EXCEEDED.
Do not overfill defects.
The maximum amount per site of the implantation is 10cm³.
The device should be sealed to prevent motion and migration, used in areas where the graft can be adequately contained.
Do not leave defect open (e.g. muscle well sutured and use of appropriate membrane).
As with any surgical procedure, care should be exercised in handling patients with preexisting conditions that may affect the success of the surgical procedure.
This device does not adhere after implantation until healing. This device is a void filler used in cavity with respected cortices or associated with an osteosynthesis device.
Do not compromise blood supply to the defect area.
Maximum contact between the product and the recipient bone must be established.
This device has to be in direct contact of at least 3 surfaces.
Adequate post-operative immobilization is necessary to prevent movement that may lead to soft tissue ingrowth.
The Synthetic Bone Substitute should be manipulated with precaution in order to preserve the porous structure. It must be out of contact before implantation.
The delivery anology is made of glass. A care must therefore be taken: the syringe cannot be used for other ancillary purposes.
The Synthetic Bone Substitute is intended for use by physician familiar with bone void filling and rigid fixation techniques.
The Synthetic Bone Substitute is to be handled or implanted by trained qualified physicians having read this instruction for use.
Caution: Federal law in USA and Canada restricts this device to sale by or on the order of a physician.

PREPARATION
Before use, the product must be hydrated without excess with physiological saline to prevent osmotic damage. If desired, the product may be mixed with patient's whole blood following hydration with a solution saline.
Tubes
Remove the vial from sterile package using sterile technique. Dispense desired quantity of granules into a sterile dish. Granules should be wetted with physiological saline to facilitate delivery to the surgical site using conventional instruments such as curettes or plastic carriers.
Syringe
Remove the syringe from sterile packaging using sterile technique. Withdraw the plunger to facilitate delivery of granules which will be wetted. Draw physiological saline or use patient's own blood above level of granules. Express excess liquid by depressing plunger slightly, without undue force.

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