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

Preparation according to EN ISO 17664

Parts which belong to the Image Plate- and Film Holder System Set COLOR and require preparation:
Aiming rings without lead windows, Bite blocks, Indicator arms

Aiming rings with lead windows only requires cleaning with a disinfectant wipe.

[MD] Medical device. Parts are medical instruments which come into contact with mucous membranes and skin or tissue which has been affected by illness. Remember to observe all current local and national regulations concerning the cleaning, disinfection, sterilisation and storage of the parts! In the Federal Republic of Germany, for example, a risk evaluation must be carried out and medical products are to be categorised according to the Recommendations of the Commission on Hospital Hygiene and the Prevention of Infections. The medical device has been categorised as: „semikritisch B with higher levels of preparation“.

Transport/Storage

After use parts must be placed in a suitable closed container which is filled with an instrument disinfectant solution collected and transported in a securely sealed vessel to the site of cleaning.

Preparation for cleaning

No special requirements are necessary. Disassemble the parts. Where necessary, remove larger dirt particles etc. using a disinfectant disposable wipe.

Cleaning and disinfection, automatic

[T] Thermally disinfectable
Use cleaning and disinfection devices strictly according to manufacturer's instructions, e. g. Typ Model Miele (Vario TD)

1. Place the parts into the unit in such a way that the interior surfaces will be rinsed thoroughly and the water can drain off.
2. Set up the cleaning and disinfection cycles according to the manufacturer's instructions.
3. After the cleaning and disinfection cycles have finished check the parts carefully f. signs of residual dirt and repeat the cycles if necessary.

Cleaning and disinfection, manual

Cleaning takes place in a wet chemical process using a VAH/DGHM listed instrument disinfection agent:

- Do not use ultrasonic cleaning or chemical autoclave.
- Do not use glutaraldehyde with a phenol base.
- Place the parts fully submerged in an instrument-disinfectant and cleaning agent. Observe manufacturer's instructions carefully.
- After the cleaning and disinfection cycles have finished check the parts carefully for signs of residual dirt and repeat the cycles if necessary.
- Rinse the parts thoroughly with water.
- Drying should be carried out using hygienically clean, disposable wipes, or under the application of compressed air.

Check and test of function

Frequent application and cleaning has little effect on the parts. The end of the product's life cycle is determined by the amount of wear and tear or damage during usage. Check the parts visually and where necessary, damaged parts must be disposed of. Before disposal the parts must be disinfected and sterilised.

Packaging

Parts can be sterilised either unpacked or packed in a container or in sterile product packaging. Use a suitable sterile product packaging!

Steam Sterilisation

! WARNUNG

The wrong sterilisation can reduce effectiveness and can lead to damage to the product

- Only steam sterilisation is permissible. • Do not exceed 138 °C. • Observe the process parameters.
- Do not use any other process. • Place filter paper on sterilisation tray before use.
- Always place objects for sterilisation together with instruments in order to achieve a balanced temperature distribution during sterilisation.
- Avoid metal contact between any of the objects for sterilisation and any instruments or with the tray edge itself.
- Place the objects to be sterilised in such a way that they do not come too close to the heating elements and do not come into contact with the internal walls and place on the middle tray.
- Only sterilise the product in a disassembled state.

[T] Minimum 3 min (max. 20 min) at 135 °C -1°/+2°, sterilisation process Model B (with pre-vacuum)

- Sterilised implements must be sterilised in a small steam steriliser before use according to European guidelines EN 13060 in sterilisation cycle B and subsequently dried.
- Observe the manufacturer's instructions concerning usage of the steam sterilisation unit and the correct placement of objects to be sterilised.
- Observe all national and local guidelines and regulations concerning sterilisation of medical products.

Wipe- & spray disinfectable

[W] Instruments not suitable for disinfection in washer disinfectors or for steam sterilization (Aiming rings RWT window x-ray)

Storage

Store the parts so they are protected against any possible contamination.

Important information

The instructions described above have been validated by Kentzler-Kaschner Dental as suitable and are recommended for the preparation of the parts before and after use. The surgery manager is responsible for ensuring that the correct and safe preparation of the instruments is carried out using suitable appliances, materials and trained personnel in order to produce the desired results. Validation and routine monitoring of the preparation procedures are necessary. Any and every deviation from the procedures described above by preparation staff is solely the responsibility of the personnel in terms of effectivity and the possible negative repercussions.