de	Gebrauchsanweisung	2 – 3
en	Directions for Use.....	4 – 5
fr	Mode d'emploi.....	6 – 7
it	Manuale d'uso.....	8 – 9
es	Manual de usuario	10 – 11
nl	Gebruiksaanwijzing	12 – 13
cs	Návod k použití	14 – 15
pl	Instrukcja użytkowania.....	16 – 17
ro	Instrucțiuni de utilizare	18 – 19
tr	Kullanım Talimatı	20 – 21
el	Οδηγίες χρήσης.....	22 – 23
pt	Manual de instruções.....	24 – 25
ru	Инструкция по применению	26 – 27
zh	使用说明	28 – 29
ko	사용 설명서.....	30 – 31
sr	Uputstvo za upotrebu	32 – 33
sv	Bruksanvisning	34 – 35
da	Brugsanvisning	36 – 37
hu	Használati útmutató.....	38 – 39
id	Petunjuk Penggunaan.....	40 – 41
ja	取扱説明書	42 – 43
uk	Інструкція з використання	44 – 45
ar	لامعة سرال لى لى	46 – 47

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VW000394 Rev. 1/11.01.2017

For dental use only, Rx only

RECIPROC® *blue* instruments for root canal preparation

1. Composition

The working part of the RECIPROC® *blue* instruments is made of nickel-titanium alloy.

RECIPROC® *blue* undergoes a thermal treatment process for more flexibility and a greater cyclic fatigue resistance compared to traditional nickel-titanium. Due to this processing, RECIPROC® *blue* instruments may appear slightly curved. This is not a manufacturing defect and does not weaken the instrument performance. It is not necessary to straighten the instrument prior to use in the root canal.

2. Contents

RECIPROC® *blue* blister package, 4 or 6 sterile instruments each

3. Indications for Use

- Preparation of the root canal system.
- Retreatment of the root canal system (removal of gutta-percha filling material and carrier-based obturators).

Due to the thermal-treatment process RECIPROC® *blue* instruments can be slightly pre-curved (max. 10°) to have an easier access to the orifice of canals which are not easily accessible.

RECIPROC® *blue* may only be used in reciprocation with a designated drive system with original VDW RECIPROC® settings, e.g. VDW.CONNECT Drive® or VDW.GOLD® RECIPROC®. The use of reciprocating drive systems without original VDW RECIPROC® settings can lead to misuse.

Due to its specific design, the instrument cuts counterclockwise in reciprocation.

4. Contraindications

If in rare cases, the hand file used for the working length determination (after the RECIPROC® *blue* instrument has reached 2/3 of the working length) needs to be pre-curved in order to reach the working length, it is necessary to create a glide path up to ISO size 15. If the ISO size 15 hand instrument has to be pre-curved in order to reach full working length after the creation of a glide path, there is an abrupt apical curvature in the apical region. The use of RECIPROC® *blue* or possibly pre-curved RECIPROC® *blue* instruments is contra-indicated in this instance. In these cases the canal preparation has to be finished with hand files. This limitation also applies to continuous rotary instruments.

5. ⚠ Warnings

None known.



6. ▽ Precautions

- The instruments may only be used by qualified persons in dentists' offices or clinics.
- RECIPROC® *blue* instruments are sterile, single use instruments, designed for the preparation and the retreatment of no more than one molar and cannot be reused. The instrument cannot be re-sterilized.
- Examine the instrument during treatment for signs of wear such as untwisting, bending or damaged cutting edges. In these cases, instruments should be discarded.
- Clean the flutes from debris in an interim stand every three pecks and irrigate the canal according to the appropriate irrigation protocol.
- Only light pressure should be applied on RECIPROC® *blue* instruments. If resistance is encountered or if the instrument does not advance, clean the instrument, irrigate the canal according to the appropriate irrigation protocol and check the canal for possible blockage with a C-PILOT® file ISO size 10. Then continue with the treatment briefly using a lateral brushing motion. If it still advances with difficulty or is unable to advance, it is necessary to create a glide path with C-PILOT® files up to ISO size 15.
- Use a rubber dam.

7. Adverse Reactions

None known.

Symbols used

	Explanation
	Manufacturer
	Packaging unit
	Catalogue number
	Batch code
	Use-by date
	Sterilized using irradiation
	Do not use if packaging is damaged
	Consult instructions for use
	Do not re-use
	CE marking
	Warnings
	Precautions
	RECIPROC®
	Nickel-Titanium

8. Preparation with RECIPROC® blue

Step-by-Step Instructions:

- Create a straight line access to the root canal entrance
- Select the correct RECIPROC® blue instrument:

In most cases, the R25 will be suitable in size for the root canal treatment. Consult the pre-operative radiograph to see if the canal is likely to be considered as narrow, medium or wide.

If the canal is partially or completely invisible on the radiograph:

The canal is considered narrow = use an R25.

If the canal is completely visible on the radiograph:

1. Take an ISO size 30 hand instrument; insert it passively into the canal. If it reaches working length, the canal is considered large = use the R50.

2. If an ISO size 30 hand instrument does not go passively to working length, try an ISO size 20 hand instrument. If this goes passively to length, the canal is considered medium = use an R40.

3. If a size 20 hand instrument does not go passively to length, the canal is considered narrow = use an R25.

Passively means that the instrument goes directly to working length without any watch-winding and can be moved up-and-down without filing action.

The amplitude of the up-and-down movement should be 4-5 mm.

Preparation with the R25 – Getting started:

Prior to preparation with an R25, estimate the working length with the help of a pre-operative radiograph. Set the stopper at 2/3 of that length.

Preparation with the R40 and R50 – Getting started:

Before selecting the appropriate RECIPROC® blue instrument, check the working length electronically by using a C-PILOT® or K-file. Set the stopper at that length.

Further steps:

1. Place irrigant in the access cavity of the root canal.

2. Introduce the RECIPROC® blue instrument into the canal. Activate the motor when the instrument is at the root canal orifice.

3. Move the instrument in a slow in-and-out pecking motion. The amplitude of the in-and-out movements should not exceed 3 mm. Only very light pressure should be applied. The instrument will advance easily in the canal.

One in-and-out movement = one peck. If necessary, the instrument can be used in a brushing motion to enlarge the root canal entrance.

4. Clean the instrument in the interim stand after 3 pecks.

5. Irrigate the canal.

6. Make sure the canal is free with an ISO size 10 C-PILOT® file.

7. Repeat step 3-6 until approx. 2/3 of the working length has been reached (indicated by the silicone stopper).

8. When using an R25:

Electronic length determination of the working length.

When using a R40 or R50:

After preparation of 2/3 of the working length, the working length should be rechecked with an apex locator.

9. Afterwards continue with step 3-6 until full working length has been reached.

10. Determine the apical preparation size – apical gauging by inserting a hand instrument which is one ISO size larger than the RECIPROC® blue instrument that reached full working length. If the hand instrument can be inserted (light tug-back) to 0.7 mm short of full working length and advances no further, the apical region does not need further enlargement. If the hand instrument reaches full working length, preparation of the canal should be completed with a larger RECIPROC® blue instrument.

9. Retreatment with RECIPROC® blue

Retreatment of gutta-percha root canal fillings and carrier-based obturators using an R25 – Instructions step-by-step:

1. Remove the bulk of the gutta-percha in the coronal third of the canal with an appropriate instrument e.g. Gates Enlarger or with the help of an ultrasonic device such as the VDW.ULTRA®.

2. If necessary, soften the gutta-percha with a little solvent, such as eucalyptus oil.

3. Use the R25 as described until working length has been reached. If resistance is encountered, do not apply pressure. Remove the instrument from the canal, apply a little solvent and repeat the process. The carrier may be removed in one piece during the use of the RECIPROC® blue instrument; otherwise, it will be removed in small pieces with the gutta-percha.

4. Use a brushing motion against lateral walls to remove residual obturation material.

5. After reaching working length with the R25, use an R40 or R50 for an increased apical enlargement, if necessary.

10. Disinfection, cleaning and sterilization

RECIPROC® blue is intended for single-use only.

Disinfection, cleaning and sterilization procedures are not recommended.

Re-use of RECIPROC® blue can increase the risk of cross contamination and breakage.

Further languages can be found on www.reciproc.com

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