DE

Gebrauchsanweisung
RESOLON® (PA 6-6/6)
Steriles, nicht-resorbierbares
Nahtmaterial

EN

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= RN = Stumpfe Rundkörpernadel

= R = Rundkörpernadel

The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.

Polyamide 6-6/6

Sterile, nonabsorbable
suture material

DESCRIPTION

RESOLON® is a monofilament, nonabsorbable, surgical suture material made from polyamide 6-6/6, extruded from a copolymer of polyamide 6 and polyamide 6/6. RESOLON® is available dyed blue (copper phthalocyanine blue, Colour Index 74160). RESOLON® is available in various thicknesses and lengths and also as a needle-thread combination. RESOLON® satisfies the requirements of the European Pharmacopoeia for sterile nonabsorbable suture material (0324).

INDICATIONS

RESOLON® (Polyamide) is used for coaptation of soft tissues and/ or for ligatures in the surgical disciplines dermatologic surgery, general surgery and microsurgery.

MODE OF APPLICATION

The suture material to be used and its size should be chosen in accordance with the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound. RESOLON® may only be used by qualified medical personnel.

ACTION

RESOLON® causes only a slight inflammatory tissue reaction initially and is then slowly and progressively encapsulated by connective tissue. Although RESOLON® is not absorbed, gradual hydrolysis of the polyamide in vivo can eventually lead to a progressive loss of tensile strength.

CONTRAINDICATIONS

As RESOLON® eventually loses some of its tensile strength, it is not suitable for indications in which lasting stability of the suture material is required.

WARNINGS / PRECAUTIONS / INTERACTIONS

Care should be taken by known hypersensitivity against one of the constituents of the nickel-chromium containing austenitic alloy of needles.

As with all suture material, care must be taken to ensure that the thread is not damaged during handling. In particular, it must not be

kinked or crushed by surgical instruments such as needle holders. When tightening the suture always pull on the thread between the needle and the puncture channel. Do not pull the thread too strongly or over sharp objects. When stretching the thread avoid friction with the rubber glove, as this can damage the thread.

In order not to damage the needle, always grasp it 1/3 – 1/2 of the distance from the end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard needles in "Sharps" containers.

The risk of wound dehiscence varies with the location of the wound and the suture material used, therefore the user should be familiar with the surgical techniques in which RESOLON® is to be used.

Post-operative care shall ensure that the suture is not subjected to considerable tension. RESOLON® must not be reused, there is a risk of infection for the patient and/or user and the functionality of the product cannot be guaranteed.

ADVERSE REACTIONS

As with all foreign bodies, prolonged contact between RESOLON® and saline solutions can lead to the formation of calculi (urinary and biliary tracts). Inflammatory tissue reactions and transient irritation in the area of the wound may occur initially. The foreign-body effect can exacerbate pre-existing tissue inflammation initially.

STERILITY

RESOLON® is sterilised by ethylene oxide. Already opened or damaged packs must not be used! RESOLON® is for single use only and should not be resterilised. Discard open, unused sutures! This is also applicable for Multiple-suture packs/ multipacks.

STORAGE CONDITIONS

Store at not more than 25 °C protected from moisture and direct heat/ solar radiation.

Do not use after the expiry date!

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

REF = Reference number

LOT = Batch number

 = Use by year - month - day

 = Consult instructions for use

 = Do not reuse

 = Do not resterilise

 = Do not use if package is damaged

STERILEEO = Sterilised using ethylene oxide

 = Upper limit of temperature

 = Dyed, monofilament, nonabsorbable

PA = Polyamide

 = HIBC code

 = Manufacturer  = Loop suture

 = Ligature pack: sterile individual suture on a small roll


R_{only} = Prescription only (only for USA)

 = Package contents

 = Folding out the suture carrier to access the needle

CE 0197 = CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

Needle code

Needle code

Example of needle size and needle curvature

 = P = Spatula needle

 = PM = Spatula needle PREMIUM-cut

 = S = Reverse cutting needle

 = SI = Inner cutting needle

 = RT = Round-bodied cutting needle

 = SM = Reverse cutting needle PREMIUM-cut

 = SMI = Inner cutting needle PREMIUM-cut

 = RTM = Round-bodied cutting needle PREMIUM-cut

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