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GB

Instructions for use

Precision impression material, silicone based, condensation curing, kneadable consistency

Fig. 1

1. Tray preparation

Select a tray type suited for desired application. For optimal adhesion of impression material, we recommend use of **sili adhesive**. Apply a thin layer of adhesive film to tray, then let dry for 60 seconds (Fig. 1).

Fig. 2

2. Mixing

Fill measuring scoop evenly to top with **silaplast FUTUR** (Fig. 2). Smooth **silaplast FUTUR** onto mixing pad with a spatula and make a series of crosswise grooves (Fig. 3). Add recommended amount of catalyst (Fig. 4). Mix thoroughly, first using spatula and then by kneading between fingers until uniform colour is achieved (approx. 30 sec.).


Fig. 3

3. Cleaning and disinfecting

After removing tray from mouth, rinse impression with lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaral-dehyde for 15 min.

4. Pouring the model

We recommend waiting at least 30 minutes, but no more than 72 hours before pouring the model.

All class III and IV dental plasters and standard model-casting plastics can be used.

5. Electro plating

The impression can be electro plated using all commercially available copper or silver plating baths.

Fig. 4

6. Cleaning the tray

After setting, material can be dislodged with a blunt instrument. Remaining adhesive film can easily be removed with **sili Solvent** spray (FCHC-free; use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.

Important hints

- We recommend using **silasoft® Normal**, **silasoft® Special** or **silasoft® direct** for corrective impression after making a preliminary impression with **silaplast FUTUR**.
- Overdosing of catalyst and increased temperatures accelerate, underdosing and decreased temperatures retard the setting time.
- Catalyst bottles or tubes should be closed tightly immediately after use.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- Avoid contact of catalyst with eyes, may cause irritation. In case of accidental contact rinse eyes with plenty of water and seek medical care immediately.
- Do not leave any material residues in the mouth.
- Please follow the instructions of the safety data sheets!


Avoid contact with eyes and skin!
Wear suitable eye protection and gloves!

Further information:


Silicone based materials are proven a million times. On condition of a proper application, **undesired effects** are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test **before the application** of the material.

Only for dental use by qualified personnel.

silaplast FUTUR / cat[®] universal FUTUR

silaplast FUTUR / cat^f silaplast FUTUR

Indications for use:

- Preliminary impressions for two-stage putty and wash techniques
- Situation impressions
- Check bites

Technical data:

DIN EN ISO 4823 - Type 0

- **Mixing ratio:**
silaplast FUTUR / cat^f silaplast FUTUR
1 scoop = 13.5 g /
1 scoop = 0.50 g
silaplast FUTUR / cat[®] universal FUTUR
1 scoop = 13.5 g /
6 cm paste = 0.30 g
- **Colour code:**
silaplast FUTUR: white
cat[®] universal FUTUR: blue
cat^f silaplast FUTUR: light blue
- **Mixing time:**
approx. 30 sec.
- **Working time:**
approx. 1 min. 15 sec.*
- **Setting time:**
cat[®] universal FUTUR:
approx. 4 min. 45 sec.*
cat^f silaplast FUTUR:
approx. 3 min. 45 sec.*
- **Time in the mouth:**
cat[®] universal FUTUR:
approx. 3 min. 30 sec.
cat^f silaplast FUTUR:
approx. 2 min. 30 sec.
- **Strain in compression:**
approx. 3.8 %
- **Recovery from deformation:**
approx. 98.6 %
- **Linear dimensional change:** approx. 0.3 %
- **Application:**
At 23 °C ± 2 °C /
73 °F ± 4 °F,
50 ± 5 % rel. humidity
- **Storage:**

* from initiation of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. increased temperatures accelerate, decreased temperatures retard these times.

Ordering information:

- | | |
|---|--------------|
| silaplast FUTUR
Standard packing
Jar of 900 ml | 02001 |
| Refill pack
Bag of 900 ml | 02681 |
| 4-pack
Jars of 4 x 900 ml | 02002 |
| Clinic-packing
Bucket of 5400 ml | 02003 |

cat[®] universal FUTUR 02015
paste, tubes of 5 x 35 ml

cat^f silaplast FUTUR 02016
liquid, bottles of 5 x 50 ml

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 07243/510-0 · Fax: 07243/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personell).

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