<b>Venus</b> <sup>®</sup>			
<b>DE</b>	<b>Gebruiksaanwijzing</b>	<b>SE</b>	<b>Bruksanvisning</b>
<b>GB</b>	<b>Instructions for use</b>	<b>DK</b>	<b>Brugervejledning</b>
<b>FR</b>	<b>Mode d'emploi</b>	<b>NO</b>	<b>Bruksanvisning</b>
<b>ES</b>	<b>Instrucciones de uso</b>	<b>FI</b>	<b>Käyttöohjeet</b>
<b>IT</b>	<b>Istruzioni per l'uso</b>	<b>GR</b>	<b>Οδηγίες Χρήσης</b>
<b>PT</b>	<b>Instruções de uso</b>	<b>CZ</b>	<b>Návod k použití</b>
<b>NL</b>	<b>Gebruiksaanwijzing</b>	<b>HU</b>	<b>Használati utasítás</b>

**Venus**<sup>®</sup> Lichthärtendes Universalkomposit, radioopak

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## Venus<sup>®</sup> Instructions for use (GB)

Universal light curing composite, radioopaque

**Venus** is a lightcuring, radioopaque composite. It is indicated for adhesively bonded, tooth colored restorations. It is suitable for anterior and posterior application according to Microglass filler system, durable esthetic restorations and high gloss can be realized easily.

**Composition:**

**Venus** is based on a BIS-GMA matrix and contains 58.7% filler (by volume), which is:

- Barium Aluminium Fluoride glass (0.7 µm; max. < 2µm) • Highly dispersive Silicon Dioxide (0.04 µm)

**Clinical Indications:**

• Direct restoration of Class I-IV cavities (according to G.V. Black) • Direct composite veneers • Shape corrections (incls. diastemas, congenital defects) • Temporary splinting of teeth loosened by trauma or periodontal disease • Indirect restorations (inlays, veneers) • Milk teeth restorations • Stump build-up • Temporary repairs of porcelain, composite (in combination with an adequate repair-system)

**Clinical treatment protocol:**

Select an appropriate shade or shade combination from the **Venus** Shade guide prior to placing the rubber dam. This is important because teeth will significantly change colour when desiccated. If the tooth surface being treated is not to be instrumented by a diamond stone, carbide bur, drill or other abrasive instrument, the surface should be cleaned with pumice or a fluoridated polishing paste. A non-slicking plastic, Teflon coated or gold plated plastic instrument (such as the Kutzer Plasmacoat Instrument) is recommended for use in clinically placing, sculpting and shaping the restoration.
1. For prior application, the use of an Ivory rubber dam (or similar) is recommended to control contamination from moisture, blood and saliva.
2. Prepare the cavity according to generally accepted roles for composite application. Then clean and dry the cavity.
3. For ideal results, **Venus** restorative material to the cavity preparation, pre-treatment of the tooth structure with an adequate enamel/dentin adhesive is mandatory (e.g. IBOND<sup>®</sup> Self Etch, IBOND<sup>®</sup> Total Etch, GLUMA<sup>®</sup> Solid Bond of Kutzer). Once treated with the adhesive, avoid contamination with saliva or blood.

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4. Apply **Venus** Restorative material in thin layers (max. 2 mm) and adapt carefully to the tooth surface or cavity walls. In case of severe discoloration of the tooth surface or surface to be treated, it may be necessary to mask this by using a thin layer of one of the recommended **Venus** dentin shades first.

5. To avoid excess translucency and achieve an ideal esthetic result in the case of large Class IV cavities and anterior tooth reconstruction, it is necessary to utilize a **Venus** dentin shade in the deepest layers, which can then be covered by the enamel and incisal edge Halogen shades.

6. Polymerize **Venus** in thin layers (max. 2 mm) according to the table below using an acceptable QHT (Quartz-Halogen-Tungsten) curing light such as the Translux Curing Unit (>400 mW/cm²). During polymerization of the bonding resin, avoid contact with the surface of the composite. This layer allows the chemical bonding of the composite masses to each other. This layer should not be touched, contaminated or removed during the build-up phase of the restoration. It serves as chemical compensation for any subsequent composite layers that are added.

Recommended **Venus** curing times with Translux Curing Unit or other QHT light curing units:

A1, A2, A3, A3.5, B1, B2, C2, D2, T1, T2, SBl, HKA2.5	Curing time with QHT curing light	20 s
A4, B3, C3, C4, D3, D42, O42, O43, O43.5, O82, S8D, HKA5	Curing time with QHT curing light	40 s
All shades	Maximum layer thickness	2 mm

For Plasma Arc light curing units, the curing time can usually be reduced by 25%. For LED (e.g. Translux Curing unit or any unit with a similar spectrum) curing units, use the same curing times recommended with a QHT light. If the light source can be held in close approximation to the composite material, however, due to the relatively low light output of currently available LED curing units, it is recommended to increase curing time to at least 150%, if the light source will be positioned away from the surface of the material. Note: Curing lights vary greatly in intensity and spectral output. All curing times are minimum times. You should check with the light unit manufacturer and verify its efficacy non-clinically before using it on a patient.

After final build-up and polymerization is completed, **Venus** may immediately be buffed and polished. For gross finishing, use diamond burs. Fine finishing and shaping is best achieved with fine-grained diamond stones and multifluted carbide burs. Proximal areas may be finished with abrasive strips. A high final polish can be achieved with flexible, silicone rubber polishing points, cups, discs, or similar. Composite polishing pastes can be applied with polishing brushes and small buffing wheels if desired.

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