SHOFU VINTAGE LD PRESS

en Instructions for use

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en Lithium Disilicate Glass Ceramic

VINTAGE LD PRESS is Lithium disilicate glass ceramic press ingot to fabricate inlays, onlays, veneers, crowns and bridges

SHADES

T: T-1,T-2,T-3

MT: MT-W1, MT-W2, MT-W3, MT-A1, MT-A2, MT-A3, MT-A3.5, MT-A4, MT-B1, MT-B2, MT-B3, MT-B4, MT-C1, MT-C2, MT-C3, MT-C4, MT-D2, MT-D3, MT-D4

LO: LO-1, LO-2, LO-3, LO-4, LO-5
MO: MO-1, MO-2, MO-3

INDICATIONS

Fabrication of all ceramic crowns, inlays, onlays, veneers and bridges (up to 3 units excluding molar)

PRECAUTIONS

- If any inflammation or other allergic reactions occur on the operator while using this product, immediately discontinue use and seek medical advice.
- Wear protective glasses etc. while grinding and polishing this product to avoid damages to eyes.
- Avoid getting the debris of this product into the eyes. If they get into eyes while grinding or polishing, immediately flush the eyes with plenty of water and seek medical advice.
- Use local dust extractor, dust protective mask, etc. while grinding this product to avoid harmful influence of the dust on the human body.
- This product should not be used for the cases where residual abutment tooth is too small, or for the patients with bruxism or occlusal abnormality etc.
- This product is intended for use by dental professionals only.

DIRECTIONS FOR USE

- Wax up

Wax restoration to full contour in a conventional manner.

 - Maintain a minimum of 0.8 mm wax thickness over the copings except for the cervical area, where at least 0.3 mm wax thickness should be ensured.
 - Select investment ring size according to the weight of wax pattern. Prepare appropriate number of ingots.

Weight of the wax pattern	Number of ingots	Ring size
0.4 g or less	1	For mixing 100 g investment material
1.2 g or less	2	For mixing 200 g investment material

- Spruing

Attach wax sprues (3.3 mm in diameter) to the broadest place on each waxed up restoration in such a way that insures unrestricted flow of the pressable ceramic to all areas of the restoration(s).

 - Always keep a minimum distance of 5 mm between wax patterns. Maintain a distance of 10 mm or more to the outer wall and position all restorations in the same horizontal plane.
 - Connect one sprue to each individual pontic for bridges.
- Attach/Wax the sprue and position

Sprue should be 3-5 mm long. Attach/Wax the sprue to the edge of the ring base at an angle of between 45° and 60°.
- Investment

The investment should be processed according to the investment material manufacturer's instructions for use (Ceravety Press and Cast is recommended). Avoid entrapping air bubbles. Leave undisturbed until it completely hardens.

Remove the ring and ring base with a turning movement and check the 90° angle of the investment ring.

en Manufacturer
SHOFU INC.
11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan
www.shofu.com

SHOFU DENTAL CORPORATION
1225 Stone Drive, San Marcos, California 92078-4059, USA

SHOFU DENTAL ASIA-PACIFIC PTE.LTD.
10 Science Park Road, #03-12, The Alpha, Science Park II, Singapore 117684

EC **REPL**
SHOFU DENTAL GmbH
An der Pöntz 70, 40885 Ratingen, Germany

CE 0123

en Lithium Disilicate Glass Ceramic

READ CAREFULLY BEFORE USE

5.Burnout
Place the investment ring in the burnout furnace together with the disposable plunger for preheating.

6.Pressing
Take the investment ring and disposable plunger out of the burnout furnace. Insert the ceramic ingots and press disposable plunger into the ring, then center the hot investment ring on the pressing platform of the press furnace. Follow the pressing procedure indicated below.

Ring size (g)	Entry Temp. (°C) <°F>	Incr. Temp. (°C/min) <°F/min>	Press Temp. (°C) <°F>	Hold Time (min)	Pressure (MPa) (Maximum)	Press time (min)
100	700	60	910 <1,670>	15	0.5 (Maximum)	3
200	<1,292>	<108>	925 <1,697>	25		

* Vacuum should start with the entry temperature and should be released after complete pressing.

- This product should not be preheated.
- Do not reuse the remaining pressed material.
- For pressing procedure, follow the press furnace manufacturer's instructions for use.
- Optimum press schedule is different depending on the type and model of the press furnace. Try to determine the optimal conditions at which every restoration is pressed completely.

- Divesting

Divest the restoration once the ring has cooled to room temperature. First, carefully cut the sprue using a diamond disc. Then, divest the ring with glass beads at a pressure of 0.4 MPa (4 bar) for rough divestment and of 0.2 MPa (2 bar) for fine divestment followed by Al₂O₃ blasting at a pressure of 0.2 MPa (2 bar) to remove the investment material completely.

8.Grinding and contouring
Contour and grind the restoration at a low speed with light pressure while avoid overheating the restoration.

- Staining and Layering
 - Staining technique

After contouring and adjusting occlusion, apply VINTAGE Art LF stains on the cleaned surface. Then, fire following the instructions for use of the VINTAGE Art LF stains.
 - Layering technique

Contour the area where additional porcelain built up is required. Build up VINTAGE LD porcelains on the cleaned surface and fire following the instructions for use of the porcelain.

 - Refer to combination chart

NOTE ON USE

- Follow the instructions for use of each dental material or instrument to be used in conjunction with this product.
- Do not use this product in conjunction with PFMs, metal frames or ceramic frames (alumina and zirconia).
- Do not use this product for any purpose other than specifically outlined in INDICATIONS in these instructions for use.
- Use this product within the expiration date indicated on the package and container.

(Example ☞ YYYY-MM-DD → Year-Month-Date of the expiration date)

STORAGE
Store at room temperature (1-30 °C/34-86 °F). Avoid high temperature and high humidity. Keep away from direct sunlight.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

Location of the symbols glossary:
https://www.shofu.com/en/symbols_glossary/

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