VITA AKZENT® Plus FLUOGLAZE LT SPRAY

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VITA AKZENT® Plus FLUOGLAZE LT SPRAY

Quickstart Working Instructions

CTE-independent, highly fluorescent glaze spray

VITA – perfect match.

VITA

VITA AKZENT® Plus FLUOGLAZE LT SPRAY – FACTS & FEATURES

VITA AKZENT Plus FLUOGLAZE LT SPRAY is:

- a CTE-independent, highly fluorescent glaze spray for all standard veneering ceramics, monolithic zirconia, silicate or feldspar ceramics
- suitable for achieving a natural fluorescence, especially for monolithic, all-ceramic restorations without auto-fluorescence
- suitable for sealing the ceramic surface

When using VITA AKZENT Plus FLUOGLAZE LT SPRAY, you will benefit from:

- simple enhanced fluorescence
- fast and reliable glazing
- highly esthetic results through natural fluorescence

SYSTEM COMPONENTS

VITA AKZENT® Plus FLUOGLAZE LT SPRAY

Using VITA AKZENT Plus FLUOGLAZE LT SPRAY

✓ VITA AKZENT Plus FLUOGLAZE LT SPRAY application:

- Clean the restoration
- Shake spray can well immediately before use
- A distance of approx. 10 cm between the spray head and the surface to be sprayed is required
- Spray from all sides to apply a uniform layer
- Place the restoration on a firing tray for stains firing and glaze firing, and fire according to the firing parameters

> Optimal layer thickness

> Not enough applied

> Too much applied

! Please note:

- If the gloss is insufficient, the process can be repeated.

✓ Firing parameters


	Predry. °C	→ min.	↗ °C/min.	approx. temp. °C	→ min.	VAC min.	↘ °C
VITA AKZENT PLUS GLAZE LT paste	400	6.00	50	740-850*	1.00	-	500*

*Firing temperature and/or long-term cooling depend on the veneering ceramic or substructure material used. Observe the manufacturer's instructions.


1 Creation of natural fluorescence

2 Example of a natural crown and monolithic zirconia front tooth crowns

natural crowns


monolithic
zirconia crowns

> no auto-
fluorescence

> fluorescence
too weak

> streaky due to
use of paste

with competitor's fluorescent

agent

> fluorescence
too high

with VITA AKZENT **Plus** FLUOGLAZE LT

**Perfectly matched and homogeneously
applied fluorescent and glaze agent**

> VITA AKZENT **Plus** FLUOGLAZE LT

VITA AKZENT® Plus – Product portfolio

VITA AKZENT Plus EFFECT STAINS

- highly chromatic and opaque stains
- to reproduce individual shades and to adapt the shade to the characteristics of a natural shade
- when applied thickly, they cover the underlying base shade completely

VITA AKZENT Plus BODY STAINS

- translucent finishing materials that have the effect of a thin color filter on the base material
- for color changes and color shifts towards reddish, yellowish, orange, brownish or gray-greenish

VITA AKZENT Plus CHROMA STAINS

- glazing stains for modifying the chromaticity of the base material within a shade group

VITA AKZENT Plus GLAZE

- glaze material for a brilliant, glasslike, homogeneous and dense surface after firing

VITA AKZENT Plus GLAZE LT

- LT = Low Temperature: low-melting glaze material for ceramics with firing temperatures below 800°C


VITA AKZENT Plus FLUOGLAZE LT

- CTE-independent, highly fluorescent glaze spray to achieve natural fluorescence

VITA AKZENT Plus Finishing Agent

- for a natural surface texture with a matte finish

VITA AKZENT® Plus – Available shades

TECHNICAL HOTLINE

Questions about VITA product solutions?

The **dental technician and technical consultant**

Ralf Mehlin (DT) can answer your questions here:

- ▶ **Phone +49(0) 7761 / 56 22 22**
- ▶ **Email info@vita-zahnfabrik.com**
- ▶ **8 a.m. to 5 p.m CET**

Additional international contact information can be found

- ▶ at **www.vita-zahnfabrik.com/contacts**

ONLINE SUPPORT

Questions about processing and product solutions?

You can find competent answers online in the form of **FAQs, videos and electronic tutorials.**

More detailed information can also be found in the comprehensive VITA AKZENT Plus Working Instructions.

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VITA AKZENT® Plus FLUOGLAZE LT SPRAY

Utilisation de VITA AKZENT Plus FLUOGLAZE LT SPRAY

Application de VITA AKZENT Plus FLUOGLAZE LT SPRAY

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VITA Zahnfabrik H. Rauter GmbH & Co.KG