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5. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide). Note: els composite does not contain TEGDMA or HEMA.

6. Interactions

Avoid underfilling materials which may hinder polymerisation owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

7. Processing stages

7.1. Tooth cleaning

Brush the tooth being treated and the teeth next to it with fluoride-free toothpaste. If necessary, clean interdental spaces using strips and dental floss.

7.2. Shade selection

SAREMCO shades relate to VITA standards.

7.3. Drying

Dry sufficiently, apply rubber dam.

7.4. Cavity preparation

Prepare the cavity as usual. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

7.5. Underfilling

In case of near-to-pulp preparations, it is advisable to protect the pulp using an underfilling (e.g. calcium hydroxide). Do not use zinc oxide eugenol!

7.6. Adhesion

The adhesive has to be conditioned and applied in correspondence to instruction for use of the product. SAREMCO recommends cmf adhesive system, els dubond and els unibond.

7.7. Restoration with els composite

Application of the syringe: After each usage, turn the syringe screw back by one turn and re-attach the screw cap.

Application of the tip: Remove tip cap. Push the tip into the corresponding recess of the injector and fix and click it into place. (SAREMCO recommends SAREMCO applicator). Apply material with constant pressure of the lever.

Gradually adapt els composite up to a maximum layer thickness of 2 mm. Light-cure every layer for 40 seconds. When using matrices, remove them and expose filling to light again from lingual and buccal. Note: els composite produces an extremely low shrinkage stress and shrinkage.

Warning: Protect light-curing products from strong sources of light.

The specified exposure time refers to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm² and a wavelength of 400 - 500 nm. They apply to a maximal layer thickness of 2 mm. The necessary exposure time may vary depending on the light source and its instructions for use. In case of doubt, check the lamp's light output and the necessary exposure time before operation *in vitro*. During polymerisation, an inhibition layer forms on the surface; this must not be touched or removed, if other composite layers are going to be applied.

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Instruction for use els composite

1. Product description

Light-curing, highly-filled, radio-opaque micro hybrid composite with extremely low shrinkage stress, for anterior and posterior restorations. DIN EN ISO 4049-compliant.

2. Composition

inorganic filler (barium glass and silica 74 % wt, 49 % by volume, particle size between 0.04 - 3.00 µm, median 0.7 µm), BisGMA, BisEMA, catalysts, inhibitors, pigments

3. Indication

- Restoration of class I, II, III, IV and V of front and side teeth
- Extraorally made inlays and onlays (direct and indirect), veneers
- Aesthetic corrections of interdental spaces, enamel hypoplasia, discolorations etc.
- Blockages, fractures

4. Contra-indication

Protect light-curing products from strong sources of light and heat: els composite was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze!

5. Batch number and expiry date

The batch number should be specified to identify products in case of inquiries. Products should no longer be used once the expiry date has elapsed.

6. Finishing, polishing

Prepare the filling with 40 µm and 12 µm diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers. The filling can be prepared and polished immediately after polymerisation.

7. Storage

Protect light-curing products from strong sources of light and heat: els composite was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze!

8. Finishing, Polishing

Prepare the filling with 40 µm and 12 µm diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers. The filling can be prepared and polished immediately after polymerisation.

9. Batch number and expiry date

The batch number should be specified to identify products in case of inquiries. Products should no longer be used once the expiry date has elapsed.

10. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. When using tips, make sure, that before usage the cap has been removed. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap and put on a new glove. In case of an allergic reaction, seek medical advice. When polishing or removing composites, it is recommended to always use a water cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

11. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, seek thorough with water. Consult an eye specialist.

12. Hygiene

Use application instruments for one patient only. Dispose products away from patients to avoid contamination.

13. Hygiene

Our liability is restricted to the quality of our products. In case of a product being of defective quality, only its value will be replaced. For damages, notably that occur due to non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental.

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