VOCO

Ionofil

GB	In accordance with ISO 29917/1993	CE ⁰⁴⁸²
DE	E	
FR	E	
ES	E	
PT	E	
IT	E	
GR	E	
NL	E	
DK	E	
FI	E	
NO	E	
SE	E	

Instructions for use

Product description:

Ionofil is a glass ionomer restorative material with chemical bonding on enamel and dentine.

Fields of application:

- For restoration of deciduous teeth (especially class I)
- Repair of V-shaped defects and enamel or root erosions in the cervical area incl. class V fillings
- Cavities caused by carious lesions
- Fillings of class III and small fillings of class I (incl. extended fissure sealing) as well as build-up fillings for crowns, linings
- Attachment of posts and orthodontic bands

Method of use:

Creation of retentions generally is not necessary in view of the good adhesion of **Ionofil** to the dentine. Undercuts, however, provide additional hold to the filling. Layer thickness of the fillings should not be less than 1 mm.

In case of carious defects, clean cavity accordingly (for instance with a 3 % hydrogen peroxide) and then rinse carefully with water and dry (air jet). Do not overdry dentine. In case of a remaining dentine thickness of more than 1 - 1.5 mm, an additional lining is not necessary. A pulp protection with calcium hydroxide is recommended if the dentine thickness is less than 1 mm.

Shake powder well before use.

Place powder and liquid on a glass plate or a suitable mixing pad and work the powder into the liquid a little at a time to begin with, then in larger quantities. Use 1 drop of liquid for 1 level measuring scoop of powder (mixing ratio: 4.75 g of powder : 1 g of liquid). When dispensing the liquid, hold the dropper bottle absolutely in vertical position! A strong plastic spatula should be used for mixing. Mixing time should not exceed 1 minute.

Working time from end of mixing procedure is approximately 2 - 3 minutes. The temperature should be 15 °C - 23 °C.

Setting time, after insertion into the cavity, is approximately 2.5 minutes.

During filling procedure ensure that working area is kept dry.

Use of a rubberdam is recommended.

Class V cavities should be at least 1 mm in depth. Avoid feather-shaped margins. In case of doubt, prepare in a conventional way.

After application of **Ionofil**, cover filling immediately with a strip or a matrix.

During setting phase **Ionofil** is susceptible to water contamination. Immediately upon removal of the strips or the matrix, the filling therefore has to be covered completely with a varnish (see applicable instructions for use).

After setting of **Ionofil**, no sooner than 4 - 6 minutes after application (depending on mixing temperature and mixing ratio), excess material may be removed with a sharp instrument and contouring can be carried out (flexible disks). Then cover again with varnish (see applicable instructions for use). The best results for final polishing can be obtained about 24 hours after insertion of the filling.

Precautions:

In order to avoid irritations in close proximity to the pulp, use pulp protection (less than 1 - 1.5 mm). Close powder and liquid thoroughly after use in order to prevent humidity from entering and thus reducing the high strength of the material.

Avoid contact with mucosa and eyes.

Ionofil contains Na-Ca fluoroasilicate, polyacrylic acid and parabens.

Ionofil should not be used in case of known hypersensitivities (allergies) to any of these ingredients. Do not use insulating agents for the dentine in order to maintain full adhesion (exception: pulp protection).

Storage:

Store **Ionofil** powder and liquid containers tightly closed, at room temperature (4 °C - 23 °C).

Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

DE	Entspricht ISO 29917/1993
E	
FR	E
ES	E
PT	E
IT	E
GR	E
NL	E
DK	E
FI	E
NO	E
SE	E

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VC 60 BB1110 E1 0114 99 © by VOCO
VOCO GmbH