PERform / Inkotherm 85

CE 0297

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performd.odt Stand 16.02.2021

PERform / Inkotherm 85 CE 0297
INKOVAC-SYSTEM
Acrylic Resin

Perform resin, is a liquid resin especially made for the PERform InkoVac-System. It is stable in colour and cadmium free. Mixing proportion: 10 parts powder 5,5 parts liquid (use special measuring cups). Pressure polymerisation 30 minutes, 3-6 bar in 45 °C temperate warm water. For model casts dentures, rebases and repairs. Mixing proportion: 10 parts powder, 10 parts liquid, pressure polymerisation 15 minutes. 2-6 bar in 55 °C warm water. Perform is very economical because of its individual use. Auto-curing resin for dentures in the PERform - INKOVAC - System and the model casting technique, long handling time.

Technical dates: DIN 13907

Fracture strength: 70-75 N/mm²
Modules of rupture: 2100-2200 N/mm²
Typ II Klasse 1 gem. ISO 1567,
gem. Medical device Class. IIa CE
Residual monomer 24 hours 2-3%,
after 7 days below 1%

For Use in the PERform-InkoVac-System read the separate instruction manual of the system .

Perform System

Prepare the Perform flask for pouring as instructed in the Perform-InkoVac-System. Remove the wax up from the model by using hot water and insulate the warm model with Perform separating medium for acrylic by using vacuum.
Mix the Perform resin in the ratio 10 : 5.5 of powder and liquid (use Perform measuring cups), avoiding the formation of bubbles. The large marks on the measuring cups indicate the respective amounts required for a prosthesis of average size.

Note: The powder should be stored in a refrigerator. Mixing can be effected on a small vibrator, when air bubbles will rise. The mixed plastic material should be poured quickly into the sprue, via the funnel, until the plastic rises in the second opening. meanwhile the pouring of resin , the flask must be installed on the flask connector by using vacuum. After approximately 2 - 3 minutes, the vacuum is released, and the and the flask can be removed from the flask connector.

The vacuum membrane is checked to see whether

the vacuum in the flask is maintained. Place the flask into the polymerisation unit that the flask is totally covered with water. The funnel have to be out of the water level. Fill the pouring funnels with cold water or ice (use perform icecube former) that the pouring channel will be cooled (delay of polymerisation to compensate the shrinkage). Close the polymerisation unit and polymerise for 30 minutes with 45 °C and 5.5 bar pressure.

Model casting prosthesis, repairs and relinings use the material as follow.

Indication:

For completion of model casting prosthesis, for relinings, borders and repairs of dentures.
Material: Cold polymerising material based on methylmethacrylate. Available powder and liquid. Colour stability due to a catalyst system, without tertiary amines.

Characteristics:

Easy processing, extremely long plastic and pouring phase for modelling and at room temperature delayed start of polymerising. Perform can be applied economically. The high mechanical properties correspond to those of a good hot polymerising material. Natural shading, safe and quick curing.

Preparation:

For model casting dentures fix the teeth, which are standing in wax, by a plaster or silicone precast. The precasts, be carefully cleaned by scalding off the wax remainders. The plaster model must be well watered, so that during pressure polymerising in the Hedent PPA, no air is being pressed from the plaster into the resin which would form bubbles.

Isolating: Apply a thin layer of Inkolant-Alginate-Insulation to the model and plaster precasts.

Mixing ratio: Mix 10g Powder : 7ml liquid.
Free dosage is possible.

Working time:

Setting time: 1/2 minute. Flowing phase: up to 3rd minute, thereafter rigid.
Plastic modelling phase up to 20th minute. Start of polymerisation at least after 25th minute. Polymerising can also be done immediately after the flowing phase. Working times refers to a material and room temperature of 23 °C. Higher temperatures shorten, lower temperatures extend the indicated working time. The working time of Perform is extremely long and does not dry out.

Mixing:

Place liquid into mixing bowl and quickly stir in the powder until a creamy consistency, capable of dripping, is reached.

Processing:

The mixed resin is to be filled into the precasts during the flowing phase. At the beginning of the plastic phase the resin is rigid, does not longer flow out of the precast and can be modelled.

Polymerisation:

After a minimum of 25 minutes, polymerisation is done in a pressure polymerising unit - Hedent PPA - under a pressure of 3 - 6 bar and a water temperature of 45 °C.

Important:

The alginate isolation must be dried well, in order to avoid white discolouring at the basis areas. Application of activator is not necessary with Perform .

Packages:

500 g Powder - 250 ml Liquid
1000 g Powder - 500 ml Liquid
5000 g Powder - 1000 ml Liquid

Important

Use only original Hedent PERform powder and liquid. Proceed only with that instruction. In case of different handling the polymerisation cannot be sufficient and the patient could notice irritation of skin and taste. Also the mechanical values can be lost. To trace back, always the LOT number has to be noted for each work. (European law of medical devices)

Caution Perform liquid monomer contains: Methylmethacrylate. Highly flammable. Harmful by inhalation and in contact with skin and if swallowed. Keep bottle always closed. Keep bottle in a well ventilated place. Keep away from children and also keep away from sources of ignition. No smoking. Bottle has to be closed always. After contact with skin, wash immediately with plenty of water and soap. Do not apply into drains. Do not throw to refuse if bottles are not empty. Pay attention for the rules of toxic waste Storage temp max 25 °C / 77 °F, Flashpoint 10 °C / 40 °F

Material will not bond to synthetic polymer teeth. For mechanical retention to synthetic polymer teeth use special PERform retention burr ref. no. DA-CU 1

Hedent GmbH,
Obere Zeil 6 – 8,
D – 61440 Oberursel, Germany

Tel. +49 6171 52036
Fax. +49 6171 52090
Email: info@hedent.de

performd.odt Stand 16.02.2021