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GB **KKD® KONDISIL - Condensation Type Silicone Impression System**

KKD® KONDISIL V-1 Putty - Silicone Tray Material
High-viscous Silicone Tray Material for preliminary impressions. Also for bite registrations and for laboratory procedures (as a space maintainer, as an insulating material etc.) very suitable. V-1 Putty material achieves in a few minutes an optimal shore hardness value and simultaneously develops exceptional elasticity. - Colour: white

V-1 Putty Mixing rate:
1 measuring spoon (21 g) V-1 Putty to 2,5 cm (0,5 g) paste accelerator
Mixing time: 1 minute
Working time: 2 minutes
Setting time in mouth: 3 ½ minutes
Complete setting time: 5 ½ minutes

Physical dates:
Shore hardness immediately: 70
Shore hardness after 2 hours: 79
Max. strain in compression: 1,3%
Memory after deformation: 98,6%
Linear dimensional change: 0,2%

Test conducted at room temperature of 23° C and relative humidity of 50 % (+/- 5%).
EN ISO 4823:2000 Type 0

Presentation: Pck.: 1,6 kg with hardener
REF: 18000S

Individuals:
Jar: 1,6 kg without hardener 17998
Pail: 10 kg without hardener 18001
(please order approx. 7-10 units hardener to it)
35 ml paste accelerator (red) 18003

as model material we recommend to use, according to final use, class III (e.g. Topdent Hartgips, KKD Cal) and class IV (e.g. KKD Satin Roc, KKD Veradur 03) gypsum

KKD® KONDISIL V-3 light - Silicone Impression Material
This low-viscous Silicone Impression Material for refined impressions has an excellent fluidity and a precise detail reproduction. It can be used with an impression tray or with a syringe. When used in combination with the V-1 putty material an optimal compound is guaranteed. - Colour: blue

V-3 low-viscous Mixing rate:
8 cm (6 g) V-3 base material to 5 cm (1 g) paste accelerator
Mixing time: 45 to 60 sec.
Working time: 2 minutes
Setting time in mouth: 3 ½ minutes
Complete setting time: 5 ½ minutes

Physical dates:
Shore hardness immediately: 33
Shore hardness after 2 hours: 39
Max. strain in compression: 12,2%
Memory after deformation: 99,2%
Linear dimensional change: 1,1%

Test conducted at room temperature of 23° C and relative humidity of 50 % (+/- 5%).
EN ISO 4823 Type 3

Presentation: Pkg.: 200 ml with hardener and mixing pad
REF: 18008S

Individuals:
200 ml V-3 base material 18009
35 ml paste accelerator (red) 18003

as model material we recommend to use, according to final use, class III (e.g. Topdent Hartgips, KKD Cal) and class IV (e.g. KKD Satin Roc, KKD Veradur 03) gypsum

- smooth consistency
- easy to mix and to handle
- precise reproductions of details and borders
- great tensile strength
- elastic
- dimensionally stable
- mint-scented
- long shelf life

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