**DENOMINAZIONE / PRODUCT
NAME / NOM DU PRODUIT /
PRODUKTBEZEICHNUNG /
DENOMINACIÓN / NOME DO
PRODUTO****DESCRIZIONE / PRODUCT
DESCRIPTION / DESCRIPTION /
PRODUKTBESCHREIBUNG /
DESCRIPCIÓN / DESCRIÇÃO DO
PRODUTO****OMNIFOL**

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OMNIFOL sterilization rolls are used for the sterilization of medical-surgical devices and of materials for medical use which are especially suitable for steam or EtO GAS. They are made by 2 layers: one made of high resistance medical paper and the other made of blue polyester/polypropylene. Sealing temperature between 170 °C and 200 °C. The seal is made by two parallel strings which are separated by a groove for rolls with a width of 5 cm. The seal is made by three parallel strings which are separated by two grooves for rolls with other dimensions. On the paper layer two chemical indicators of process are printed: one for steam and the other for ethylene-oxyde, both indicating that the sterilization process has been carried out.

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Storage: sterilized pouches and rolls must be handled as less as possible, kept possibly in closed cabinets in a cleaned environment, at a temperature and humidity under control. If maintained as described above, at a temperature between 8 and 35 °C and a relative humidity of between 30 and 70%, OMNIFOL packaging systems will keep the contents sterile for a maximum of 6 months and, in any case, no later than the indicated expired date (table 1 DIN 58953-8:2010).

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Rolls 200 metres long available in the following widths: 5 – 7,5 – 10 cm.
Roll 100 metres, 20 x 5 cm.

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