
astron® Labtec Pro

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MMA-free denture material with
pleasant wearing comfort

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- easy to handle
- no need for special devices
- **licence-free**
- natural colour
- excellent wearing comfort
- very break-resistant
- easy to reline
- suitable for all manufactured teeth

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astron Labtec Pro is an unique denture resin, free of MMA, flexible and due to its translucency very aesthetic. It can be processed in the pouring technique or as press packing. The best indication for dentures, e.g. for patients with allergies.

Economy Kit
320 g Pulver, 240 g Flüssigkeit
320 g powder, 240 g liquid REF 26210S

Super Economy Kit
960 g Pulver, 720 g (3 x 240 g) Flüssigkeit
960 g powder, 720 g (3 x 240 g) liquid REF 26215S

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Store powder and liquid in the refrigerator at 4 – 8 °C. (Practical tips: prior to use put mixing bowl into the refrigerator for 30 minutes to extend processing time – we recommend the use of ceramic – or glass bowls to avoid any chemical bond with the material). In order to avoid uncontrollable chemical reactions please use only mixing bowls and - spatulas which do not get in contact with other resins. For RELININGS and REPAIRS we recommend to use Labtec Reline.

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1. prepare model and wax up denture in regular/usual manner (use natural stone class III)
2. prepare protection wall and let it cure
3. remove it and prepare on the tubers, resp. on both sides of the molar triangle two overflow canals
4. boil plaster model out and water it
5. isolate the still warm plaster model
6. remove the isolation immediately under running tap water and dry it
7. sandblast the teeth, make retentions and reset the teeth back in the protection wall
8. fix the model with the protection wall and put it into the refrigerator for 15 minutes
9. powder and liquid have to be precisely-weighed and then put the appropriate liquid- and powder amount for approx. 15 minutes in the refrigerator
10. **Mixing ratio: ideally, use the measuring spoons included in the package**, or 1.0 g powder to 0.85 g liquid or a multiple, using a cold ceramic bowl (for one denture 20 g powder to 17 g liquid are required), spatulate slowly and bubble-free
11. pour resin immediately in a thin stream into the protection wall, immerse in pressure pot and cure for 60 minutes at a temperature of 52°C and 2 bar
12. remove denture from pressure pot and prior to open the unit, cool it for 15 minutes under cold running tap water
13. after reocclusion remove denture from plaster model and finish and polish the denture as usual
Attention: Do not steam off the denture, otherwise the high-gloss of the surface will be damaged and becomes matt. Ultrasonic cleaning and rinsing with fresh water and detergent are sufficient.
Hint: for easier removal of the denture from the plaster model, put model without „protection wall“ for approx. 30 seconds in warm water and remove it carefully

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1. prepare model and wax up denture in regular/usual manner (use natural stone class III)
2. invest model into flask
3. prepare counter part and after curing heat flask in water bath
4. open the flask, remove wax, boil out model and counter part, water and isolate them
5. sandblast the teeth, make retentions and reset the teeth back in the flask
6. put the flasks into the refrigerator for 15 minutes
7. remove powder and liquid from the refrigerator just prior to mixing
8. **Mixing ratio: ideally, use the measuring spoons included in the package**, or 1.0 g powder to 0.85 g liquid or a multiple, using a cold ceramic bowl (for one denture ,30 g powder to 25,5 g liquid are required), spatulate slowly and bubble-free
9. important: as soon as mix reaches a doughy state (5 – 8 minutes) pack into regular denture mold using normal press-packing procedure
10. close mold slowly to permit resin to flow. Place clamp and remove it together with flask
11. immerse flask in 52°C warm water and cure for 2 hours
12. cool flask for minimum 15 minutes under cold running tap water before removing denture from mold
13. after reocclusion remove denture from plaster model and finish and polish the denture as usual
Attention: Do not steam off the denture, otherwise the high-gloss of the surface will be damaged and becomes matt. Ultrasonic cleaning and rinsing with fresh water and detergent are sufficient.
Hint: for easier removal of the denture from the plaster model, put model for approx. 30 seconds in warm water and remove it carefully

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Care of the astron® Labtec Pro denture: We recommend astron® DuraClean® – REF 26170S

HINT: We recommend to water the finished denture for approx. 3 – 4 hours, then to rinse it under running water with a soft brush and to disinfect it afterwards with an alcohol-free disinfectant.

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CH REP RODENT AG
Astron Dental Corp.