# **FREEPRINT®**

## **3D PRINTING MATERIAL**

### **INDIVIDUAL FUNCTIONAL TRY-INS**

# **DETAX**

|           |                             |    |
|-----------|-----------------------------|----|
| <b>DE</b> | Gebrauchsanweisung.....     | 3  |
| <b>EN</b> | Instructions for use .....  | 6  |
| <b>FR</b> | Mode d'emploi .....         | 9  |
| <b>ES</b> | Modo de empleo.....         | 12 |
| <b>IT</b> | Istruzioni per l'uso.....   | 15 |
| <b>DA</b> | Brugsanvisning.....         | 18 |
| <b>EL</b> | Οδηγίες χρήσης.....         | 21 |
| <b>HR</b> | Upute za uporabu .....      | 24 |
| <b>LV</b> | Lietošanas instrukcija..... | 27 |
| <b>LT</b> | Naudojimo instrukcija ..... | 30 |
| <b>NL</b> | Gebruiksaanwijzing.....     | 33 |
| <b>NO</b> | Brugsanvisning.....         | 36 |
| <b>PL</b> | Instrukcja stosowania ..... | 39 |

|           |                                 |    |
|-----------|---------------------------------|----|
| <b>PT</b> | Manual de instruções.....       | 42 |
| <b>RO</b> | Instrucțiuni de utilizare ..... | 45 |
| <b>SV</b> | Brugsanvisning .....            | 48 |
| <b>SK</b> | Návod na použitie.....          | 51 |
| <b>SL</b> | Navodila za uporabo .....       | 54 |
| <b>CS</b> | Návod k použití.....            | 57 |
| <b>HU</b> | Használati útmutató.....        | 60 |
| <b>RU</b> | Инструкция по применению .....  | 63 |
| <b>TR</b> | Kullanım kılavuzu .....         | 66 |
| <b>FI</b> | Käyttöohje.....                 | 69 |
| <b>ET</b> | Kasutusjuhend .....             | 72 |
| <b>BG</b> | Инструкция за употреба.....     | 75 |

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Resin for dental 3D printing

**INDICATION**

Individual functional try-ins of complete dentures

**PATIENT TARGET GROUP**

Persons being treated in the context of a dental procedure.

**INTENDED USERS**

Dentist, dental technician

**SUITABLE FOR THE FOLLOWING DLP PRINTERS**

"Annex 1" (enclosed separately)

**PROCESSING**

- ▶ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the process description).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ A layer thickness of 100 µm is suitable for the printing process.
- ▶ Polish the surface mechanically if necessary. Pre-polishing is done with rotating brushes and pre-polishing paste, high-gloss polishing with buffing and high-gloss polishing agents for resin.
- ▶ Processing temperature 23 °C ± 2 °C.

**SAFETY INFORMATION**

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ Biocompatibility is guaranteed only after complete polymerization.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

**NOTES**

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Read and understand the safety data sheet

**For users and/or patients:**

Any serious incident that has occurred in connection with this product must be reported immediately to [incident@detax.com](mailto:incident@detax.com) and to the competent public authority of the Member State in which the user and/or patient resides.

## STORAGE

- Store **FREEPRINT® tryin** in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

## CONTRAINDICATION

Contains (meth)acrylates and phosphine oxides.

Ingredients of **FREEPRINT® tryin** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply **FREEPRINT® tryin** intraorally only when fully polymerized.

## SIDE EFFECTS

Product may cause allergic reactions.

## DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

### Storage:

### Processing:

at 23 °C ± 2 °C

# Manufacturing process

Data processing and creation of the support structure according to the specifications of the CAD software manufacturer

## Building process

Creation of a print job in compliance with the machine and material parameters.

## Finishing process

After raising the platform, a dripping time of approx. 10 minutes is recommended. If possible, the finishing should take place immediately after the building process.

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**DETAX GmbH**