## KariesScreenTest +P („KST +P“) order n° MD2427 Instructions for Use (IFU; English language)

### Introduction

The KST after Laurisch et al. is used for the identification and semiquantitative count of *streptococcus mutans* and *lactobacillus species* from human saliva specimen applying selective culture media *in vitro*.

A modified Rogosa-Agar (light color) is used to select for *lactobacillus species*.

A modified Mitis Salivarius-Agar (blue color) is used to select for *streptococcus mutans*.

The +P Salivatest after Ericson et al. is used for (a) the quantification of salivation over time and (b) measurement of saliva pH in the absence and presence of acidic buffer.

### Indication

KST +P is used to identify bacteria species that are strongly associated with caries development, and for determination of the saliva buffering capacity which is lowered in a cariogenic environment.

KST +P is an *in vitro* diagnostic to be used by skilled medical personnel, only.

### Warnings

For medical professional use only.

Store out of the reach of children. Do not eat and not drink. Read these instructions for use and the safety data sheets before using KST +P.

### Contraindications

Patient sample for the performance of KST+P may not be taken (1) during local/oral or systemic antibiotic therapy and two weeks thereafter, and not (2) within 12 hours after an antibiotic mouth rinse.

Patients are requested to not do one of the following things within one hour prior to specimen collection: -no drinking, -no eating; -no chewing gum; -no smoking; -no tooth brushing.

### Test Procedure

Take out one vial of KST media and one vial of +P Salivatest per patient from the box in the fridge and check the agar surfaces excluding contamination or bacterial growth. Warm up the media slide to room temperature before use (30 min. at RT or 5 min at 37°C). KST +P should be performed by skilled medical personnel, only, in accordance with these instructions for use and good medical practice.

(1) Quantify the Saliva collected using the scale on the collection tube or a measuring pastette.

(2) Measure the Saliva pH using one of the two pH strips comparing against the scale at the end of this IFU.

(3) Using the enclosed pastette transfer 1 ml of saliva into the buffering solution vial. Reclose the vial and let the saliva react with the buffering solution over 5 minutes. Do not shake!

(4) Reopen the vial and briefly introduce the other pH strip into the mix in the vial wetting as 3 reaction fields, and measure the pH of the mix representing the saliva buffering capacity comparing against the scale at the end of this IFU.

(5) Sample collection: For patients 6 years and older, let the patient chew one of the enclosed paraffin pellets over 5 minutes to stimulate salivation. Collect the saliva in the tube of the +P Saliva test after removing the buffering solution and the pH strips.

For patients 5 years or younger, collect saliva and/or plaque specimen using a fine brush or swab.

(6) Remove the agar carrier from the KST test vial.

(7) Place a NaHCO<sub>3</sub>-tablet at the bottom of the vial to generate appropriate culture conditions.

If test is performed with plaque only or very little saliva, add two drops tap water onto the tablet.

(8) Carefully remove both foils protecting the growth

media without touching the agar.

(9) Thoroughly, inoculate both agar media with saliva / the collected specimen wetting the entire agar surfaces. Do not scratch or punch the media surface.

(10) Slide the agar carrier back into the vial and screw-close the vial tightly.

(11) Mark the vial or the lid of the agar carrier with the patient name / code ensuring unique patient identification.

(12) Incubate the vial upright (lid up) over 48h at 37°C (99°F). For appropriate incubators, please refer to [www.aurosan-shop.de](http://www.aurosan-shop.de) or contact us for advice.

Perform the +P Saliva test as below (steps 10-13).

(13) Immediately after removing the vial from the incubator analyze the bacterial growth by comparing the density of the colony forming units (CFU) with the respective reference pictures below.

### Interpretation of results

Presence of both streptococcus mutans and lactobacillus in high numbers (>10<sup>5</sup> CFU) are associated with a high risk to develop caries.

saliva pH without stimulation: [6,5-6,9]

saliva pH after stimulation: [6,5-7,2]

buffering capacity of the saliva (pH after stimulation and buffering):

> 6.0: excellent; > 5.0: sufficient;

< 5.0: insufficient; < 4.0: poor.

The lower the buffering capacity the better the growth condition for procariogenic bacteria! Therapeutic and prophylactic measures should be discussed between dentist and patient. A communication aid that can be individualized is available at [www.revois-dental.com](http://www.revois-dental.com)

### Test performance and quality control

The growth media of KST +P are highly selective. The test may be interpreted only against the respective reference picture. For quality control the following control strains may be used:

Streptococcus mutans ATCC™ 25175 (e.g. AUROSAN order n° MD0266P)

Lactobacillus acidophilus ATCC™ 4356 (AUROSAN order n° MD0243P)

For the performance of quality controls, the necessary equipment and training, please refer to the AUROSAN mediathek at [www.aurosan.de](http://www.aurosan.de) or contact us, directly.

### Test disposal

Single use only! Contaminated or used vials need to be disposed in accordance with applicable national or local laws.

Some regulations may permit the disposal in regular trash after autoclavation. Appropriate autoclaving bags are available at AUROSAN (order n° MD-AL-00501).

### Storage conditions

KST +P should be stored at 2-8°C (36-46°F) – do not freeze. Protect from light and humidity. Do not use after the expiration dates printed on the outer box.

### Package Content and Order Information

Manufacturer: AUROSAN GmbH, Frankenstr. 231, D-45134 Essen, Germany

Tel +49 (0)201 21961 701

Fax +49 (0)201 21961 731

E-Mail [service@aurasan.de](mailto:service@aurasan.de)

[www.aurosan.de](http://www.aurosan.de) | [www.revois-dental.com](http://www.revois-dental.com)

Content MD2427 – KariesScreenTest +P:

5 vials with 2 media, 10 NaHCO<sub>3</sub> tablets, 10 paraffin pellets, 5 +P Saliva tests incl. collection tube, buffering solution, 2 pH-strips, 5 transfer pastettes, IFU

Additional product information incl. how to use – videos is available online at [www.kariesscreentest.de](http://www.kariesscreentest.de) [www.revois-dental.com](http://www.revois-dental.com) [www.revois-shop.de](http://www.revois-shop.de)

### Order information:

Please email your orders to [service@aurasan.de](mailto:service@aurasan.de) / [info@r-mt.de](mailto:info@r-mt.de) or to your local distributor

### Liability disclaimer

This *in-vitro* diagnostic is to be used by skilled medical personnel, only, and within the indication designated in this IFU. Processing, application and interpretation should strictly follow this IFU. Therapeutic or prophylactic measures recommended or initiated by the dentist as an interpretation of KST +P are the sole responsibility of the dentist. The user of KST +P is solely responsible for the suitability of material not delivered with the test and for its use for any purpose not explicitly stated in this IFU. Description, explanation, data or citation do not constitute any warranty of attributes and are not binding.

KST +P is made in Germany.

AUROSAN GmbH is certified in accordance with ISO13485.

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Hersteller: AUROSAN GmbH, Frankenstr. 231, D-45134 Essen, Germany