## KariesScreenTest („KST“) MD2426 Instructions for Use (IFU; English language)

### Introduction

The KST *after Laurisch et al.* is used for the identification and semiquantitative count of *streptococcus mutans* and *lactobacillus species* from human saliva specimen applying selective culture media *in vitro*.

A modified Rogosa-Agar (light color) is used to select for *lactobacillus species*.

A modified Mitis Salivarius-Agar (blue color) is used to select for *streptococcus mutans*.

### Indication

KST is used to identify bacteria species that are strongly associated with karies development.

KST is an *in vitro* diagnostic to be used by skilled medical personnel, only.

### Warnings

For medical professional use only.

Store out of the reach of children.

Read these instructions for use and the safety data sheets before using KST.

### Contraindications

Patient sample for the performance of KST may not be taken (1) during local/oral or systemic antibiotic therapy and two weeks thereafter, and not (2) within 12 hours after an antibiotic mouth rinse.

Patients must not do one of the following things within one hour prior to the specimen collection: -no drinking; -no eating; -no chewing gum; -no smoking; -no teeth brushing.

### Test Procedure

Take out one vial per patient from the box in the fridge and check the agar surfaces excluding contamination or bacterial growth. Warm up to RT before use (30 min. at room temperature or 5 min at 37°C).

KST should be performed by skilled medical personal, only, in accordance with these instructions for use and good medical practice.

(1) Sample collection: *For patients 6 years and older*, let the patient chew one of the enclosed paraffin pellets over 5 minutes to stimulate salivation. Collect the saliva in a suitable container. *For patients 5 years or younger*, collect saliva and/or plaque specimen using a fine brush or swab.

(2) Remove the agar carrier from the test vial.

(3) Place a NaHCO<sub>3</sub>-tablet at the bottom of the vial to generate appropriate culture conditions. If test is performed with plaque only or very little saliva, add two drops tap water onto the tablet.

(4) Carefully remove both foils protecting the growth media without touching the agar.

(5) Thoroughly, inoculate both agar media with saliva / the collected specimen wetting the entire agar surfaces. Do not scratch or punch the media surface.

(6) Slide the agar carrier back into the vial and screw-close the vial tightly.

(7) Mark the vial or the lid of the agar carrier with

the patient name / code ensuring unique patient identification.

(8) Incubate the vial upright (lid up) over 48h at 36 ± 1°C (approx. 99°F). For appropriate incubators, please refer to [www.aurosan-shop.de](http://www.aurosan-shop.de) or contact us for advise.

(9) Immediately after removing the vial from the incubator analyze the bacterial growth by comparing the density of the colony forming units (CFU) with the respective reference pictures below.

### Interpretation of results

Presence of both streptococcus mutans and lactobacillus in high numbers (>10<sup>5</sup> CFU) are associated with a high risk to develop caries.

Therapeutic and prophylactic measures should be discussed between dentist and patient. A communication aid that can be individualized is available at [www.revois-dental.com](http://www.revois-dental.com)

### Test performance and quality control

The growth media of KST are highly selective. The test may be interpreted only against the respective reference picture. For internal quality assurance the following control strains may be used: *Streptococcus mutans* ATCC™ 25175 (e.g. AUROSAN order n° MD0266P); *Lactobacillus acidophilus* ATCC™ 4356 (e.g. AUROSAN order n° MD0243P). For the performance of quality controls, the necessary equipment and training, please refer to the AUROSAN mediathek at [www.aurosan.de](http://www.aurosan.de) or contact us, directly.

### Test disposal

Single use only! Contaminated or used vials need to be disposed as biohazards in accordance with applicable national or local laws. Some regulations may permit the disposal in regular trash after autoclavation. Appropriate autoclaving bags are available (e.g. AUROSAN order n° MD-AL-00501).

### Storage conditions

KST should be stored at 2-8°C (36-46°F) – do not freeze. Protect from light and humidity. Do not use after the expiration date printed on the agar carrier.

### Package Content and Order Information

Manufacturer: AUROSAN GmbH, Frankenstr. 231, D-45134 Essen, Germany  
Tel +49 (0)201 21961 701  
Fax +49 (0)201 21961 731  
E-Mail [service@aurasan.de](mailto:service@aurasan.de)  
[www.aurosan.de](http://www.aurosan.de) | [www.revois-dental.com](http://www.revois-dental.com)

Content MD2426 – KariesScreenest:  
10 vials with 2 media, 10 NaHCO<sub>3</sub> tables, 10 paraffin pellets, IFU  
*NaHCO<sub>3</sub> tables and paraffin pellets can be purchased separately, as well.*

Extensive product information online available at [www.kariesscreentest.de](http://www.kariesscreentest.de)  
[www.revois-dental.com](http://www.revois-dental.com)  
[www.revois-shop.de](http://www.revois-shop.de)

Purchase Orders: Please e-mail your orders to [service@aurasan.de](mailto:service@aurasan.de) / [info@r-mt.de](mailto:info@r-mt.de) or to your local dealer / dental depot.

### Liability disclaimer

This *in-vitro* diagnostic is to be used by skilled medical personell, only, and within the indication designated in this IFU. Processing, application and interpretation should strictly follow this IFU.

Therapeutic or prophylactic measures recommended or initiated by the dentist as an interpretation of KST are the sole responsibility of the dentist. The user of KST is solely responsible for the suitability of testing material not delivered with the test and its use for any purpose not explicitly stated in this IFU. Description, explanation, data or citation do not constitute any warranty of attributes and are not binding.

### KST is made in Germany.

AUROSAN is certified in accordance with ISO13485.

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Hersteller: AUROSAN GmbH