





## UNIGLOVES COMFORT®

SIZE	ITEM-NO
XS	1201
S	1202
M	1203
L	1204
XL	1205



## EIGENSCHAFTEN

<b>MATERIAL</b>	<i>Latex</i>
<b>COLOUR</b>	<i>natural latex</i>
<b>LÄNGE</b>	<i>240 mm</i>
<b>SIZE</b>	<i>XS (5.5-6.0), S (6.5-7.0), M (7.5-8.0), L (8.5-9), XL (9.5-10)</i>
<b>AQL</b>	<i>1.5</i>
<b>PRODUCT FEATURES</b>	<i>beaded cuff, chlorinated on the inside &amp; outside, non-sterile, powder-free, smooth palm</i>
<b>NORM GLOVES</b>	<i>CE, EN 374, EN 455, T. 1-4</i>

**VERPACKUNGSIHALT**

box of 100 pieces

**VERSANDEINHEIT**

10 boxes

**Dieses Dokument wird Ihnen zur Verfügung gestellt von:**

Med-Dent24 Handelsgesellschaft mbH  
Gotenweg 7  
63128 Dietzenbach

Telefon: +49 (0)6074 803 3615  
Telefax: +49 (0)6074 803 4092  
E-Mail: [info@med-dent24.com](mailto:info@med-dent24.com)

Registergericht: Offenbach am Main  
Registernummer: HRB 56518  
Umsatzsteuer-Identifikation-No.: DE307869540

Geschäftsführerin: Ana Ban