

D600

Ultrasonic Scaler

Pleasant & Efficient




Intelligent power output
Experience painless treatment

Ultrasonic Endodontic Expert

Rapid preparation, root canal orifice locating, removal of foreign body in root canal, endodontic irrigation, root canal retreatment, supragingival scaling and periodontal treatment.


Compatible with more than 70 DTE tips.
Metal handpiece with LED can be autoclaved.



Sine-wave Drive

Avoid stimulation caused by instantaneous current output. Bring patients with comfortable treatment. Embody comfortable scaling.

D600 mild sine wave power output



High configuration

GD5 x2 GD12 x2 PD4 x2 PD7 x2 PD12 x2 ED60 x2 ED66 x2

Standard configuration

GD5 x1 GD12 x1 PD4 x1 PD7 x1 PD12 x1 ED60 x1 ED66 x1

Note: Woodpecker reserves the right of modifying the tip configuration without further notice.




Intelligent degenerative feedback

According to the different surface state of teeth, automatically select correct power so as to realize stable and continuous scaling effect.


Smart touch system

Clear control panel display. Ultrasensitive touch control. Waterproof design.




Automatic water supply system

Solutions such as sodium hypochlorite, chlorhexidine, hydrogen peroxide, and saline can be used. Easy to change water; double waterway design eliminates redundant steps. Clear 50 ml scale mark makes it convenient to formulate solutions.



Detachable handpiece with LED

Handpiece with LED. Handpiece adopts latest sealing technique and can be autoclaved. Patients feel more comfortable during treatment.




D600 LED with light / Without light


Supragingival scaling





Implant maintenance



Subgingival scaling



Root canal treatment



Cavity preparation



GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.

Registration Add & Manufactory Add:
Information Industrial Park, Guilin National High-Tech
Information Park 541004 Guilin, Guangxi China
Tel:
Europe Sales Dept.: +86-773-5873196, +86-773-2125222
North America, South America &
Oceania Sales Dept.: +86-773-5873198, +86-773-2125123
Asia & Africa Sales Dept.: +86-773-5855350, +86-773-2125896
Fax: +86-773-5822450
E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com
Website: <http://www.glwoodpecker.com>
ZMM/W1-10-098 V1.0

Distributor:



Scan and login website
for more information



Dieses Dokument wird Ihnen zur Verfügung gestellt von:

Med-Dent24 Handelsgesellschaft mbH
Gotenweg 7
63128 Dietzenbach

Telefon: +49 (0)6074 803 3615
Telefax: +49 (0)6074 803 4092
E-Mail: info@med-dent24.com

Registergericht: Offenbach am Main
Registernummer: HRB 56518
Umsatzsteuer-Identifikation-No.: DE307869540

Geschäftsführerin: Ana Ban