

Robust. Powerful. Efficient.

KaVo **expert** series

KaVo EXPERTmatic series

- **Low head dimensions** for easy and efficient access to the preparation area
- **Always ready for use** thanks to replaceable water filters (EXPERTmatic E25 L / E25 C)
- **Easy Instrument change** – quickly ready for use thanks to the flexible snap-in nose
- **2 high-speed attachments** – transmission ratio 1:5, E25 L with light, E25 C without light



KaVo EXPERTmatic
LUX E25 L



KaVo EXPERTmatic
E25 C



KaVo EXPERTmatic handpieces and contra-angles

KaVo EXPERTmatic LUX E25 L / E25 C

- High-speed attachments, transmission ratio 1:5
- Fast and precise work is possible using a small head, snap-in nose and ceramic ball bearings
- Maximum reliability through internal spray media and replaceable water filter
- E25 L with glass rod light conductor (25,000 LUX), E25 C without light



Special EXPERTmatic features:

- EXPERTgrip surface for effective gripping
- Exchangeable water filter for consistently good spray
- Light output up to 25,000 lux
- Chuck retention force of up to 30 N
- Carbide guide bushing for perfect drill concentricity
- Original KaVo ceramic ball bearing
- 18-month warranty
- Made in Germany

KaVo EXPERTmatic LUX E20 L / E20 C

- Transmission ratio 1:1
- Easy access due to a small head
- Efficient design with single spray and internal spray media
- E20 L with glass rod light conductor, E20 C without light



KaVo EXPERTmatic LUX E15 L / E15 C

- Reduction 5.4:1
- Easy access due to a small head
- Efficient design with single spray and internal spray media
- E15 L with glass rod light conductor, E15 C without light



KaVo EXPERTmatic E10 C

- Handpiece, transmission ratio 1:1
- Efficient design with single spray and internal spray media
- Suitable for handpiece and contra-angle drill



KaVo EXPERTmatic E31 C

- Polishing contra-angle reduction 5.4:1
- Multiple seals protect against ingress of polishing paste
- With thread and insert
- For screw-in and snap-on cup



KaVo EXPERTmatic E61 C

- EVA-PROPHY contra-angle
- Ideal for polishing interdental spaces, etc.
- 36 lock settings for precise locking of files



KaVo expert series

Article	Description	Material number
E25 L	EXPERTmatic high-speed attachment LUX, 1:5	1.007.5550
E25 C	EXPERTmatic high-speed attachment (without light), 1:5	1.007.5551
E20 L	EXPERTmatic LUX, 1:1	1.007.5540
E20 C	EXPERTmatic contra-angle (without light), 1:1	1.007.5541
E15 L	EXPERTmatic LUX, 5.4:1	1.007.5530
E15 C	EXPERTmatic contra-angle (without light), 5.1:1	1.007.5531
E10 C	EXPERTmatic handpiece, 1:1	1.007.5560
E31 C	EXPERTmatic polishing contra-angle, 5.4:1	1.008.4998
E61 C	EXPERTmatic EVA-PROPHY contra-angle	1.008.4999



KaVo. Dental Excellence.

Dieses Dokument wird Ihnen zur Verfügung gestellt von:

Med-Dent24 Handelsgesellschaft mbH
Gotenweg 7
63128 Dietzenbach

Telefon: +49 (0)6074 803 3615
Telefax: +49 (0)6074 803 4092
E-Mail: info@med-dent24.com

Registergericht: Offenbach am Main
Registernummer: HRB 56518
Umsatzsteuer-Identifikation-No.: DE307869540

Geschäftsführerin: Ana Ban