

KaVo MASTER series

Enjoy maximum power with innovative KaVo comfort.

1. Best treatment quality.
2. Maximum treatment comfort.
3. Highest patient satisfaction.

Handpieces and contra-angle handpieces

MASTERmatic™ series

Features of the MASTERmatic series:

- Plasmatic coating for perfect grip
- Replaceable spray microfilter
- Glass rod light conductor
- Externally guided, separate cooling media
- M25 L/M25 L: More space due to the 100°/19° head/knee angle combination and smaller heads
- 24 months warranty for MASTERmatic and INTRA™ heads
- *except M10 L

MASTERmatic LUX speed increasing contra-angle handpiece 1:5 transmission ratio
 Maximum treatment flexibility and options provided by 14 exchangeable heads
 M20 L | Mat. No. 1.009.3580

MASTERmatic LUX LUX 2:1
 Maximum treatment flexibility and options provided by 14 exchangeable heads
 M07 L | Mat. No. 1.009.3670

MASTERmatic LUX LUX 2:1
 Maximum treatment flexibility and options provided by 14 exchangeable heads
 M07 L | Mat. No. 1.009.3670

COMFORTdrive™
 High-speed contra-angle handpiece – very light, weight and handy – integrated brushless motor – glass rod light conductor – ceramic ball bearings – replaceable spray microfilter – 3-port spray – for direct connection only to KaVo CDPH/RTBase – max. speed: 200,000 rpm – 24 months warranty
 200 XD | Mat. No. 1.000.5500

MASTERmatic LUX speed increasing contra-angle handpiece 1:5 transmission ratio
 Miniature head specifically for short-shank bars – 3-port spray – push-button chuck – ceramic ball bearing – head 4:3 & smaller than GENTLEFLOW LUX 25 LP – very low noise level (55 dBA) | Mat. No. 1.009.3640

MASTERmatic LUX reducing shank 7:4:1
 Maximum treatment flexibility and options provided by 14 exchangeable heads
 M20 L | Mat. No. 1.009.3580

MASTERmatic LUX LUX handpiece 1:1
 For both handpiece and contra-angle handpiece bars
 M10 L | Mat. No. 1.009.3570

COMFORTdrive contra-angle handpiece
 High-speed contra-angle handpiece – very light, weight and handy – integrated brushless motor – glass rod light conductor – ceramic ball bearings – replaceable spray microfilter – 3-port spray – for direct connection only to KaVo CDPH/RTBase – max. speed: 200,000 rpm – 24 months warranty
 200 XD | Mat. No. 1.000.5500

INTRA™ Heads

INTRA head 1:1
 Push-button chuck – for contra-angle handpiece bars – includes debris deflector disc and reinforced chuck system – max. speed: 40,000 rpm
 L68 | Mat. No. 1.008.1833

INTRA head 2:1
 Push-button chuck – for contra-angle handpiece bars – includes debris deflector disc and reinforced chuck system – max. speed: 20,000 rpm
 L67 | Mat. No. 1.008.1832

INTRA head 1:1
 Push-button chuck – ball bearings – for contra-angle handpiece bars and root canal instruments – suitable for function of the ESTETICA E70/E80 Vision in combination with MASTERmatic M20 L
 L68 B | Mat. No. 1.008.1834

INTRA head 3:1
 Push-button chuck – ball bearings – for contra-angle handpiece bars and root canal instruments – suitable for function of the ESTETICA E70/E80 Vision in combination with MASTERmatic M20 L
 L66 B | Mat. No. 1.008.1831

INTRA Mini head 1:1
 Push-button chuck – for miniature contra-angle handpiece bars (6 mm Ø/8 in) – max. speed: 40,000 rpm
 L52 | Mat. No. 1.008.1825

INTRA head 1:1
 Push-button chuck – for FG-burs Ø 1.6 mm (1/16 in) – max. speed: 40,000 rpm
 L80 | Mat. No. 1.004.4456

INTRA Micro head 1:1
 Round bur, size 2 and U-inverted cone bur size 0 – for microsurgery – max. speed: 20,000 rpm
 L22 | Mat. No. 1.008.1835

INTRA Prophy head 2:1
 With thread and knob – for screw-type and snap-on polishing attachments – with multiple seals to prevent penetration of polishing paste
 max. speed: 15,000 rpm
 L31 | Mat. No. 1.008.1824

INTRA EVA head
 Files rotate freely – with buffered 0.8 mm (1/32 in) stroke – for files by Intensiv – max. speed: 20,000 rpm
 L61 G | Mat. No. 1.008.1828

INTRA EVA head
 With 36 separate locking positions for locking files – with buffered 0.8 mm (1/32 in) stroke – for files by Intensiv – max. speed: 20,000 rpm
 L61 G | Mat. No. 1.008.1829

INTRA Endo head
 Push-button chuck – for manual root canal instruments with grip diameter 3.5–4.0 mm – lift motion 0.4 mm – instrument reciprocates 300° – max. speed: 3,000 rpm
 L3 | Mat. No. 1.008.1822

INTRA Endo head
 Push-button chuck – for manual root canal instruments with grip diameter 3.5–4.0 mm – lift motion 0.4 mm – instrument reciprocates 300° – max. speed: 3,000 rpm
 L3 | Mat. No. 1.008.1822

INTRA Endo/Prophy head
 Push-button chuck – for root canal instruments and polishing attachments with 2.35 mm shaft diameter (later) – reciprocating rotation through 60° – max. speed: 3,000 rpm
 L53 | Mat. No. 1.008.1826

Electric Motors / Air Motors / ELECTROmatic

Electric motors

INTRA LUX motor KL 703 LED
 Brushless electric motor – speed range from 100 to 40,000 rpm in connection with KaVo electronic – controlled torque – max. torque 2.7 Nm
 (24 in) for use with KaVo MULTiflex LED – interchangeable by user – 36 months warranty
 KL 703 | Mat. No. 1.007.9150

INTRA LUX motor KL 701
 Brushless electric motor – speed range from 2,000 to 40,000 rpm in connection with KaVo electronic – controlled torque – max. torque 2.7 Nm
 (24 in) for use with KaVo MULTiflex LED – interchangeable by user – 36 months warranty
 KL 701 | Mat. No. 1.001.8776

INTRAmatic™ L-LUX motor 181 L
 Speed range from 5,000 to 20,000 rpm – forward and backward running – steerable – 4-hole connection – KaVo high pressure bulb – interchangeable by user – 12 months warranty
 181 L | Mat. No. 0.535.5210

INTRAmatic™ L-motor 181 H
 Speed range from 5,000 to 20,000 rpm – forward and backward running – steerable – 4-hole connection – INTRAmatic 2-coupling – 12 months warranty
 181 H | Mat. No. 0.535.5200

INTRAmatic™ L-motor 181 M
 Speed range up to 20,000 rpm – forward and backward running – 4-hole connection – no optic – no water – no chip air – to be used with SMARTmatic K program – 12 months warranty
 181 K | Mat. No. 0.535.5200

ELECTROmatic™ – The conversion kit from air to electric

Air electric motor conversion kit (controller, tubing, brushless electric motor and fixtures) for any treatment unit. Easy retrofitting of existing air-driven treatment units. Allows the dentist to upgrade to an electric motor and all its advantages, within minutes. Speed control via foot pedal of central unit.
 Optional with INTRA LUX Motor (C-5, S-5) Versions or with COMFORTdrive (C-5, S-5) Version. With standard tube (750mm) or special long tube (2200mm) with different Ref. see price list.

ELECTROmatic Premium / Premium Dual
 Highest power and fast, precise operation – faster due to 3-Speed Programs – with 100 Hz Invert. Database – Reciprocating Movement – Color Touch Screen – Fast & easy set up
 TH-S | Mat. No. 1.012.8545
 TH-M-S | Mat. No. 1.012.8547
 TH-C-S | Mat. No. 1.012.8549

ELECTROmatic Plus
 Highest power and fast, precise operation – Speed control via foot pedal of central unit.
 PM-S | Mat. No. 1.012.8541
 PC-S | Mat. No. 1.012.8543

ELECTROmatic
 Highest power and fast, precise operation – Speed control via foot pedal of central unit.
 M-S | Mat. No. 1.012.8537
 C-S | Mat. No. 1.012.8539

Turbines

MASTERtorque™ series

Features of the MASTERtorque series:

- Plasmatic coating for perfect grip
- Replaceable spray microfilter
- Carbide retention force up to 32 N
- Carbide guiding bush
- Ceramic ball bearing
- Glass rod light conductor
- Perfect light and good access, small head and 100°/19° head/knee angle combination
- 24 months warranty

MASTERtorque LUX Turbine
 Maximum power 31 Watt – low noise level up to 57 dBA – head height 13.1 mm – 4-port spray
 M9000 L | Mat. No. 1.008.7900
 M9000 L51 | Mat. No. 1.008.5400
 M9000 L | Mat. No. 1.007.9900
 M9000 L | chocolate brown | Mat. No. 1.007.7100
 M9000 L | anthracite | Mat. No. 1.007.7100

MASTERtorque Mini LUX Turbine
 Extra small head, head height 12.1 mm – low noise level – up to 59 dBA – maximum power 30 Watt – 3-port spray – for standard and miniature burs
 M7000 L | Mat. No. 3.001.0000
 M7000 L51 | Mat. No. 3.001.0060

* This instrument fits a non KaVo coupling system (L5 = KaVo turbine fits to Sirona Click&Go-coupling)
 Maximum power data are measured at max. pressure settings for MULTiflex versions.

KaVo EXPERT series

Rely on advanced KaVo performance.

1. Maximum longevity.
2. Efficient operations.
3. Reliable features.

Handpieces and contra-angle handpieces

EXPERTmatic™ series

Features of the EXPERTmatic series:

- KaVo CDPH/RT for better grip
- Internally guided, separate cooling media
- E25 L/E20 L/ETSL: Glass rod light conductor
- 36 months warranty

EXPERTmatic LUX speed increasing contra-angle handpiece 1:5
 3-port spray – push-button chuck – ceramic ball bearings – replaceable spray microfilter
 E25 L | Mat. No. 1.007.5550

EXPERTmatic LUX contra-angle handpiece 1:1
 Push-button chuck – for contra-angle handpiece bars – debris deflector disc and reinforced chuck system
 E20 L | Mat. No. 1.007.5540

EXPERTmatic LUX reducing contra-angle handpiece 5:4:1
 Push-button chuck – for contra-angle handpiece bars – debris deflector disc and reinforced chuck system
 E15 L | Mat. No. 1.007.5530

EXPERTmatic speed increasing contra-angle handpiece with 1:5 transmission ratio
 3-port spray – push-button chuck – ceramic ball bearings – replaceable spray microfilter
 E25 C (without light) | Mat. No. 1.007.5551

EXPERTmatic contra-angle handpiece 1:1
 Push-button chuck – for contra-angle handpiece bars – debris deflector disc and reinforced chuck system
 E20 C (without light) | Mat. No. 1.007.5541

EXPERTmatic reducing contra-angle handpiece 5:4:1
 Push-button chuck – for contra-angle handpiece bars – debris deflector disc and reinforced chuck system
 E15 C (without light) | Mat. No. 1.007.5531

EXPERTmatic handpiece 1:1
 For both handpiece and contra-angle handpiece bars
 E70 C (without light) | Mat. No. 1.007.5560

EXPERTmatic prophylaxis contra-angle handpiece 2:7:1
 Prophylactic contra-angle handpiece for screw-type and snap-on polishing attachments
 E37 C (without light) | Mat. No. 1.008.4998

EXPERTmatic prophylaxis contra-angle handpiece 2:7:1
 Prep control head with 36 locking positions for polishing inserts – with buffered 0.8 mm (1/32 in) stroke – for files by Intensiv
 E61 C (without light) | Mat. No. 1.008.4999

EXPERTmatic handpiece 1:1
 For both handpiece and contra-angle handpiece bars
 E70 C (without light) | Mat. No. 1.007.5560

EXPERTmatic prophylaxis contra-angle handpiece 2:7:1
 Prophylactic contra-angle handpiece for screw-type and snap-on polishing attachments
 E37 C (without light) | Mat. No. 1.008.4998

EXPERTmatic prophylaxis contra-angle handpiece 2:7:1
 Prep control head with 36 locking positions for polishing inserts – with buffered 0.8 mm (1/32 in) stroke – for files by Intensiv
 E61 C (without light) | Mat. No. 1.008.4999

ELECTROmatic™ – The conversion kit from air to electric

Air electric motor conversion kit (controller, tubing, brushless electric motor and fixtures) for any treatment unit. Easy retrofitting of existing air-driven treatment units. Allows the dentist to upgrade to an electric motor and all its advantages, within minutes. Speed control via foot pedal of central unit.
 Optional with INTRA LUX Motor (C-5, S-5) Versions or with COMFORTdrive (C-5, S-5) Version. With standard tube (750mm) or special long tube (2200mm) with different Ref. see price list.

ELECTROmatic Premium / Premium Dual
 Highest power and fast, precise operation – faster due to 3-Speed Programs – with 100 Hz Invert. Database – Reciprocating Movement – Color Touch Screen – Fast & easy set up
 TH-S | Mat. No. 1.012.8545
 TH-M-S | Mat. No. 1.012.8547
 TH-C-S | Mat. No. 1.012.8549

ELECTROmatic Plus
 Highest power and fast, precise operation – Speed control via foot pedal of central unit.
 PM-S | Mat. No. 1.012.8541
 PC-S | Mat. No. 1.012.8543

ELECTROmatic
 Highest power and fast, precise operation – Speed control via foot pedal of central unit.
 M-S | Mat. No. 1.012.8537
 C-S | Mat. No. 1.012.8539

Turbines

EXPERTtorque™ series

Features of the EXPERTtorque series:

- Replaceable spray microfilter
- Push-button chuck
- Chuck retention force up to 30 N
- Carbide guiding bush
- Ceramic ball bearing
- Perfect light, small head and 100°/19° head/knee angle combination
- 18 months warranty

EXPERTtorque LUX Turbine
 Maximum power 29 Watt – glass rod light conductor (25,000 Lux) – 4-port spray – head height 13.6 mm – noise level up to 62 dBA
 E680 L | Mat. No. 1.006.6700
 E680 L51 | Mat. No. 1.006.4600
 E680 L51 | Mat. No. 1.006.4700

EXPERTtorque Turbine
 Maximum power 29 Watt – 4-port spray – three-port spray – glass rod light conductor – maximum power 27 Watt – noise level up to 65 dBA
 E677 L | Mat. No. 1.007.3600
 E677 L51 | Mat. No. 1.006.0100

EXPERTtorque Mini LUX Turbine
 Extra small head, head height 12.2 mm – three-port spray – glass rod light conductor – maximum power 27 Watt – noise level up to 65 dBA
 E677 L | Mat. No. 1.007.3600
 E677 L51 | Mat. No. 1.006.0100

* This instrument fits a non KaVo coupling system
 L5 = KaVo turbine fits to Sirona Click&Go-coupling, LN = KaVo turbine fits to NSK FlexiQuick-coupling
 Maximum power data are measured at max. pressure settings for MULTiflex versions.

KaVo SMART series

Start safe with affordable KaVo quality.

1. Made in Germany.
2. Premium materials.
3. Entry level price for the world of KaVo.

Handpieces and contra-angle handpieces

SMARTmatic™ series

Features of the SMARTmatic series:

- Smaller head
- Ceramic sliders within the clamping system
- Rust-free stainless steel
- Extremely precise ball bearings
- Lightweight
- Optimum balance
- 12 months warranty

SMARTmatic S20 5 1:1 with external spray
 For bars with shaft Ø 2.35 mm – internal spray transmission to the SprayCip® – max. speed: 40,000 rpm
 S20 S1 | Mat. No. 1.011.6752
 Also available without spray
 S20 | Mat. No. 1.011.6750

SMARTmatic S20 5 1:1 with external spray
 For bars with shaft Ø 2.35 mm – internal spray transmission to the SprayCip® – max. speed: 40,000 rpm – for straight and contra-angle handpiece drills
 S20 S1 | Mat. No. 1.011.6752
 Also available without spray
 S20 | Mat. No. 1.011.6750

SMARTmatic ENDO S81 B:1
 For endodontic files with shaft Ø 2.35 mm – max. speed: 40,000 rpm
 S81 | Mat. No. 1.011.6780

SMARTmatic ENDO S321 3:2:1
 For endodontic files with shaft Ø 2.35 mm – max. speed: 40,000 rpm
 S321 | Mat. No. 1.011.6790

SMARTmatic PROPHY S31 B:1
 Screw-in and snap-on – max. speed: 20,000 rpm – double-threaded for screw-type polisher/brushes – unscrewable holder for snap-on polishing applications
 S31 | Mat. No. 1.011.6760

SMARTmatic PROPHY S19 B:1 B:1
 Doriol – max. speed: 20,000 rpm – for disposable prophylaxis angles
 S19 | Mat. No. 1.011.6740

SMARTmatic PROPHY S33 B:1 / 70° oscillating rotation direction
 Snap-on – max. speed: 20,000 rpm – with holder for snap-on polishing applications
 S33 | Mat. No. 1.011.6800

SMARTmatic S53 / Prophy, Endo B:1 / 70° oscillating rotation direction
 For bars with shaft Ø 2.35 mm – max. speed: 20,000 rpm – usable with polishing and root canal instruments
 S53 | Mat. No. 1.011.6810

MULTiflex™ Couplings

MULTiflex LED coupling
 4-hole LUX connection – spray ring regulation with lock settings – no lock back due to integrated self-extinguishing valve – self-exchangeable KaVo MULTiflex LED – 12 months warranty
 465 LED | Mat. No. 1.007.3197

MULTiflex LUX coupling
 4-hole LUX connection – spray ring regulation with lock settings – no lock back due to integrated self-extinguishing valve – self-exchangeable KaVo high pressure bulb – 12 months warranty
 465 LKN | Mat. No. 0.533.1550

MULTiflex LUX coupling
 F151 Housed cold light supply – integrated glass rod light conductor – 12 months warranty
 Housed | Mat. No. 0.533.1930

MULTiflex coupling
 4-hole standard connection – spray ring regulation with lock settings – anti-retrofitting valve – 12 months warranty
 465 RH | Mat. No. 0.533.1600

MULTiflex coupling
 4-hole standard connection – without spray regulation – anti-retrofitting valve – 12 months warranty
 467 | Mat. No. 0.533.1210

MULTiflex LUX coupling
 4-hole standard connection without spray adjustment ring – anti-retrofitting valve – 12 months warranty
 460 LED | Mat. No. 1.007.3201
 460 LE | Mat. No. 1.001.7599

MULTiflex coupling
 4-hole standard connection without spray adjustment ring – 12 months warranty
 460 L | Mat. No. 1.001.7600

Prophylaxis and Minimally invasive applications

SONICflex™ 2008 quick
 For all applications ranging from prophylaxis to surgery – more than 50 different tips – rapid tip switch – frequency: 6,000 Hz – 3 adjustable power levels – 2006 L/L5: Glass rod light conductor – for KaVo or Sirona couplings – direction of motion: circular – includes torque wrench – SONICflex quick Set includes scaler tips no. SA, BA, GA – 24 months warranty
 SONICflex quick Set 2008 L | Mat. No. 1.007.8664
 SONICflex quick Set 2008 L51 | Mat. No. 1.007.8665
 SONICflex quick Set 2008 L51 | Mat. No. 1.007.8665
 SONICflex quick handpiece 2008 L51 | Mat. No. 1.007.5440
 SONICflex quick handpiece 2008 (without light) | Mat. No. 1.005.9311

* This instrument fits a non KaVo coupling system (L5 = KaVo instrument fits to Sirona Click&Go-coupling)

Prophylaxis and Minimally invasive applications

SONICflex™ 2003
 For all applications ranging from prophylaxis to surgery – 50 different tips – frequency: 6,000 Hz – 3 adjustable power levels – direction of motion: circular – includes torque wrench – SONICflex Set includes scaler tips no. 5, 6, 7 – 2003 L: glass rod light conductor – 12 months warranty
 SONICflex LUX Set 2003 L | Mat. No. 1.008.8333
 SONICflex LUX handpiece 2003 L | Mat. No. 1.000.4243
 SONICflex Set 2003 (without light) | Mat. No. 1.008.8334
 SONICflex handpiece 2003 (without light) | Mat. No. 1.000.4246

Prophylaxis

SONICflex™ 2000 N
 Frequency: 6,000 Hz – amplitude: 100µm – direction of movement: circular – internal water cooling – suitable for MULTiflex coupling – includes torque wrench and SONICflex scaler tips no. 5 "Invert" – 12 months warranty
 2000 N | Mat. No. 0.570.5020

SMARTmatic PROPHY S31 B:1
 Screw-in and snap-on – max. speed: 20,000 rpm – double-threaded for screw-type polisher/brushes – unscrewable holder for snap-on polishing applications
 S31 | Mat. No. 1.011.6760

SMARTmatic PROPHY S19 B:1
 Doriol – max. speed: 20,000 rpm – for disposable prophylaxis angles
 S19 | Mat. No. 1.011.6740

SMARTmatic PROPHY S33 B:1 / 70° oscillating rotation direction
 Snap-on – max. speed: 20,000 rpm – with holder for snap-on polishing applications
 S33 | Mat. No. 1.011.6800

Prophylaxis and Minimally invasive applications

PROPHYflex™ 4
 An air-polishing device for the gentle cleaning of tooth surfaces – all parts reproducible – perfect ergonomics – clean treatment – prevents disturbing powder mist – flexible working – for KaVo and Sirona coupling – 12 months warranty
 PROPHYflex 4 Standard | Mat. No. 3.002.8000
 PROPHYflex 4 S | Mat. No. 3.002.8000
 PROPHYflex 4 S | Flamingo | Mat. No. 3.004.5950
 PROPHYflex 4 S | Wave | Mat. No. 3.002.8000
 PROPHYflex 4 S | Wave | Mat. No. 3.002.8000
 PROPHYflex 4 S | Wave | Mat. No. 3.004.5900
 PROPHYflex 4 S | Wave | Mat. No. 3.002.2200
 PROPHYflex 4 S | Lime | Mat. No. 3.004.5950

* This instrument fits a non KaVo coupling system (S = KaVo instrument fits to Sirona Click&Go-coupling)

PROPHYflexarts™
 Supragingival – spherical powder for safe removal of soft plaque – calcium carbonate neutral – 30 sticks (Mat. No. 1.010.2020 neutral – 4 bottles (250 g each) | Mat. No. 1.010.1790
 berry | Mat. No. 1.010.0015
 cherry | Mat. No. 1.010.0016
 mint | Mat. No. 1.010.0017
 PROPHYflex 4-Perio Kit | Mat. No. 1.011.9403
 * 80 sticks (15 g each)

Perio Kit
 Used for subgingival treatments with the PROPHYflex – removes biofilm in periodontal pockets – cleans gingiva – improved soft soft cleaning
 Perio Kit | Mat. No. 1.000.3721
 Perio Set | Mat. No. 1.010.2237
 PROPHYflex 4-Perio Kit | Mat. No. 1.011.9403
 * Tips can be reproduced up to 10 times

PROPHYflex Perio Powder
 Sub-/supragingival – implant cleaning – extra smooth – Tasteless – glycerine powder
 Perio | Mat. No. 1.009.3732

DIAGNOcam™ Vision Full HD
 Three instant diagnostic images from one shot. Combining an intraoral camera, a transilluminating device for caries detection and a fluorescence instrument for root caries check into one device. All captured in a single shot and instantly displayed in Full HD resolution – 12 months warranty
 DIAGNOcam – Vision Full HD | Mat. No. 1.011.9213

DIAGNOcam™
 Radiation-free device for early detection of caries – ideal for patient education and excellent monitoring – 12 months warranty
 DIAGNOcam™ | Mat. No. 1.001.9000

DIAGNOdent™ pen and DIAGNOdent™ Display
 Device for early detection of hidden fissure and interproximal caries – detection of concretions in periodontal pockets – wireless instrument with optional external display – pen/Display 12 months warranty
 DIAGNOdent™ pen | Mat. No. 1.002.7000
 DIAGNOdent™ Display | Mat. No. 1.004.5400

Surgery portfolio

MASTERlung LUX Wireless Surgical Unit
 Large colour touchscreen (5.7") – wireless surgical foot control – data documentation (SD memory card) – multi-program – integrated cooling pump – hose set – max. torque 80 Ncm – motor INTRA LUX 5600 LED brushless motor – speed range: 0 – 40,000 rpm, motor and motor cable: auto-storable & disinfectable – 12 months warranty
 MASTERlung LUX Wireless | Mat. No. 1.009.1200
 INTRA LUX 5600 LED | Mat. No. 1.008.8000

SURGmatic S10 L Pro speed increasing contra-angle handpiece 1:5
 Transmission ratio 1:5 with up to 200,000 rpm – all-rounder with high cutting performance for surgical and standard applications – 3-port spray – external cooling connection – Cool-Head technology – EXPERTrip – 12 months warranty
 S10 L Pro (with light) | Mat. No. 1.014.4500

SURGmatic S201 L Pro Surgical contra-angle handpiece 2:0:1
 Reduction ratio 2:0:1 – Hexagon chuck system for a long service life and a max. torque of 80 Ncm – optimized clamping system – removable head, easy to clean – EXPERTrip – small head with internal cooling system – external cooling connection – optimised toothing at the head and improved ball bearings – 12 months warranty
 S201 L Pro (with hexagon light) | Mat. No. 1.013.7541

SURGmatic S201 L Pro Surgical contra-angle handpiece 2:0:1
 Reduction ratio 2:0:1 – max. torque of 55 Ncm – optimized clamping system – removable head, easy to clean – EXPERTrip – small head with internal cooling system – external cooling connection – optimised toothing at the head and improved ball bearings – 12 months warranty
 S201 L Pro (with light) | Mat. No. 1.013.7540

INTRA Surgical handpiece 1:1
 Glass rod light conductor – for straight bars with 44.5 mm (2.34 in) standard length – external cooling connection – easy to dismantle – reinforced chuck system – 12 months warranty
 3610 N1 | Mat. No. 0.524.5600

INTRA Surgical handpiece 1:1
 Glass rod light conductor – for straight bars with 44.5 mm (2.34 in) standard length – external cooling connection – easy to dismantle – reinforced chuck system – 12 months warranty
 3610 N2 | Mat. No. 0.524.5610

INTRA Surgical handpiece 1:1
 For straight bars with 62.5 mm (2.34 in) special length – external cooling connection – easy to dismantle – reinforced chuck system – 12 months warranty
 3610 N2 | Mat. No. 0.524.5610

INTRA Surgical handpiece 1:1
 For straight bars with 70 mm (2.76 in) special length – external cooling connection – easy to dismantle – reinforced chuck system – 12 months warranty
 3610 N3 | Mat. No. 0.524.5620

Surgery portfolio

EXPERTlung LUX Surgical Unit
 Large colour display (5.7") – multifunctional foot control (cable length 2.5 m) – single-user-operation for saving data – integrated cooling pump – hose set – max. torque 80 Ncm – motor INTRA LUX 5600 LED brushless motor – speed range: 0 – 40,000 rpm, motor and motor cable: auto-storable & disinfectable – 12 months warranty
 EXPERTlung LUX Surgical Unit | Mat. No. 1.008.3500
 EXPERTlung LUX LED | Mat. No. 1.008.3500

INTRA Surgical handpiece 1:1
 Glass rod light conductor – for straight bars with 44.5 mm (2.34 in) standard length – external cooling connection – easy to dismantle – reinforced chuck system – 12 months warranty
 3610 N1 | Mat. No. 0.524.5600

INTRA Surgical handpiece 1:1
 Glass rod light conductor – for straight bars with 44.5 mm (2.34 in) standard length – external cooling connection – easy to dismantle – reinforced chuck system – 12 months warranty
 3610 N2 | Mat. No. 0.524.5610

INTRA Surgical handpiece 1:1
 For straight bars with 62.5 mm (2.34 in) special length – external cooling connection – easy to dismantle – reinforced chuck system – 12 months warranty
 3610 N2 | Mat. No. 0.524.5610

INTRA Surgical handpiece 1:1
 For straight bars with 70 mm (2.76 in) special length – external cooling connection – easy to dismantle – reinforced chuck system – 12 months warranty
 3610 N3 | Mat. No. 0.524.5620

Surgery portfolio

SURtorque 459 L
 45 degree head angle – for 25 mm bars – 19 W – no air into bur direction, were out in 2 side vents – ceramic ball bearing – fiber optic – standard head size – push-button chuck – glass rod light conductor – single water jet – suitable for MULTiflex coupling
 459 L | Mat. No. 3.000.4808

SURtorque 459 C
 45 degree head angle – for 25 mm bars – 19 W – no air into bur direction, were out in 2 side vents – ceramic ball bearing – fiber optic – standard head size – push-button chuck – glass rod light conductor – single water jet – suitable for MULTiflex coupling
 459 C | Mat. No. 3.000.5063

KaVo CLEANspray™/DRYSpray™
 Cleaning and drying spray for intensive manual, validated internal cleaning of instruments and fixtures – with certification of the removal of residual proteins

Black: KaVo CLEANspray (4 cans)
 2110 P1 | Mat. No. 1.007.5719
 Black: KaVo DRYSpray (4 cans)

Dieses Dokument wird Ihnen zur Verfügung gestellt von:

Med-Dent24 Handelsgesellschaft mbH
Gotenweg 7
63128 Dietzenbach

Telefon: +49 (0)6074 803 3615
Telefax: +49 (0)6074 803 4092
E-Mail: info@med-dent24.com

Registergericht: Offenbach am Main
Registernummer: HRB 56518
Umsatzsteuer-Identifikation-No.: DE307869540

Geschäftsführerin: Ana Ban