

KaVo MASTER series

Enjoy maximum power with innovative KaVo comfort.

1. Best treatment quality.
2. Maximum treatment comfort.
3. Highest patient satisfaction.

Handpieces and contra-angle handpieces

MASTERmatic™ series

Features of the MASTERmatic series:

- Plasma-coating for perfect grip
- Replaceable spray microfilter
- Glass rod light conductor
- Externally guided, separate cooling media
- M25 L/M25 L: More space due to the 100°/19° head/knee angle combination and smaller heads
- 24 months warranty for MASTERmatic and INTRA™ heads
- *except M10 L

MASTERmatic LUX speed increasing contra-angle handpiece 1:5 transmission ratio
 Maximum treatment flexibility and options provided by 14 exchangeable heads
 M20 L | Mat. No. 1.009.3580

MASTERmatic LUX LUX 2:1
 Maximum treatment flexibility and options provided by 14 exchangeable heads
 M07 L | Mat. No. 1.009.3670

MASTERmatic LUX LUX 2:1
 Maximum treatment flexibility and options provided by 14 exchangeable heads
 M07 L | Mat. No. 1.009.3670

MASTERmatic LUX speed increasing contra-angle handpiece 1:5 transmission ratio
 Minimum head specifically for short-shank bars - 3-port spray - push-button chuck - ceramic ball bearing - head 4:3 & smaller than GENTLEFLOW LUX 2:1 LP - 25 LP - very low noise level (55 dBA) | Mat. No. 1.009.3640

MASTERmatic LUX reducing shank 7:4:1
 Maximum treatment flexibility and options provided by 14 exchangeable heads
 M20 L | Mat. No. 1.009.3580

MASTERmatic LUX LUX handpiece 1:1
 For both handpiece and contra-angle handpiece bars
 M10 L | Mat. No. 1.009.3570

COMFORTdrive™ contra-angle handpiece
 High-speed contra-angle handpiece - very light, weight and handy - integrated brushless motor - glass rod light conductor - ceramic ball bearing - replaceable spray microfilter - 3-port spray - for direct connection only to KaVo COMFORTbase - max. speed: 200,000 rpm - 24 months warranty
 200 XD | Mat. No. 1.000.5500

INTRA™ Heads

INTRA head 1:1
 Push-button chuck - for contra-angle handpiece bars - includes debris deflector disc and reinforced chuck system - max. speed: 40,000 rpm
 L68 | Mat. No. 1.008.1833

INTRA head 2:1
 Push-button chuck - for contra-angle handpiece bars - includes debris deflector disc and reinforced chuck system - max. speed: 20,000 rpm
 L67 | Mat. No. 1.008.1832

INTRA head 1:1
 Push-button chuck - ball bearings - for contra-angle handpiece bars and root canal instruments - suitable for function of the ESTETICA E70/E80 Vision in combination with MASTERmatic M20 L
 L68 B | Mat. No. 1.008.1834

INTRA head 3:1
 Push-button chuck - ball bearings - for contra-angle handpiece bars and root canal instruments - suitable for function of the ESTETICA E70/E80 Vision in combination with MASTERmatic M20 L
 L66 B | Mat. No. 1.008.1831

INTRA Mini head 1:1
 Push-button chuck - for miniature contra-angle handpiece bars (6 mm Ø/8 in) - max. speed: 40,000 rpm
 L52 | Mat. No. 1.008.1825

INTRA head 1:1
 Push-button chuck - for FG-burs Ø 1.6 mm (1/16 in) - max. speed: 40,000 rpm
 L80 | Mat. No. 1.004.4456

INTRA Micro head 1:1
 Round bur, size 2 and U-inverted cone bur size 0 - for microsurgery - max. speed: 20,000 rpm
 L22 | Mat. No. 1.008.1835

INTRA Prophy head 2:1
 With thread and knob - for screw-type and snap-on polishing attachments - with multiple seals to prevent penetration of polishing paste
 max. speed: 15,000 rpm
 L31 | Mat. No. 1.008.1824

INTRA EVA head
 Files rotate freely - with buffered 0.8 mm (1/32 in) stroke - for files by Intensiv - max. speed: 20,000 rpm
 L61 G | Mat. No. 1.008.1828

INTRA EVA head
 With 36 separate locking positions for locking files - with buffered 0.8 mm (1/32 in) stroke - for files by Intensiv - max. speed: 20,000 rpm
 L61 H | Mat. No. 1.008.1829

INTRA Endo head
 Push-button chuck - for manual root canal instruments with grip diameter 3.5-4.0 mm - lift motion 0.4 mm - instrument reciprocates 300° - max. speed: 3,000 rpm
 L3 | Mat. No. 1.008.1822

INTRA Endo head
 Push-button chuck - for manual root canal instruments with grip diameter 3.5-4.0 mm - lift motion 0.4 mm - instrument reciprocates 300° - max. speed: 3,000 rpm
 L3 | Mat. No. 1.008.1822

INTRA Endo/Prophy head
 Push-button chuck - for root canal instruments and polishing attachments with 2.35 mm shaft diameter (later) - reciprocating rotation through 60° - max. speed: 3,000 rpm
 L53 | Mat. No. 1.008.1826

Electric Motors / Air Motors / ELECTROmatic

INTRA LUX motor KL 703 LED
 Brushless electric motor - speed range from 100 to 40,000 rpm in connection with KaVo electronic - controlled torque - max. torque 2.7 Nm
 (2.4 in) - fully replaceable - KaVo MULTI LED - exchangeable by user - 36 months warranty
 KL 703 | Mat. No. 1.007.9150

INTRA LUX motor KL 701
 Brushless electric motor - speed range from 100 to 40,000 rpm in connection with KaVo electronic - controlled torque - max. torque 2.7 Nm
 (2.4 in) - fully replaceable - KaVo MULTI LED - exchangeable by user - 36 months warranty
 KL 701 | Mat. No. 1.001.8776

INTRAmatic™ L-LUX motor 181 L
 Speed range from 5,000 to 20,000 rpm - forward and backward running - steerable - 4-hole connection - KaVo high pressure bulb exchangeable by user - 12 months warranty
 181 L | Mat. No. 0.535.5210

INTRAmatic™ L-motor 181 H
 Speed range from 5,000 to 20,000 rpm - forward and backward running - steerable - 4-hole connection - INTRAmatic 2-coupling - 12 months warranty
 181 H | Mat. No. 0.535.5200

INTRAmatic™ L-motor 181 M
 Speed range from 5,000 to 20,000 rpm - forward and backward running - steerable - 4-hole connection - INTRAmatic 2-coupling - 12 months warranty
 181 M | Mat. No. 0.535.5200

INTRAmatic™ K-motor 181 K
 Speed range up to 20,000 rpm - forward and backward running - 4-hole connection - non optic - no water - no chip air - to be used with SMARTmatic K program - 12 months warranty
 181 K | Mat. No. 0.535.5200

KaVo EXPERT series

Rely on advanced KaVo performance.

1. Maximum longevity.
2. Efficient operations.
3. Reliable features.

Handpieces and contra-angle handpieces

EXPERTmatic™ series

Features of the EXPERTmatic series:

- KaVo COMFORTgrip for better grip
- Internally guided, separate cooling media
- E25 L/E20 L/ E15 L: Glass rod light conductor
- 18 months warranty

EXPERTmatic LUX speed increasing contra-angle handpiece 1:5
 3-port spray - push-button chuck - ceramic ball bearing - replaceable spray microfilter
 E25 L | Mat. No. 1.007.5550

EXPERTmatic LUX contra-angle handpiece 1:1
 Push-button chuck - for contra-angle handpiece bars - debris deflector disc and reinforced chuck system
 E20 L | Mat. No. 1.007.5540

EXPERTmatic LUX reducing contra-angle handpiece 5:4:1
 Push-button chuck - for contra-angle handpiece bars - debris deflector disc and reinforced chuck system
 E15 L | Mat. No. 1.007.5530

EXPERTmatic speed increasing contra-angle handpiece with 1:5 transmission ratio
 3-port spray - push-button chuck - ceramic ball bearings - replaceable spray microfilter
 E25 C (without light) | Mat. No. 1.007.5551

EXPERTmatic contra-angle handpiece 1:1
 Push-button chuck - for contra-angle handpiece bars - debris deflector disc and reinforced chuck system
 E20 C (without light) | Mat. No. 1.007.5541

EXPERTmatic reducing contra-angle handpiece 5:4:1
 Push-button chuck - for contra-angle handpiece bars - debris deflector disc and reinforced chuck system
 E15 C (without light) | Mat. No. 1.007.5531

EXPERTmatic handpiece 1:1
 For both handpiece and contra-angle handpiece bars
 E70 C (without light) | Mat. No. 1.007.5560

EXPERTmatic prophylaxis contra-angle handpiece 2:7:1
 Prophylactic contra-angle handpiece for screw-type and snap-on polishing attachments
 E37 C (without light) | Mat. No. 1.008.4998

EXPERTmatic prophylaxis contra-angle handpiece 2:7:1
 Prep control head with 36 locking positions for polishing inserts - with buffered 0.8 mm (1/32 in) stroke - for files by Intensiv
 E61 C (without light) | Mat. No. 1.008.4999

KaVo SMART series

Start safe with affordable KaVo quality.

1. Made in Germany.
2. Premium materials.
3. Entry level price for the world of KaVo.

Handpieces and contra-angle handpieces

SMARTmatic™ series

Features of the SMARTmatic series:

- Smaller head
- Ceramic sliders within the clamping system
- Rust-free stainless steel
- Extremely precise ball bearings
- Lightweight
- Optimum balance
- 12 months warranty

SMARTmatic S20 5 1:1 with external spray
 For bars with shaft Ø 2.35 mm - internal spray transmission to the SprayCip® - max. speed: 40,000 rpm
 S20 S1 | Mat. No. 1.011.6752
 Also available without spray
 S20 | Mat. No. 1.011.6750

SMARTmatic S20 5 1:1 with external spray
 For bars with shaft Ø 2.35 mm - internal spray transmission to the SprayCip® - max. speed: 40,000 rpm - for straight and contra-angle handpiece drills
 S10 S1 | Mat. No. 1.011.6752
 Also available without spray
 S10 | Mat. No. 1.011.6730

SMARTmatic ENDO S81 B:1
 For endodontic files with shaft Ø 2.35 mm - max. speed: 40,000 rpm
 S81 | Mat. No. 1.011.6780

SMARTmatic ENDO S321 3:2:1
 For endodontic files with shaft Ø 2.35 mm - max. speed: 40,000 rpm
 S321 | Mat. No. 1.011.6790

SMARTmatic PROPHY S31 B:1
 Screw-in and snap-on - max. speed: 20,000 rpm - double-threaded for screw-type polishers/brushes - unscrewable holder for snap-on polishing applications
 S31 | Mat. No. 1.011.6760

SMARTmatic PROPHY S19 B:1 B:1
 Doriol - max. speed: 20,000 rpm - for disposable prophylaxis angles
 S19 | Mat. No. 1.011.6740

SMARTmatic PROPHY S33 B:1 / 70° oscillating rotation direction
 Snap-on - max. speed: 20,000 rpm - with holder for snap-on polishing applications
 S33 | Mat. No. 1.011.6800

SMARTmatic S53 / Prophy, Endo B:1 / 70° oscillating rotation direction
 For bars with shaft Ø 2.35 mm - max. speed: 20,000 rpm - usable with polishing and root canal instruments
 S53 | Mat. No. 1.011.6810

Electric Motors / Air Motors / ELECTROmatic

INTRA LUX motor KL 703 LED
 Brushless electric motor - speed range from 100 to 40,000 rpm in connection with KaVo electronic - controlled torque - max. torque 2.7 Nm
 (2.4 in) - fully replaceable - KaVo MULTI LED - exchangeable by user - 36 months warranty
 KL 703 | Mat. No. 1.007.9150

INTRA LUX motor KL 701
 Brushless electric motor - speed range from 100 to 40,000 rpm in connection with KaVo electronic - controlled torque - max. torque 2.7 Nm
 (2.4 in) - fully replaceable - KaVo MULTI LED - exchangeable by user - 36 months warranty
 KL 701 | Mat. No. 1.001.8776

INTRAmatic™ L-LUX motor 181 L
 Speed range from 5,000 to 20,000 rpm - forward and backward running - steerable - 4-hole connection - KaVo high pressure bulb exchangeable by user - 12 months warranty
 181 L | Mat. No. 0.535.5210

INTRAmatic™ L-motor 181 H
 Speed range from 5,000 to 20,000 rpm - forward and backward running - steerable - 4-hole connection - INTRAmatic 2-coupling - 12 months warranty
 181 H | Mat. No. 0.535.5200

INTRAmatic™ L-motor 181 M
 Speed range from 5,000 to 20,000 rpm - forward and backward running - steerable - 4-hole connection - INTRAmatic 2-coupling - 12 months warranty
 181 M | Mat. No. 0.535.5200

INTRAmatic™ K-motor 181 K
 Speed range up to 20,000 rpm - forward and backward running - 4-hole connection - non optic - no water - no chip air - to be used with SMARTmatic K program - 12 months warranty
 181 K | Mat. No. 0.535.5200

ELECTROmatic™ - The conversion kit from air to electric

Air electric motor conversion kit (controller, tubing, brushless electric motor and fixtures) for any treatment unit. Easy retrofitting of existing air-driven treatment units. Allows the dentist to upgrade to an electric motor and all its advantages, within minutes. Speed control via foot pedal of central unit.

ELECTROmatic Premium / Premium Dual
 Highest power and fast, precise operation - faster due to 3-Speed Program - with 100 Hz Invert. Database - Reciprocating Movement - Color Touch Screen - Fast & easy set up
 TH-S | Mat. No. 1.012.8545
 TH-M-S | Mat. No. 1.012.8547
 TH-C-S | Mat. No. 1.012.8549

ELECTROmatic Plus
 Highest power and fast, precise operation - Speed control via foot pedal of central unit
 PM-S | Mat. No. 1.012.8541
 PC-S | Mat. No. 1.012.8543

ELECTROmatic
 Highest power and fast, precise operation - Speed control via foot pedal of central unit
 M-S | Mat. No. 1.012.8537
 C-S | Mat. No. 1.012.8539

Turbines

MASTERtorque™ series

Features of the MASTERtorque series:

- Replaceable spray microfilter
- Plasma-coating for perfect grip
- Carbide retention force up to 32 N
- Ceramic guiding bush
- Ceramic ball bearing
- Glass rod light conductor
- Perfect light and good access, small head and 100°/19° head/knee angle combination
- 24 months warranty

MASTERtorque LUX Turbine
 Maximum power 31 Watt - low noise level up to 57 dBA - head height 13.1 mm - 4-port spray
 M9000 L | Mat. No. 1.008.7900
 M9000 L51 | Mat. No. 1.008.5400
 M9000 L (chocolate brown) | Mat. No. 1.007.9900
 M9000 L (anthracite) | Mat. No. 1.007.7100

MASTERtorque Mini LUX Turbine
 Extra small head - head height 12.1 mm - low noise level - up to 59 dBA - maximum power 30 Watt - 3-port spray - for standard and miniature burs
 M8700 L | Mat. No. 3.001.0000
 M8700 L51 | Mat. No. 3.001.0060

* This instrument fits a non KaVo coupling system (L5 = KaVo turbine fits to Sirona Click&Go-coupling)
 Maximum power data are measured at max. pressure settings for MULTIflex versions.

Turbines

EXPERTtorque™ series

Features of the EXPERTtorque series:

- Replaceable spray microfilter
- Push-button chuck
- Chuck retention force up to 30 N
- Carbide guiding bush
- Ceramic ball bearing
- Perfect light, small head and 100°/19° head/knee angle combination
- 18 months warranty

EXPERTtorque LUX Turbine
 Maximum power 29 Watt - glass rod light conductor (25,000 Lux) - 4-port spray - head height 13.6 mm - noise level up to 62 dBA
 E680 L | Mat. No. 1.006.6700
 E680 L51 | Mat. No. 1.006.4600
 E680 L51 | Mat. No. 1.006.4700

EXPERTtorque Turbine
 Maximum power 29 Watt - 4-port spray - three-port spray - glass rod light conductor - maximum power 27 Watt - noise level up to 65 dBA
 E677 L | Mat. No. 1.007.3600
 E677 L51 | Mat. No. 1.006.0100

EXPERTtorque Mini LUX Turbine
 Extra small head - head height 12.2 mm - three-port spray - glass rod light conductor - maximum power 27 Watt - noise level up to 65 dBA
 E677 L | Mat. No. 1.007.3600
 E677 L51 | Mat. No. 1.006.0100

* This instrument fits a non KaVo coupling system
 L5 = KaVo turbine fits to Sirona Click&Go-coupling, LN = KaVo turbine fits to NSK FlexiQuick-coupling
 Maximum power data are measured at max. pressure settings for MULTIflex versions.

Turbines

SMARTtorque™ series

SMARTtorque LUX Turbine
 Standard head - Glass rod light conductor - 26 Watt Power - ceramic bearings - made in Germany - push-button chuck - four-port spray
 S619 L | Mat. No. 1.008.1641
 S619 C (without light) | Mat. No. 1.008.1642

Fix coupling turbines

S609 C (without light) - 4-hole standard connector - 26 Watt Power - Midwest
 Mat. No. 1.008.1644

S605 C (without light) - 4-hole standard connector - 22 Watt Power - Midwest
 Mat. No. 1.008.1645

Maximum power data are measured at max. pressure settings.

MULTIflex™ Couplings

MULTIflex LED coupling
 4-hole LUX connection - spray ring regulation with lock settings - no lock back due to integrated self-extinguishing valve - self-exchangeable KaVo MULTI LED - 12 months warranty
 465 LED | Mat. No. 1.007.3172

MULTIflex LUX coupling
 4-hole LUX connection - spray ring regulation with lock settings - no lock back due to integrated self-extinguishing valve - self-exchangeable KaVo high pressure bulb - 12 months warranty
 465 LKN | Mat. No. 0.553.1550

MULTIflex LUX coupling
 F151 Housed cold light supply - integrated glass rod light conductor - 12 months warranty
 Housed | Mat. No. 0.553.1930

MULTIflex coupling
 4-hole standard connection - spray ring regulation with lock settings - anti-rotation valve - 12 months warranty
 465 RH | Mat. No. 0.553.1600

MULTIflex coupling
 4-hole standard connection - without spray regulation - anti-rotation valve - 12 months warranty
 460 LED | Mat. No. 1.003.1210
 457 | Mat. No. 0.553.1210

MULTIflex LUX coupling
 4-hole standard connection without spray adjustment ring - anti-rotation valve - 12 months warranty
 460 LED | Mat. No. 1.007.2021
 460 L | Mat. No. 1.001.7599

MULTIflex coupling
 4-hole standard connection without spray adjustment ring - 12 months warranty
 460 L | Mat. No. 1.001.7600

Prophylaxis and Minimally invasive applications

SONICflex™ 2008 quick
 For all applications ranging from prophylaxis to surgery - more than 50 different tips - rapid tip switch - frequency: 6,000 Hz - 3 adjustable power levels - 2006 L/L5: Glass rod light conductor - for KaVo or Sirona couplings - direction of motion: circular - includes torque wrench - SONICflex quick Set includes scaler tips no. SA, BA, GA - 24 months warranty
 SONICflex quick Set 2008 L | Mat. No. 1.007.8664
 SONICflex quick Set 2008 L51 | Mat. No. 1.007.8665
 SONICflex quick Set 2008 L51 | Mat. No. 1.007.8665
 SONICflex quick handpiece 2008 L51 | Mat. No. 1.007.5440
 SONICflex quick handpiece 2008 (without light) | Mat. No. 1.005.9311

* This instrument fits a non KaVo coupling system (L5 = KaVo instrument fits to Sirona Click&Go-coupling)

Prophylaxis and Minimally invasive applications

SONICflex™ 2003
 For all applications ranging from prophylaxis to surgery - 50 different tips - frequency: 6,000 Hz - 3 adjustable power levels - direction of motion: circular - includes torque wrench - SONICflex Set includes scaler tips no. 5,6,7 - 2003 L: glass rod light conductor - 12 months warranty
 SONICflex LUX Set 2003 L | Mat. No. 1.008.8333
 SONICflex LUX handpiece 2003 L | Mat. No. 1.000.4243
 SONICflex Set 2003 (without light) | Mat. No. 1.008.8334
 SONICflex handpiece 2003 (without light) | Mat. No. 1.000.4246

Prophylaxis

SONICflex™ 2000 N
 Frequency: 6,000 Hz - amplitude: 100µm - direction of movement: circular - internal water cooling - suitable for MULTIflex coupling - includes torque wrench and SONICflex scaler tips, 5° Inverted - 12 months warranty
 2000 N | Mat. No. 0.570.5020

Prophylaxis and Minimally invasive applications

PROPHYflex™ 4
 An air-polishing device for the gentle cleaning of tooth surfaces - all parts reproducible - perfect ergonomics - clean treatment - prevents disturbing powder mist - flexible working - for KaVo and Sirona couplings - 12 months warranty
 PROPHYflex 4 Standard | Mat. No. 3.002.8000
 PROPHYflex 4 S | Mat. No. 3.002.8000
 PROPHYflex 4 S Flamingo | Mat. No. 3.004.5950
 PROPHYflex 4 Wave | Mat. No. 3.002.8000
 PROPHYflex 4 S Wave | Mat. No. 3.002.8000
 PROPHYflex 4 S Wave | Mat. No. 3.004.5900
 PROPHYflex 4 L (me) | Mat. No. 3.002.2200
 PROPHYflex 4 S (me) | Mat. No. 3.004.5950

* This instrument fits a non KaVo coupling system (S = KaVo instrument fits to Sirona Click&Go-coupling)

PROPHYflexarts™
 Supragingival - spherical powder for safe removal of soft plaque - calcium carbonate neutral - 30 sticks (Mat. No. 1.000.2000)
 neutral - 4 bottles (250 g each) | Mat. No. 1.007.0014
 berry | Mat. No. 1.007.0015
 cherry | Mat. No. 1.007.0016
 mint | Mat. No. 1.007.0017

PROPHYflex Powder
 Supragingival - removal of soft plaque - sodium bicarbonate
 neutral - 4 bottles (250 g each) | Mat. No. 1.007.0014
 berry | Mat. No. 1.007.0015
 cherry | Mat. No. 1.007.0016
 mint | Mat. No. 1.007.0017

* 80 sticks (15 g each)

Perio Kit
 Used for subgingival treatments with the PROPHYflex - removes biofilm in periodontal pockets - cleans gingiva - improved soft rest cleaning | Mat. No. 1.000.3721
 Perio Kit | Mat. No. 1.000.3721
 PROPHYflex 4-Perio Kit | Mat. No. 1.011.9403
 * Tips can be reproduced up to 10 times

PROPHYflex Perio Powder
 Sub-/supragingival - implant cleaning - extra smooth - Tasteless - glycine powder
 PROPHYflex Perio Powder | Mat. No. 1.009.3732

DIAGNOcam™ Vision Full HD
 Three instant diagnostic images from one shot. Combining an intraoral camera, a transilluminating device for caries detection and a fluorescence instrument for root caries check into one device. All captured in a single shot and instantly displayed in Full HD resolution - 12 months warranty
 DIAGNOcam - Vision Full HD | Mat. No. 1.011.9213

DIAGNOcam™
 Radiation-free device for early detection of caries - ideal for patient education and excellent monitoring - 12 months warranty
 DIAGNOcam | Mat. No. 1.001.9000

DIAGNOdent™ pen and DIAGNOdent™ Display
 Device for early detection of hidden fissure and interproximal caries - detection of concretions in periodontal pockets - wireless instrument with optional external display - pen/Display 12 months warranty
 DIAGNOdent pen 1790 | Mat. No. 1.002.7700
 DIAGNOdent Display 2101 | Mat. No. 1.004.5400

Surgery portfolio

MASTERlung LUX Wireless Surgical Unit
 Large colour touchscreen (31" - wireless surgical foot control) - data documentation (SD memory card) - multi-program - integrated cooling pump - hose set - max. torque 80 Ncm - motor INTRA LUX 5600 LED brushless motor - speed range: 0 - 40,000 rpm, motor and motor cable: auto-storable & disinfectable - 12 months warranty
 MASTERlung LUX Wireless | Mat. No. 1.009.1200
 INTRA LUX 5600 LED | Mat. No. 1.008.8000

SURGmatic S16 L Pro speed increasing contra-angle handpiece 1:5
 Transmission ratio 1:5 with up to 200,000 rpm - all-rounder with high cutting performance for surgical and standard applications - 3-port spray - external cooling connection - Cool-Head technology - EXPERTgrip - 12 months warranty
 S16 L Pro (with light) | Mat. No. 1.014.4500

SURGmatic S201 L Pro Surgical contra-angle handpiece 2:0:1
 Reduction ratio 2:0:1 - Hexagon chuck system for a long service life and a max. torque of 80 Ncm - optimized clamping system - removable head, easy to clean - EXPERTgrip - small head with internal cooling system - external cooling connection - optimised toothing at the head and improved ball bearings - 12 months warranty
 S201 L Pro (with hexagon - light) | Mat. No. 1.013.7541

SURGmatic S201 L Pro Surgical contra-angle handpiece 2:0:1
 Reduction ratio 2:0:1 - max. torque of 55 Ncm - optimized clamping system - removable head, easy to clean - EXPERTgrip - small head with internal cooling system - external cooling connection - optimised toothing at the head and improved ball bearings - 12 months warranty
 S201 L Pro (with light) | Mat. No. 1.013.7540

INTRA Surgical handpiece 1:1
 Glass rod light conductor - for straight bars with 44.5 mm (2.34 in) standard length - external cooling connection - easy to dismantle - reinforced chuck system - 12 months warranty
 3610 N1 | Mat. No. 0.524.5600

INTRA Surgical handpiece 1:1
 Glass rod light conductor - for straight bars with 44.5 mm (2.34 in) standard length - external cooling connection - easy to dismantle - reinforced chuck system - 12 months warranty
 3610 N2 | Mat. No. 0.524.5610

INTRA Surgical handpiece 1:1
 For straight bars with 62.5 mm (2.5 in) special length - external cooling connection - easy to dismantle - reinforced chuck system - 12 months warranty
 3610 N2 | Mat. No. 0.524.5610

INTRA Surgical handpiece 1:1
 For straight bars with 70 mm (2.8 in) special length - external cooling connection - easy to dismantle - reinforced chuck system - 12 months warranty
 3610 N3 | Mat. No. 0.524.5620

Surgery portfolio

EXPERTlung LUX Surgical Unit
 Large colour display (31" - multifunctional foot control) (cable length 2.5 m) - single-user-operation for saving data - integrated cooling pump - hose set - max. torque 80 Ncm - motor INTRA LUX 5600 LED brushless motor - speed range: 0 - 40,000 rpm, motor and motor cable: auto-storable & disinfectable - 12 months warranty
 EXPERTlung LUX Surgical Unit | Mat. No. 1.008.3500
 DIAGNOcam 5600 LED | Mat. No. 1.008.3500

INTRA Surgical handpiece 1:1
 For straight bars with 44.5 mm (2.34 in) standard length - external cooling connection - easy to dismantle - reinforced chuck system - 12 months warranty
 3610 N1 | Mat. No. 0.524.5600

INTRA Surgical handpiece 1:1
 For straight bars with 44.5 mm (2.34 in) standard length - external cooling connection - easy to dismantle - reinforced chuck system - 12 months warranty
 3610 N2 | Mat. No. 0.524.5610

INTRA Surgical handpiece 1:1
 For straight bars with 62.5 mm (2.5 in) special length - external cooling connection - easy to dismantle - reinforced chuck system - 12 months warranty
 3610 N2 | Mat. No. 0.524.5610

INTRA Surgical handpiece 1:1
 For straight bars with 70 mm (2.8 in) special length - external cooling connection - easy to dismantle - reinforced chuck system - 12 months warranty
 3610 N3 | Mat. No. 0.524.5620

Surgery portfolio

SURtorque 459 L
 45 degree head angle - for 25 mm bars - 19 W - no air into bur direction, were out in 2 side vents - ceramic ball bearing - fiber optic - standard head size - push-button chuck - glass rod light conductor - single water jet - suitable for MULTIflex coupling
 459 L | Mat. No. 3.000.4808

SURtorque 459 C
 45 degree head angle - for 25 mm bars - 19 W - no air into bur direction, were out in 2 side vents - ceramic ball bearing - fiber optic - standard head size - push-button chuck - glass rod light conductor - single water jet - suitable for MULTIflex coupling
 459 C | Mat. No

Dieses Dokument wird Ihnen zur Verfügung gestellt von:

Med-Dent24 Handelsgesellschaft mbH
Gotenweg 7
63128 Dietzenbach

Telefon: +49 (0)6074 803 3615
Telefax: +49 (0)6074 803 4092
E-Mail: info@med-dent24.com

Registergericht: Offenbach am Main
Registernummer: HRB 56518
Umsatzsteuer-Identifikation-No.: DE307869540

Geschäftsführerin: Ana Ban