




1
Kugel
Round
Boule
6




Shank Size	L mm	ISO	Order No.
Size		ISO	Order No.
Ø ₁₀ mm			US No.
RA		310 204 001 001...	1.204...
		↻ _{max}	

005	006	007	008	009	010	012	014	016	018
¼	½		1		2	3	4	5	6
005	006	007	008	009	010	012	014	016	018
50.000									


1
Kugel
Round
Boule
6

Shank Size	L mm	ISO	Order No.
Size		ISO	Order No.
Ø ₁₀ mm			US No.
RA		310 204 001 001...	1.204...
		↻ _{max}	

021	023	025	027	029	031	033	035	037	040
7	8	9			11				
021	023	025	027	029	031	033	035	037	040
30.000			25.000			20.000			


2
Umgekehrter Kegel
Inverted Cone
Cône renversé
6

Shank Size	L mm	ISO	Order No.
Size		ISO	Order No.
Ø ₁₀ mm			US No.
RA		310 204 010 001...	2.204...
		↻ _{max}	

006	007	008	009	010	012	014	016	018	021	023
33½		34		35	36	37	38	39	40	41
006	007	008	009	010	012	014	016	018	021	023
50.000								30.000		


3
Rad
Wheel
Roue
6




Shank Size	L mm	ISO	Order No.
Size		ISO	Order No.
Ø ₁₀ mm			US No.
RA		310 204 040 001...	3.204...
		↻ _{max}	

006	007	008	009	010	012	014	016	018	021	023
11½		12			14		16			
006	007	008	009	010	012	014	016	018	021	023
50.000								30.000		

21
Zylinder
Cylinder
Cylindre
6




Shank Size	L mm	ISO	Order No.
Size		ISO	Order No.
Ø ₁₀ mm			US No.
RA		310 204 107 006...	21.204...
		↻ _{max}	

3,6	3,9	4,2	4,5	4,8	5,1	5,4	5,7	6,0
008	009	010	012	014	016	018	021	023
56		57	58	59	60	61	62	63
008	009	010	012	014	016	018	021	023
50.000					30.000			